




BEDROOM SAFETY & SUPPORT

The Aspire bedroom safety & support range has been designed with user support and safety in mind. Bed poles, cradles and rails ensure carers can provide safety and support for residents without compromising comfort in their bedroom environment.



Drop Down Bed Rail
BEA013700



Bed Pole With Adjustable Crook Handle
BEA014150



Free Standing Self Help Pole
BEA013100



Fixed Bed Rail
BEA013600



Clamp On Bed Pole
BEA009300 & BEA009310



Adjustable Height Bed Cradle
BEA014610



Fixed Height Bed Cradle
BEA014605



Clamp On Bed Rail
BEA829105

*Only fits community beds with a side frame thickness of at least 50mm



FREE STANDING SELF HELP POLE

BEA013100

- Free standing support handle provides assistance for the user to reposition or transfer in and out of bed
- Handle height is easily adjusted to suit the user



BED POLE WITH ADJUSTABLE CROOK HANDLE

BEA014150

- Adjustable crook handle can be located in multiple positions
- Base is positioned under the mattress



CLAMP ON BED POLE

BEA009300 (RIGHT) & BEA009310 (LEFT)

- Folding crook handle allows the pole to swing down and out of the way as required
- Clamps on to the side frame of community beds with a side frame thickness of at least 50mm



FIXED BED RAIL

BEA013600

- Support and safety rail provides repositioning and transfer support and assists to prevent falls from bed



DROP DOWN BED RAIL

BEA013700

- May be lowered down and out of the way for patient transfers or repositioning
- Base is positioned under the mattress



FIXED HEIGHT BED CRADLE

BEA014605

- Base is positioned under the mattress



ADJUSTABLE HEIGHT BED CRADLE

BEA014610

- Adjustable height suitable for a range of mattress heights
- Base is positioned under the mattress



CLAMP ON BED RAIL

BEA829105

- Enhances patient safety during repositioning, transfers, and transitions to and from mobility aids like wheelchairs
- The clamp design ensures consistent setup throughout different facilities