

DAY CHAIRS, STOOLS & COMMODES

Aspire for... Comfort, stability, support



HIGH BACK CLASSIC DAY CHAIR

A high back orthopaedic chair that offers adjustability, comfort and support.





Postural Support





Contoured High Backrest

With a slight recline angle provides superior postural support and opens hip positioning for comfort



Also Available in Slate Vinyl

MODELS

CHP208050

CHP208055

High Back Classic Day Chair - Champagne Vinyl High Back Classic Day Chair – Slate Vinyl

FEATURES & BENEFITS SHARED BY THE HIGHBACK CLASSIC AND WATERFALL DAY CHAIRS



Height Adjustable Legs With easy push button adjustment for comfortable positioning and changing needs



Sturdy Lightweight Design For a generous weight capacity, yet lightweight for easy manoeuvrability



Non-Slip Tips On chair legs provide stability whilst sitting or transferring



Robust Powder Coated Steel Frames Are durable and easy to clean

Seat Depth (mm)	Overall Depth (mm)	Seat Width (mm)	Overall Width (mm)	Seat Height (mm)	Backrest Height (mm)	Overall Height (mm)	Unit Weight (kg)	SWL (kg)
			Ī					kg
490	600	490	645	425–610	720	1160–1312	16	160

WATERFALL DAY CHAIR

A high back orthopaedic chair with contoured and adjustable backrest cushions that offer improved comfort and postural support.



Ink Vinyl



Latte Vinyl



Oatmeal Fabric



Mocha Fabric



Waterfall

Tri-Pillow Backrest Provides comfort and support

Removable Cushions

Seat and Back Cushion overlays can be easily detached for cleaning

Upholstery Range

Ink and Latte vinyl plus Oatmeal and Mocha fabric options



Height Adjustable



MODELS

CHP208070

CHP208070

Waterfall Day Chair – Ink Vinyl Waterfall Day Chair – Latte Vinyl CHP208075
CHP208076

Waterfall Day Chair – Oatmeal Fabric
Waterfall Day Chair – Mocha Fabric



Softly Padded Armrests, Seat & Back

Distribute pressure whilst simultaneously offering support surfaces for prolonged sitting comfort



Transport Legs With Wheels

Available as an accessory to facilitate easier manoeuvring of the chairs (without a patient in the chair)

See back cover for details



Easy CleanVinyl Upholstery

YEAR WARRANTY

✓ DVA CONTRACTED

AS / NZS Approved

Waterfall Day Chair is registered in the ARTG, AUST R ARTG 230087

High Back Classic Day Chair is registered in the ARTG, AUST R ARTG 230087



ADJUSTABLE & PRESSURE REDUCING DAY CHAIRS

A truly versatile height, width and depth adjustable high back orthopaedic chair with built in lumbar and lateral support and knee bolster. Available in pressure reducing variety.

ADJUSTABLE DAY CHAIR

Head Cushion

Height adjustable and removable head cushion for extra comfort

Contoured Seat and Backrest Padding

Designed for pressure redistribution and user comfort

Sturdy design

For a generous weight capacity of 160kg, yet lightweight for easy manoeuvrability





Also Available in Latte Vinyl



PRESSURE REDUCING ADJUSTABLE DAY CHAIR

Pressure Redistribution

Carrflex™ PU 4-Way Stretch Fabric for superior chemical resistance, retains flexibility and breathability

High quality Australian memory foam provides a comfortable and supportive surface, improving user comfort and pressure redistribution

Maximised Immersion and Envelopment to effectively redistribute pressure, improves blood flow and promotes healing





Available in 4-way stretch fabric

MODELS

CHP208025

CHP208030

Adjustable Day Chair – Latte Vinyl Adjustable Day Chair – Ink Vinyl CHP208035

Adjustable Day Chair - Pressure Reducing Blue Stretch Material

FEATURES & BENEFITS SHARED BY THE ADJUSTABLE AND MAXI ADJUSTABLE DAY CHAIRS



Height Adjustable Legs With easy push button adjustment for comfortable positioning

and changing needs



Tool Free Adjustment Points Allow the chair to be customised to suit long term use or readily adapted to suit multiple users



Adjustable Seat Depth, Width and Height Allows customised configuration to a user's needs and circumstances



Unique Width
Adjustable Frames
Cater to the full spectrum
of users shapes and sizes

Seat Depth (mm)	Overall Depth (mm)	Seat Width (mm)	Overall Width (mm)	Seat Height (mm)	Backrest Height (mm)	Overall Height (mm)	Unit Weight (kg)	SWL (kg)	
			Ī					kg	
475–590	760–610	510-610	575–675	410–530	740–790	1070–1220	26	160	

MAXI ADJUSTABLE DAY CHAIR

Bariatric Day Chair with 310kg weight capacity is fully height, width and depth adjustable to suit an array of large user shapes and sizes

Seat with Anterior Tilt

Facilitates standing from sitting and provides greater comfort for stiff or painful hips

Reinforced Design and Weight Capacity

Caters to a MAXI weight capacity of 310kg especially for bariatric users





Push Handle & Transport Wheels



Adjustable Seat Depth & Backrest Height



Height & Width Adjustable Frame





MODELS

CHP208010

Maxi Adjustable Day Chair – Burgundy Vinyl



Softly Padded Armrests

Distribute pressure whilst simultaneously offering support surfaces for prolonged sitting comfort




Push Handle and Transport Wheels

Allow a carer to tilt the chairs backwards onto the rear leg wheels in order to manoeuvre it about (without a patient in the chair)

See back cover for details



Easy CleanVinyl Upholstery



✓ DVA CONTRACTED

AS/ NZS Approved

ARTG, AUST R ARTG 230087

Adjustable Day Chair is registered in the ARTG, AUST R ARTG 230087



LOWBACK DAY CHAIR

A versatile, low back orthopaedic chair

PERCHING STOOL

A height adjustable stool that provides relief from fatigue in the kitchen or living area



Contoured and **Reclined Backrest**

Contour and slight recline of backrest angle helps maintain posture while supporting an open and comfortable hip position



Raised Armrests

Assist with sitting and standing

Seat with Anterior Tilt

Facilitates easier standing and sitting and provides greater comfort for stiff or painful hips



MODELS

CHP208000

CHP208005

Low Back Classic Day Chair – Champagne Vinyl Low Back Classic Day Chair - Slate Vinyl

MODELS

CHP213500

Perching Stool - Champagne Vinyl

FEATURES & BENEFITS SHARED BY THE LOWBACK DAY CHAIR AND PERCHING STOOL



Height Adjustable Legs Easy push button adjustment for comfortable positioning



Soft Padding

Distributes pressure and offers support for prolonged sitting comfort



Non-Slip Tips

Provides stability whilst sitting or transferring



Robust Powder Coated Steel Frames

Are durable and easy to clean



Easy Clean Vinyl upholstery



Lightweight For easy manoeuvrability





YEAR WARRANTY



	Seat Depth (mm)	Overall Depth (mm)	Seat Width (mm)	Overall Width (mm)	Seat Height (mm)	Backrest Height (mm)	Overall Height (mm)	Unit Weight (kg)	SWL (kg)
									kg
Low Back Chair	440	500	460	580	450–605	440	865–1015	9	160
Perching Stool	310	470	375	550	640-690	440	940–1090	6.8	175

CLASSIC BEDSIDE COMMODE

Allows for toileting in the bed room for users unable to travel between the bedroom and bathroom



Padded Armrests

Are raised to provide an easily reachable support point for lowering onto the seat

Lightweight

Bedside commode with padded seat, back and arms



Also Available in

DELUXE BEDSIDE COMMODE

Allows for toileting in the bed room with increased comfort for users unable to travel between the bedroom and bathroom



Contoured and Reclined Backrest

Contour and slight recline of backrest angle helps maintain posture while supporting an open and comfortable hip position

Softly Padded Armrests, Seat and Back

Offer support surfaces for prolonged sitting comfort



Also Available in Champagne Vinyl



MODELS

BEC045005

BEC045006

Classic Bedside Commode – Champagne Vinyl Classic Bedside Commode – Slate Vinyl

MODELS

BEC045116

BEC045115

Deluxe Bedside Commode – Champagne Vinyl Deluxe Bedside Commode – Slate Vinyl

FEATURES & BENEFITS SHARED BY THE CLASSIC AND DELUXE BEDSIDE COMMODES



Height Adjustable Legs
Easy push button adjustment
for comfortable positioning



Seat Access
Folds out of the way to access the toilet seat and commode bowl



Three-Sided Shroud
Discreetly conceals bowl



LightweightFor easy manoeuvrability




Removable Bowl
With lid and carry handle



Easy Clean Vinyl upholstery









	Seat Depth (mm)	Overall Depth (mm)	Seat Width (mm)	Overall Width (mm)	Seat Height (mm)	Backrest Height (mm)	Overall Height (mm)	Unit Weight (kg)	SWL (kg)
		\vdash							kg
Classic	440	570	450	570	490–615	455	915–1040	5.5	125
Deluxe	450	520	460	600	410–565	425	865–1015	12	160

aspire DAY CHAIRS AND STOOLS ACCESSORIES

TRANSPORT LEGS WITH WHEELS

FEATURES & BENEFITS

- · Assists carers to safely and easily tilt and manoeuvre an empty chair around a room
- · Wheeled legs insert directly into the frame to replace regular rear legs
- Transport wheels come as standard with the Adjustable day chair and Maxi Adjustable Day Chair and can be fitted as options on the Waterfall and High Back Classic day chairs





Transport Wheels CHP208090

COMPATIBLE CHAIRS

- CHP208050
- CHP208055
- CHP208070

High Back Classic Day Chair – Champagne Vinyl High Back Classic Day Chair – Slate Vinyl

Waterfall Day Chair – Ink Vinyl

CHP208071

CHP208075

CHP208076

Waterfall Day Chair - Latte Vinyl

Waterfall Day Chair - Oatmeal Fabric

Waterfall Day Chair – Mocha Fabric

PADDED LEGREST



CHP208080 CHP208085 Aspire Padded Legrest - Champagne Vinyl Aspire Padded Legrest – Slate Vinyl





FEATURES & BENEFITS



Comfortable **Resting Solution** Legs can be raised when

sitting in a day chair to provide extra comfort



Height Adjustable Easy push button adjustment for comfortable positioning

and changing needs



Robust Powder Coated Steel Frames Are durable and easy to clean



Padded Surface Generous support of legs and feet