
MAXI DAY CHAIR

Aspire for... Independence, Support & Comfort







FEATURES & BENEFITS

- Sturdy, flexible design for added stability and safety for the larger user, generous 310kg weight capacity.
- Ergonomic seat with anterior tilt to improve posture and reduce fatigue.
- Height, width and depth adjustability in seat and back for superior customised comfort.
- Easily manoeuvrable with push handle and transport wheels.
- Easy clean, breathable and comfortable vinyl upholstery and powder coated frame for effortless maintenance and longevity.
- Available in Burgundy.





1 YEAR WARRANTY

CODE	SEAT DEPTH	SEAT WIDTH	BACKREST	BACKREST	SEAT HEIGHT	BACKREST	PRODUCT
	RANGE	RANGE	WIDTH	HEIGHT RANGE	RANGE	HEIGHT	WEIGHT
CHP208010	390 - 530 mm	600 - 810 mm	600 mm	1120 - 1340 mm	370 - 520 mm	740mm	45 kg