Aspire Patient Positioning Products facilitate safe patient positioning in accordance with responsible manual handling techniques.

The range is designed with high quality fabrics and materials with attributes that enable gripping, sliding or swiveling to facilitate repositioning with minimum exertion and effort for patients, users and carers alike.

## **ASPIRE SLIDE SHEET**



- Low friction nylon fabric sheet is folded over and placed underneath the patient
- The top sheet slides freely on the bottom sheet, facilitating easy repositioning of the patient

Code	Description	Width	Length
TSE673070	Aspire Slide Sheet - Small	1450 mm	1000 mm
TSE673080	Aspire Slide Sheet - Large	1450 mm	2000 mm



## **ASPIRE SWIVEL CUSHION**



- Low friction cushion layers allow the user to effortlessly swivel in their seated position to assist maneuvering into or out of a vehicle seat
- Plush velour fabric with generous padding for maximum comfort on long car trips

Code	Description	Width
TSE673056	Aspire Swivel Cushion	440 mm





**ONE WAY SLIDE** 

- Assists in the correct positioning of a person in the back of a chair
- Helps maintain posture for the user, preventing them from slipping forward
- Reduces the risk of strain or injury from a fall or slip

Code	Description	Width	Length
TSE673050	One Way Slide	450mm	440mm



ONE WAY SLIDE FOOT OPERATED

- Non-slip mesh grips securely to smooth surfaces such as vinyl upholstery
- Soft seating surface can remain under patient throughout the day and facilitate multiple repositions
- Foot strap threads through back of chair to assist the carer with positioning

Code	Description	Width	Length
TSE673052	One Way Slide - Foot	450mm	1300mn



POSITIONING SLIDE FOOT OPERATED

- For patients with a lower risk of slipping while sitting in a chair
- Manoeuvre the patient easily to the back of the chair or sideways for transfer assistance
- Foot strap threads through back of chair to assist the carer with positioning

Code	Description	Width	Length
TSE673054	Positioning Slide	460mm	1375mm





