


aspire for... Quality Convenience Efficiency







## STACKABLE IV POLE

## Minimise Storage Space

The Aspire Stackable IV Pole is a unique design with interconnecting base frames that allow for IV Poles to be stored within their own footprint, reducing storage space.


## **Durable Base and Pole**

With a strong cast iron base, the Aspire Stackable IV Pole provides improved stability and durability whilst reducing manual handling efforts due to the lightweight 7kg overall weight.





Stackable IV Pole is registered in the ARTG, AUST R ARTG 303542



Common 24" Diameter Base: 4060mm

Extra Storage Space: 2450mm

Aspire Stackable IV Poles: 1610mm





Code	Description	Overall Height	Overall Width x Depth	Product Weight
WAT820185	Stackable IV Pole	1354 - 2214mm	637 x 518mm	7kg

