



# TILT & RECLINE MOBILE SHOWER COMMUNE



## Tilt-in-space

-5° to 35° of gas-assisted tilt-in-space functionality minimises manual handling efforts during shower



View our Aspire Shower Commode Accessory Flyer

- **30° of Backrest Angle Adjustment** in 10° increments for effortless repositioning and maximum user comfort
- **Open Armrests Style** accommodates a wider sitting surface whilst maintaining a smaller base footprint for improved accessibility
- **Tool-free Height Adjustment** to 4 settings allow the commode to be easily adapted to varying needs
- **4 x TENTE Total Locking Castors** provide brake access from all angles
- **Height Adjustable, Swing Away Footrests** with calf pads for optimal positioning, support and ease of transfers
- **Open Back Straps** provide greater access for bathing needs
- **Backrest May Be Configured** as full length or 3/4 length to fit an advanced headrest

Commode is registered in the ARTG AUST R ARTG 286695

**5** YEAR WARRANTY  
Corrosion

**1** YEAR WARRANTY  
Components



Code	Description	Height	Width	Seat Width	Backrest Height	Armrest Height	Anterior Seat Height	Posterior Seat Height
BTC066030	460mm Tilt & Recline	1315 - 1485mm	560mm	460mm	900mm	275mm	460 - 630mm	460 - 620mm
BTC066040	530mm Tilt & Recline	1315 - 1485mm	630mm	530mm	900mm	275mm	460 - 630mm	460 - 620mm