

WHEELCHAIR RANGE

aspire for...

Comfort
Posture
Mobility













EVOKE 2

REHAB RX

ASSIST 2

NSIT 2 LITE TRANSIT





aspire EVOKE 2 Lightweight & Highly Adjustable

The Aspire EVOKE 2 and EVOKE 2 HD are highly adjustable, lightweight, aluminium wheelchairs designed to enable a personalised fit for every user. A comprehensive range of compatible accessories is available for increased independence and enhanced clinical outcomes.



FEATURES

Frame: Folding Lightweight Aluminium

Legrest: Swing-away & Removable

Height & Angle Adjustable **Footrest:**

Armrest: Swing-back & Removable

Quick Release, Solid 24" Tyres, Wheels:

Height & Depth Adjustable

Castors: 8", Height & Angle Adjustable

Brakes: Easy Access Push-to-Lock Brakes

Depth & Height Adjustable Seat:

Backrest: Height, Tension, Angle & Depth Adjustable

Dacron Fire Retardant Fabric **Upholstery:**


Contoured Overlay is optional, **Cushion:** Standard cushion Included

Weight: From 14.8kg (SWL 140kg - HD 180kg)

EVOKE 2 pictured with optional EVOKE 2 pictured with Standard Seat Contoured Padded Overlay & Back without Overlay Customisable backrests Height adjustable for tension, angle, height back canes 17" - 19" and depth Quick release 24" solid rear wheels Swing-back armrests. (Height and depth adjustable option available) Depth and height adjustable seating Swing away and removable legrests Rear axle adjustability for amputees and users requiring different seat angles. leight and angle adjustable footrests Height and angle adjustable front castors 8"

OPTIONAL CONTOURED PADDED OVERLAY

This offers the user increased comfort and is ideal for those who are in their wheelchair for longer periods of time. The contouring assists in positioning and pressure distribution, while the fabric allows air movement, increasing breathability. It is easily removed for washing.



CUSTOMISABLE BACKREST

Tension Adjustable

✓ Angle Adjustable

Height Adjustable













SEAT WIDTH	CODE
350 mm (14")	MWS449500
400 mm (16")	MWS449510
450 mm (18")	MWS449520
500 mm (20")	MWS449530
570 mm (22")	MWS449540 (HEAVY DUTY model)
610 mm (24")	MWS449550 (HEAVY DUTY model)

OPTIONAL ACCESSORIES

- Elevating Legrest
- Advanced Articulating **Elevating Legrest**
- Stump Supports
- Oxygen Bottle Holder
- Brake Extensions
- Lateral Supports
- Contoured Padded Overlay
- Gutter Arm Rests*
- Drum Brakes

- Anti-tip Bars
- IV Pole
- Safety Belt Camber Kit
- Calf Straps
- Headrest Assembly
- Extended Back Canes 18" - 21"







ADJUSTABLE

*For Detailed Features, Specifications and Accessories see pages 14 - 20 Note: Side transfers are not possible when gutter armrests are attached


Aspire EVOKE 2 is crash tested to AS/NZS3696.19


aspire EVOKE 2 JNR


All the adjustable features of the EVOKE 2 have been reduced in size to suit the unique needs of children, including smaller 16 inch rear wheels.

CUSTOMISABLE BACKREST

- Depth Adjustable
- (V) Angle Adjustable
- (V) Height Adjustable









Frame:

Folding Lightweight Aluminium

Swing-away & Removable Legrest:

Height & Angle Adjustable Footrest:

Height Adjustable & Armrest:

Quick Release, Solid 20" Tyres Wheels:

8", Height & Angle **Castors:**

Adjustable

Brakes: Easy Access Push-to-lock

Height & Depth Adjustable Seat:

Height, Angle & Depth **Backrest:**

Adjustable

Contoured Overlay is Optional, **Upholstery:** Standard Cushion Included

Cushion: Included Weight: From 14kg

*For Detailed Features, Specifications and Accessories see pages 14 - 20













AS / NZS Approved	
oke JNR is registered in the	
RTG, AUST R ARTG 291670	Acni

SEAT WIDTH	CODE
300 mm (12")	MWS449560

Aspire EVOKE 2 JNR is crash tested to AS/NZS3696.19

aspire REHAB RX JNR

The Aspire REHAB RX JNR is designed for advanced postural support and pressure care for children via easy repositioning and advanced adjustability.

ADVANCED COMFORT Breathable PU coated multi-stretch fabric on seat and back cushions ensure optimal skin microclimate, comfort and improve pressure care outcomes.

TILT-IN-SPACE (TIS) Functionality promotes improved positioning. As well as this, it allows body areas.

CUSHIONING Multi-layer, castellated seat and back cushions with internal postural reinforcement optimise comfort and positioning.







FEATURES

Frame:

Legrest:

Footrest:

Armrest:

Wheels:

Castors:

Cushion:









TILT-IN-SPACE















Attendant Operated Drum Brakes **Brakes:** & Push-to-lock Handbrakes Depth & Width Adjustable Seat:

6", Height Adjustable

Swing-away, Removable, Padded

Removable, Height & Depth Adjustable

Height & Angle Adjustable

Height & Depth Adjustable

Quick Release, Solid 20" Tyres,

Backrest: Height, Tension & Angle Adjustable

Height, Depth, Angle & Width **Headrest:** Adjustable, Removable

Rigid Steel

& Elevating

PU Breathable Multi-stretch **Upholstery:** Fabric Covers

> Multilayer Castellated Back & Seat Cushions in Breathable PU

Weight: From 39.5kg

*For Detailed Features, Specifications and Accessories see pages 14 - 20

SEAT WIDTH	CODE
320 - 370 mm (13" - 15")	MWS449630

Aspire REHAB RX JNR is crash tested to AS/NZS3696.19















aspire REHAB RX Advanced Tilt-In-Space

The REHAB RX is geared towards the high needs, rehabilitation environment. Highly configurable and with pressure care cushioning surfaces, the REHAB RX is designed for the user with more complex



FEATURES

Seat:

Frame: Rigid Steel Frame

Swing-away, Removable, Padded Legrest:

& Elevating

Footrest: Height & Angle Adjustable

Removable, Height & Depth Adjustable **Armrest:**

Quick Release, Solid 24" Tyres. Wheels: Height & Depth Adjustable

6", Height Adjustable **Castors:**

Attendant Operated Drum Brakes **Brakes:**

& Push-to-lock Handbrakes

Depth. Width & Height Adjustable

Height, Tension & Angle Adjustable **Backrest:**

Height, Depth, Angle & Width **Headrest:**

Adjustable, Removable

PU Breathable Multi-Stretch Fabric Covers **Upholstery:**

Multilaver Castallated Back and Seat Cushion **Cushion:**

Covered in Breathable PU Fabric

Weight: From 41kg

*For Detailed Features, Specifications and Accessories see pages 14 - 20

OPTIONAL ACCESSORIES

A wide range of accessories to meet clinical needs

- Base Frame Extension
- 20" or 16" Alternative
- Wheels
- Therapist Adjustment Kit
- Stump Supports
- Safety Belt
- · Oxygen Bottle Holder
- Meal Tray
- Anti-tip Bars
- IV pole

HEAD SUPPORT Easy-to-adjust neadrest offers the user superior support & comfort. Breathable hygiene PU surfaces, easy to clean Swing-back armrests • CUSHION (height and depth adjustable option available) Quick release wheels, height & depth removable and padded Height & angle Height adjustable 6" adjustable footrests

THE HIGHLY ADJUSTABLE FEATURES FULFILL MANY CLINICAL REQUIREMENTS GEARED TOWARDS COMFORT AND SUPPORT

SOLUTION

pressure areas.

HIGH LEVEL REHABILITATION

Tilt-In-Space (TIS) functionality

promotes increased comfort and

improved positioning. As well as

redistribution to help alleviate

this, it allows body weight

SUPPORT SURFACES:

- Multi-layer foam construction with castellated top layer for enhanced comfort, pressure redistribution and durability
- Breathable PU multi-stretch cover for improved pressure redistribution and skin micro-climate management

BACKREST - CONTOURED:

- Gentle lumbar build up for postural support
- Moderate lateral supports for assisted midline sitting

CUSHION:

- Integrated positioning skeleton for enhanced sitting posture
- Fabric tensioners provide a smooth support surface
- Gentle abductor and adductor contouring provide improved postural support

ADJUSTABLE ELEVATING LEGRESTS

Elevating legrests with height & angle adjustable footrests.

INFECTION CONTROL Easy to clean surfaces.

CUSTOM COMFORT After-market backrests

& cushions ensure flexibility.



SAFETY IN DESIGN

Padded fittings and supports keep the user safe and secure at all times.

Aspire REHAB RX is crash tested to AS/NZS3696.19







SEAT WIDTH	CODE
410 - 460 mm (16" - 18")	MWS449710
460 - 510 mm (18" - 20")	MWS449720
510 - 560 mm (20" - 22")	MWS449740









BACKREST HEIGHT

adjustable

castors

SEAT DEPTH

ADJUSTABLE

aspire REHAB RS Classic Tilt-In-Space

The REHAB RS is an ideal solution for Aged Care and Community users who requires increased postural support and bodyweight redistribution. A reclined position is easily achieved, making the REHAB RS also a great option for those who fatigue easily.



FEATURES

Rigid Steel Frame Frame:

Swing-away, Removable, Padded & Elevating Legrest:

Footrest: Height & Angle Adjustable

Removable, Height & Depth Adjustable Armrest:

Quick Release, Solid 24" Tyres, Wheels:

Height & Depth Adjustable

Castors: 6", Height Adjustable

Attendant Operated Drum Brakes **Brakes:**

& Push-to-lock Handbrakes

Seat: Depth, Width & Height Adjustable

Height, Tension & Angle Adjustable **Backrest:**

Height, Depth, Angle & Width **Headrest:** Adjustable, Removable

Upholstery: PU Breathable Multi-Stretch Fabric Covers

Cushion: Standard Foam Comfort Cushions on Back & Seat

Weight: From 40.7kg

*For Detailed Features, Specifications and Accessories see pages 14 - 20

OPTIONAL ACCESSORIES

A wide range of accessories to meet clinical needs

- Base Frame Extension
- 20" or 16" Alternative
- Wheels
- Therapist Adjustment Kit
- Stump Supports
- Safety Belt
- Oxygen Bottle Holder
- Meal Tray
- Anti-tip Bars
- IV pole



The easy-to-adjust, supportive headrest offers the user superior head support, while the gently contoured backrest facilitates posture and positioning.

SUPPORTIVE & COMFORTABLE

INCREASED HEAD SUPPORT

Padded armrests, elevating legrests & a supportive headrest ensure that users of the REHAB RS are supported and comfortable.

ADVANCED COMFORT

Breathable PU coated multi-stretch fabric on seat and back cushions ensure optimal skin microclimate, comfort and improve pressure care outcomes.



EASY TO CLEAN

Wipe clean surfaces ensure infection control compliance.



SAFETY AT ITS CORE

Padded armrests. elevating legrests & a supportive headrest ensure that users of the REHAB RS remain safe and secure.



Aspire REHAB RS is crash tested to AS/NZS3696.19







SEAT	WIDTH	

410 - 460 mm (16" - 18")

460 - 510 mm (18" - 20")

510 - 560 mm (20" - 22")

MWS449620

CODE

MWS449600

MWS449610





THE HIGHLY ADJUSTABLE FEATURES FULFILL MANY CLINICAL REQUIREMENTS GEARED TOWARDS COMFORT AND SUPPORT



BACKREST HEIGHT







LOW TO MEDIUM LEVEL REHABILITATION SOLUTION

Tilt-In-Space (TIS) functionality promotes increased support and improved positioning. This allows caregivers to position the user for optimal comfort. Adjustable, elevating legrests allow for precise leg positioning.





DVA Contracted

aspire ASSIST 2 Everyday

The Aspire ASSIST 2 is a versatile and robust self-propelling wheelchair, designed to meet the rigorous demands of Community, Hire, Hospital and Aged Care environments. The Aspire ASSIST 2 is a reliable and comfortable solution with various accessories available to suit the needs of individuals and caregivers alike.





FEATURES

Seat:

Folding Steel Frame Frame: Swing-away, Removable, Legrest: (Elevating optional accessory)

Height Adjustable Footrest:

Armrest: Swing-back, Removable

Wheels: Quick Release, Solid 24" Tyres

8", Height Adjustable **Castors:**

Brakes: Push-to-lock

Fixed, Light Padding **Backrest:** Fixed, Light Padding

Upholstery: Dacron Fire Retardant Fabric

Cushion: Optional Accessory

Weight: From 17.3kg

ROBUST & RELIABLE

With a strong steel frame, removable swing-away legrests and quick release self-propelling rear wheels, the Aspire ASSIST 2 is a practical yet robust patient transport solution.

CATERS TO A WIDE VARIETY OF ENVIRONMENTS

The large, puncture-proof wheels ensure that the Aspire ASSIST 2 handles varied terrain with ease.

LOW MAINTENANCE FOLDING

Solid rear tyres and a centre cross brace design offer a low maintenance and portable mobility solution within a facility, hire or community environment.

SEAT WIDTH	CODE
350 mm (14")	MWS449300
400 mm (16")	MWS449310
450 mm (18")	MWS449320
500 mm (20")	MWS449330









*For Detailed Features, Specifications and Accessories see pages 14 - 20



The Aspire ASSIST 2 WARD Vinyl wheelchair has all the features of the ASSIST 2, but with durable hospital grade vinyl, suitable for the hospital & clinic environment. The ASSIST 2 WARD Vinyl has been carefully designed with fully welded seams and concealed mounting of upholstery to ensure that rigorous infection control standards can be met.





FEATURES

Frame: Folding Steel Frame Swing-away, Removable Legrest: (Elevating optional accessory)

Footrest: Height Adjustable

Armrest: Swing-back, Removable

Wheels: Quick Release, Solid 24" Tyres

8", Height Adjustable **Castors:**

Push-to-lock **Brakes:**

Fixed, Padded, Fire Retardant Vinyl Seat:

Fixed, Padded, Fire Retardant Vinyl

Upholstery: Hospital Grade Fire Retardant Vinyl

Cushion: Optional Accessory

Weight: From 17.3kg

Simple Frame Fix Mod to Stop

Security: Chair Being Folded (Available as option on request)

INFECTION CONTROL FOCUS

Durable hospital vinyl that is easy to wipe clean

Fully welded seams to inhibit bacterial growth

Concealed mounting of upholstery to prevent fluid ingress


400	/	~ - /	~ I II	- 11	- > /
ASS	⊢ I '	\sim		01	ıv
~~		ューヽ	_		

The centre cross-brace of the ASSIST 2 WARD can be customised to fix the frame for Hospital Asset Security.

1300 133 120

DURABLE TRANSPORT SOLUTION

The ASSIST 2 WARD Vinyl is a versatile transport wheelchair. Quick release wheels with solid, maintenance-free tyres provide on-call reliability and are suitable for use in common areas, corridors, wards, bathrooms and outdoors.



















aspireTRANSIT 2 Transporter

The Aspire TRANSIT 2 wheelchair is a versatile and robust transport wheelchair that is propelled by an attendant or carer. The Aspire TRANSIT 2 is a reliable and comfortable solution, with various accessories available to suit the needs of individuals and carers in the Community, Hospital or Aged Care environments.





FEATURES

Legrest:

Brakes:

Frame: Folding Steel Frame

Swing-away, Removable, (Elevating optional accessory)

Height Adjustable Footrest:

Armrest: Swing-away, Removable

Wheels: Quick Release, Solid 12" Tyres

Attendant Push-to-lock

Castors: 8", Height Adjustable

Seat: Fixed, Light Padding

Backrest: Fixed, Light Padding

Dacron Fire Retardant Fabric **Upholstery:**

Cushion: Optional Accessory

Weight: From 16.4kg

SWING-AWAY ARMRESTS

The padded armrests are swing-away and removable to facilitate the user when transferring in and out of the chair.

DURABLE TRANSPORT SOLUTION

The Aspire TRANSIT 2 is a versatile transport wheelchair, with rear 12 inch puncture-proof, transit wheels to allow an attendant or carer to easily navigate a client to their destination.

FOLDING DESIGN

The centre cross-brace design allows for the Aspire TRANSIT 2 to be folded up for unoccupied transport or storage within a facility, hire or community environment.

SEAT WIDTH	CODE
400 mm (16")	MWS449410
450 mm (18")	MWS449420
500 mm (20")	MWS449430











aspire LITE TRANSIT

The Aspire LITE TRANSIT is an lightweight wheelchair ideal for carer assisted transport. The compact lightweight frame makes it easy to transfer in and out of the car. Attendant brakes and 12 inch rear wheels assist the carer in negotiating indoor and outdoor terrain smoothly and safely. While the height adjustable armrests allow the user close access to tables or quick adjustment for comfort.



The Aspire LITE is an 'ultra lightweight' wheelchair that is ideal for occasional use. The Aspire LITE is designed with the carer in mind. The compact, lightweight frame makes it easy to transfer in and out of a car. The push button, height adjustable armrests make it versatile, allowing the user to be positioned at a table and quickly adjusted for comfort.



HEIGHT

ADJUSTABLE

ARMRESTS &

S.W.L

FEATURES

SAFETY

CHAIR

TRANSIT LITE

attendant brakes

Reduced manual

handling effort.

LIGHTWEIGHT

Folding Aluminium Frame Frame: Legrest: Swing-away & Removable

Armrest: Height Adjustable Wheels: Solid 12" Tyres **Castors:** 6" Height

Attendant Brakes **Brakes:** & Push-to-lock Brakes **Upholstery:** Padded Dacron Fabric

Weight: Chair Weight 11.6kg

SEAT WIDTH

460 mm (18")

TRANSPORTABLE Safety lap belt &

Easily folded due to it's centre cross brace design. The attendant handles can also be folded down, taking up minimal car space.

Frame:

Folding Aluminium Frame Legrest: Swing-away & Removable

Height Adjustable **Armrest:** Wheels: Solid 8" Tyres

Castors: 6" Height Push-to-lock **Brakes:**

Upholstery: Padded Dacron Fabric Chair Weight 9.5kg Weight:

TRANSPORTABLE SAFETY & STORAGE

Safety lap belt & rear storage pocket

ULTRA LIGHTWEIGHT At 9.5kg, allows for


also be folded down, taking

minimal manual handling.

Easily folded due to it's centre cross brace design. The attendant handles can

up minimal car space.

DVA Contracted







CODE

MWS442300

HEIGHT

ADJUSTABLE



	SUITS								
DESCRIPTION									B.
	Evoke 2	Evoke 2 JNR	Rehab RX	Rehab RX	Rehab RS	Assist 2	Transit 2	Lite Transit	Lite
Aluminium	•	JNR	JNR						
Steel									
Folding									
SEAT									
Depth adjustable									
Width adjustable									
Cushion included									
WHEELS									
Quick release									
Solid tyres									
24" wheels									
20" wheels									
12" wheels									
8" wheels									
LEGRESTS									
Swing-away									
Removable									
Padded									
Elevating									
ARMRESTS									
Swing-away									
Removable									
Height adjustable									
Depth adjustable									
BACKREST									
Height adjustable back canes									
Tension adjustable back caries									
BRAKES									
Drum Brakes									
Attendant Brakes									
OTHER FEATURES Self-propelled									
Attendant-propelled only									
Tilts									
Reclines									
Crash tested									

TECHNICAL SPECIFICATIONS

CODE	SEAT WIDTH	OVERALL WIDTH	SEAT DEPTH	SEAT HEIGHT	BACKREST HEIGHT	ARMREST HEIGHT	MINIMUM CHAIR WEIGHT	S.W.L
EVOKE 2				\$	†	†	KG	KG
MWS449500	350 mm (14")	560 mm	400 - 475 mm	440 - 480 mm	430 - 490 mm	205 - 225 mm	From 14.8 kg	140 kg
MWS449510	400 mm (16")	610 mm	400 - 475 mm	440 - 480 mm	430 - 490 mm	205 - 225 mm	From 14.9 kg	140 kg
MWS449520	450 mm (18")	660 mm	400 - 475 mm	440 - 480 mm	430 - 490 mm	205 - 225 mm	From 15.0 kg	140 kg
MWS449530	500 mm (20")	710 mm	400 - 475 mm	440 - 480 mm	430 - 490 mm	205 - 225 mm	From 15.1 kg	140 kg
MWS449540 (HD)	570 mm (22")	780 mm	400 - 475 mm	440 - 480 mm	430 - 490 mm	205 - 225 mm	From 16.0 kg	180 kg
MWS449550 (HD)	610 mm (24")	820 mm	400 - 475 mm	440 - 480 mm	430 - 490 mm	205 - 225 mm	From 16.2 kg	180 kg
EVOKE 2 JNR			-	\$	†	†	KG	KG
MWS449560	300 mm (12")	510 mm	300 - 375 mm	360 - 410 mm	430 - 490 mm	205 - 225 mm	14 kg	140 kg
REHAB RX JNR	□ □		₩	\$	\$	†	KG	KG
MWS449630	320 - 370 mm (13"- 15")	600 mm	440 - 490 mm	Front = 410 - 440 mm Rear = 400 - 475 mm	460 - 710 mm	310 - 460 mm	39.5 kg	180 kg
REHAB RX	□	□		↑	*	†	KG	KG
MWS449710	410 - 460 mm (16" - 18")	690 mm	440 - 530 mm	Front = 410 - 440 mm Rear = 400 - 475 mm	460 - 710 mm	310 - 460 mm	41 kg	180 kg
MWS449720	460 - 510 mm (18" - 20")	740 mm	440 - 530 mm	Front = 410 - 440 mm Rear = 400 - 475 mm	460 - 710 mm	310 - 460 mm	43 kg	180 kg
MWS449740	510 - 560 mm (20" - 22")	790 mm	470 - 560 mm	Front = 410 - 440 mm Rear = 400 - 475 mm	460 - 710 mm	310 - 460 mm	47.2 kg	180 kg
REHAB RS	□		→	*	\$	\$	KG	KG
MWS449600	410 - 460 mm (16" - 18")	690 mm	440 - 530 mm	Front = 410 - 440 mm Rear = 400 - 475 mm	460 - 710 mm	310 - 460 mm	40.7 kg	180 kg
MWS449610	460 - 510 mm (18" - 20")	740 mm	440 - 530 mm	Front = 410 - 440 mm Rear = 400 - 475 mm	460 - 710 mm	310 - 460 mm	42.5 kg	180 kg
MWS449620	510 - 560 mm (20" - 22")	790 mm	470 - 560 mm	Front = 410 - 440 mm Rear = 400 - 475 mm	460 - 710 mm	310 - 460 mm	46.7 kg	180 kg

e: Seat height for REHAB RS, RX, RX JNR, EVOKE 2 & EVOKE 2 JNR is without seat cushion. Standard seat cushion height is 100 mr

www.aspirecare.com.au 1300 133 120

TECHNICAL SPECIFICATIONS

CODE	SEAT WIDTH	OVERALL WIDTH	SEAT DEPTH	SEAT HEIGHT	BACKREST HEIGHT	ARMREST HEIGHT	MINIMUM CHAIR WEIGHT	S.W.L
ASSIST 2	→	→		\$	†	†	KG	KG
MWS449300	350 mm (14")	570 mm	400 mm	460 - 500 mm	440 mm	210 mm	17.3 kg	140 kg
MWS449310	400 mm (16")	620 mm	400 mm	460 - 500 mm	440 mm	210 mm	17.5 kg	140 kg
MWS449320	450 mm (18")	670 mm	400 mm	460 - 500 mm	440 mm	210 mm	17.7 kg	140 kg
MWS449330	500 mm (20")	720 mm	400 mm	460 - 500 mm	440 mm	210 mm	17.9 kg	140 kg
MWS449340 (DELUXE VINYL)	350 mm (14")	570 mm	400 mm	460 - 500 mm	440 mm	210 mm	17.4 kg	140 kg
MWS449350 (DELUXE VINYL)	400 mm (16")	620 mm	400 mm	460 - 500 mm	440 mm	210 mm	17.6 kg	140 kg
MWS449360 (DELUXE VINYL)	450 mm (18")	670 mm	400 mm	460 - 500 mm	440 mm	210 mm	17.8 kg	140 kg
MWS449370 (DELUXE VINYL)	500 mm (20")	720 mm	400 mm	460 - 500 mm	440 mm	210 mm	18 kg	140 kg
MWS449380 (HD DELUXE VINYL)	570 mm (22")	790 mm	400 mm	460 - 500 mm	440 mm	210 mm	18.9 kg	180 kg
MWS449390 (HD DELUXE VINYL)	610 mm (24")	830 mm	400 mm	460 - 500 mm	440 mm	210 mm	19.2 kg	180 kg
TRANSIT 2	→	→		‡	†	‡	KG	KG
MWS449410	400 mm (16")	540 mm	400 mm	460 - 500 mm	440 mm	210 mm	16.4 kg	140 kg
MWS449420	450 mm (18")	590 mm	400 mm	460 - 500 mm	440 mm	210 mm	16.6 kg	140 kg
MWS449430	500 mm (20")	640 mm	400 mm	460 - 500 mm	440 mm	210 mm	16.8 kg	140 kg
LITE TRANSIT				***	†	†	KG	KG
MWS442301	460 mm (18")	600 mm	400 mm	460 mm	465 mm	235 mm	11.6 kg	120 kg
LITE				\$	‡	†	KG	KG
MWS442300	460 mm (18")	550 mm	400 mm	460 mm	465 mm	235 mm	9.5 kg	120 kg

ACCESSORIES









BASE FRAME EXTENSION KITS

STUMP SUPPORT

ELEVATING LEGRESTS










OXYGEN BOTTLE BAG









CONTOURED PADDED **OVERLAY**





GUTTER ARMRESTS







SAFETY LAP BELT

BRAKE EXTENSION LEVERS

LATERAL CHEST SUPPORT

THERAPIST-FRIENDLY, TOOL-FREE ADJUSTMENT KIT



IV POLE







1300 133 120

ACCESSORIES MATRIX

			SUITS							
CODE	DESCRIPTION		Evoke 2	Evoke 2 JNR	Rehab RX	Rehab RX	Rehab RS	Assist 2	Transit 2	
ELEVATING L	EGRESTS									
WCA848100	Elevating Legrest - Left - 300/320/350/400/410mm	A	•	•	0	0	0	•	•	
WCA848110	Elevating Legrest - Right - 300/320/350/400/410mm		•	•	0	0	0			
WCA848120	Elevating Legrest - Left - 450/460/500/510mm		•			0	0	•	•	
WCA848130	Elevating Legrest - Right - 450/460/500/510mm		•			0	0	•		
WCA848140	Elevating Legrest - Left - 570/610mm		•					•		
WCA848150	Elevating Legrest - Right - 570/610mm		•							
STUMP SUPPO	DRTS									
WCA848220	Stump Support - Left - Silver - 300/320/350/400/450mm		•	•	•	•	•	•	•	
WCA848230	Stump Support - Right - Silver - 300/320/350/400/450mm					•				
WCA848240	Stump Support - Left - Silver - 500/550/570/610mm		•			•	•		•	
WCA848250	Stump Support - Right - Silver - 500/550/570/610mm	Th	•			•				
WCA848260	CoG Extension Bracket - for Stump Support		•	•						
OXYGEN BOT	OXYGEN BOTTLE HOLDERS									
WCA848290	Oxygen Bottle Holder	1	•	•						
WAA695600	Oxygen Bottle Bag				•	•	•			
BRAKE EXTER	NSIONS									
WCA848310	Brake Extension Lever - Left		•	•						
WCA848320	Brake Extension Lever - Right	7	•	•						
LATERAL SUF	PPORTS									
WCA848430	Lateral Chest Support - Left	4	•	•						
WCA848440	Lateral Chest Support - Right	1								
WCA848450	Lateral Chest Support - Left					•				
WCA848460	Lateral Chest Support - Right									

ACCESSORIES MATRIX

			SUITS					
CODE	DESCRIPTION	Evoke 2	Evoke 2 JNR	Rehab RX JNR	Rehab RX	Rehab RS	Assist 2	Transit 2
CONTOURED	PADDED OVERLAYS							
WCA848330	Contoured Padded Overlay - 350mm Wide	•						
WCA848340	Contoured Padded Overlay - 400mm Wide							
WCA848350	Contoured Padded Overlay - 450mm Wide	•						•
WCA848360	Contoured Padded Overlay - 500mm Wide							
WCA848370	Contoured Padded Overlay - 570mm Wide	•						
WCA848380	Contoured Padded Overlay - 610mm Wide	•						
DEPTH & HEIG	HT ADJUSTABLE ARMRESTS							
WCA848390	Depth & Height Adjustable Armrest - Left	•						
WCA848400	Depth & Height Adjustable Armrest - Right	•						
GUTTER ARM	RESTS *							
WCA848410	Gutter Armrest Pad - Left	•						
WCA848420	Gutter Armrest Pad - Right							
DRUM BRAKE	s							
WCA848510	Rear Wheels w/Drum Brakes Kit - 16in		•	•	•			
WCA848520	Rear Wheels w/Drum Brakes Kit - 20in		•	0				
WCA848500	Rear Wheels w/Drum Brakes Kit - 24in	•			0	0		
DEPTH EXTEN	ISIONS							
WCA848530	Base Frame Extension Kit - 50mm - Black							
WCA848540	Base Frame Extension Kit - 100mm - Black							
OTHER								
WCA848550	Therapist-Friendly Tool-Free Adjustment Kit	•		•	•			
WCA848560	Meal Tray - 320mm to 410mm							
WCA848570	Meal Tray - 460mm to 510mm				•			
WCA848270	Anti-Tip Bars	•	•	0	0	0	•	
WCA848280	IV Pole	•	•	•	•	•		
WCA848300	Safety Lap Belt - Metal Seatbelt Buckle		•				•	
WCA848470	Extended Back Cane - Pair - 19/21in	•	•					
WCA848480	Headrest Assembly - Deluxe Style	•	•					
WCA848580	Calf Strap	•	•	•	•	•		
WCA848490	Camber Kit	•	•					
AAL848590	Security Mod - ASSIST 2 WARD Only							

O - Standard Feature Included *Note: Side transfers are not possible when gutter armrests are attached



aspire Associated Products

































Distributed by:

ASBR20.69.2 010621