


# RELAX BED

The Relax Bed has a fixed base and combines luxurious finishes with multiple positioning and therapeutic functions to support users and carers alike.





### **FEATURES & BENEFITS**

- **Night Assist:** Luminous underbed halo for safe entry and exit throughout the night.
- **Massage:** Relaxing vibration modes for personalised neck, shoulder and back relief.
- Press & Hold Safety Handset: Intuitive, large button interface enhances usability.
- Box-on-Box Design: Facilitates easier entry and exit, allowing true flush positioning of beds in companion or dual adjustable bed configurations.
- **German Engineered Motors:** For smooth, quiet adjustment.
- **Personalised Positioning:** Hand control allows users to find their perfect position to maximise comfort and relaxation.
- **USB Charge:** Connect and charge your devices conveniently.
- **Height Adjustable Legs:** Allows the bed height to be adjusted to 320, 370, 380 or 430mm tall. No tools are needed to make the adjustment.

### **FUNCTIONS**

- **Head Tilt:** Relaxes neck muscles, providing optimal support whilstreading or watching TV.
- Back Raise: Improves postural support and can ease symptoms of reflux, respiratory conditions, snoring and sleep apnoea.
- **Leg Raise:** Promotes circulation, relieving tired, swollen and/or painful limbs.
- **Zero Gravity:** Helps to relax the circulatory system and muscles of the lower back.
- Memory Function: Allows the user to set and return to their preferred positions at any time with one click of a button.











OVERALL WIDTH (mm)	OVERALL LENGTH (mm)	OVERALL HEIGHT (mm)
<b>←</b> →		<b>\$</b>
920	2030	Height adjustable to 320, 370, 380 or 430
1070		
1370		
1520		

#### **ACCESSORIES**

### ALLURE HEADBOARD



- · Stylish headboard to fit ComfiMotion bed collection.
- · Subtle square padded design.

Long single King single Long double Queen King

#### MEMORY FOAM MATTRESS



- · Integrated side walls.
- · Contoured Immersion.
- · Superior breathability.
- · Anti- microbial.

Classic - Long Single Classic – King Single Classic – Long Double Classic - Queen Incontinence – Long Single Incontinence – King Single Incontinence - Long Double Incontinence - Queen

### CHESTERFIELD HEADBOARD



- · Stylish headboard to fit ComfiMotion bed collection.
- · Chesterfield padded design.

Long Single King Single Long Double Queen King

#### **POCKET SPRING MATTRESS**



- · Weight activated support zones that contour to individual body shapes to improve comfort, support and sleep quality.
- · Integrated side walls.
- · Luxuriously soft pillowtop layer.

Long Single King Single Long Double Queen

#### CROOK HANDLE BED STICK



- · Provides multiple grip points to assist repositioning in bed, and safe entry and exit.
- · Designed to be used inconjunction with ComfiMotion bed collection.



- · Improves bed edge awareness and offers multiple grip points for repositioning and bed mobility.
- · Designed to be used inconjunction with ComfiMotion bed collection.

Plush Fibrefill sleep surface and

### **BREEZE PILLOW**



- Surface Cooling Technology helps regulate your body temperature for a more comfortable sleep experience.
- Medium Density, Immersive Memory Foam.

## PLUSH PILLOW



- Memory Foam core provides ergonomic support and comfort.
  - Medium density.