aidacare HEALTHWEIGH CHAIR SCALE







The Healthweigh Chair Scale allows for carers to measure the weight of residents in a safe and comfortable manner. The ergonomically contoured seat and folding armrests ensure patient comfort and ease of transfers whilst the large LCD screen provides accurate weight measurements for carers.

1	Movement compensation technology					
2	BMI function					
3	Fold back armrests					
4	Rear locking castors					

5	Automatic Reweigh		
6	Ergonomic contoured seating surface		
7	Auto power off		
8	Rechargeable battery		

	SEAT WIDTH	OVERALL WIDTH	SEAT DEPTH	OVERALL DEPTH	OVERALL HEIGHT	UNIT WEIGHT	S.W.L	
SS377706 ealthweigh Chair Scale	430mm	630mm	420mm	850mm	965mm	18kg	300kg	

