

PRESSURE SAFETY UNDERLAY

Specifically for use underneath an Alternating Pressure Overlay the Pressure Safety Underlay is designed to provide 'bottoming out' protection while minimising the overall mattress height on the bed.

The foam underlay can be rotated in four directions to maximise longevity and performance.







PRESSURE SAFETY UNDERLAY



Used underneath Alternating Pressure Overlays

Provides bottoming out protection for significant pressure points or in the event of an overlay pressure failure



Low Height

Keeps combined underlay and mattress surface height low, reducing the potential severity of falls and helping to maintain IGCG0G01-2-52 side rail safety compliance



CARRFlex Bi-Elastic Cover (European made)

Anti microbial, anti bacterial, fire retardant, stretch characteristics reduce skin shear



Four Way Rotation

Rotation ensures longevity of foam performance



Fully Welded Cover

Welded seems with water resistant zippers and reinforced waterfall flaps provide optimum infection control.

SPECIFICATIONS

Code	BEM047050
Description	Lifecomfort Pressure Safety Underlay
Length	2000mm
Width	900mm
Height	50mm
Warranty	5 years foam, 2 years cover
SWL	135kg

CARE

Removable and fully washable cover up to 95°c

Wipe towards mattress centre to avoid spread of bacteria

Do Not use solvents or non-approved bleaches

Do Not iron

Do Not twist or wring out



Australian Made