

# aspire LIFT RECLINE CHAIRS RANGE COMPARISON



## GOLD RANGE

	RAPHAEL				DA VINCI	REMBRANDT	MATTISSE	RENOIR			MONET
Description	4 Motor – Quattro Action				4 Motor – Quattro Action	4 Motor – Quattro Action	4 Motor – Quattro Action	4 Motor – Quattro Action			2 Motor – Dual Action Zero Gravity
Size	Small	Medium	Large	Maxi	Large	Maxi	Large	Small	Medium	Large	Large
Codes	1. CHP227290 2. CHP227291	1. CHP227300 2. CHP227301 3. CHP227302	1. CHP227310 2. CHP227311 3. CHP227312	1. CHP227319 2. CHP227322	4. CHP227500 5. CHP227510	6. CHP227550	7. CHP227560	8. CHP227575	8. CHP227580	8. CHP227585	9. CHP227540
Seat Height	470mm	470mm	500mm	510mm	480mm	495mm	500mm	450mm	470mm	470mm	495mm
Seat Depth	480mm	500mm	530mm	540mm	500mm	535mm	546mm	510mm	530mm	530mm	535mm
Seat Width	480mm	520mm	550mm	650mm	530mm	560mm	571mm	510mm	550mm	590mm	560mm
Overall Height	1060mm	1090mm	1080mm	1050mm	1050mm	1095mm	1067mm	1010mm	1060mm	1060mm	1065mm
Overall Width	740mm	780mm	810mm	930mm	870mm	1015mm	965mm	750mm	800mm	830mm	890mm
Backrest Height	630mm	680mm	700mm	650mm	630mm	600mm	730mm	660mm	700mm	700mm	570mm
SWL	130kg	130kg	130kg	180kg	159kg	159kg	159kg	130kg	130kg	130kg	159kg
Warranty	2 Years	2 Years	2 Years	2 Years	2 – 5 Years*	2 – 5 Years	2 – 5 Years	2 Years	2 Years	2 Years	2 – 5 Years



## SILVER RANGE

## BRONZE RANGE

	SIGNATURE			PORTO	MAUI	OREGON	IDAHO	MONTANA	MONTANA HD
Description	1 Motor – Single Action Space Saver		2 Motor – Dual Action	1 Motor – Single Action Zero Gravity	2 Motor – Dual Action	1 Motor – Single Action	2 Motor – Dual Action	1 Motor – Single Action	2 Motor – Dual Action
Size	Small	Medium	Large	Petite	Medium	Medium	Medium	Maxi	HD
Code	10. CHP228100 – Single 10. CHP228200 – Dual 11. CHP228110 – Single 11. CHP228210 – Dual 12. CHP228120 – Single 12. CHP228220 – Dual	10. CHP228130 – Single 10. CHP228230 – Dual 11. CHP228140 – Single 11. CHP228240 – Dual 12. CHP228150 – Single 12. CHP228250 – Dual	10. CHP228160 – Single 10. CHP228260 – Dual 11. CHP228170 – Single 11. CHP228270 – Dual 12. CHP228180 – Single 12. CHP228280 – Dual	13. CHP227200	14. CHP227110	16. CHP227150 17. CHP227155	16. CHP227160 17. CHP227165	1. CHP227400	18. HD300 – CHP227405 18. HD450 – CHP227410
Seat Height	480mm	480mm	500mm	490mm	500mm	480mm	480mm	510mm	500mm
Seat Depth	460 – 485mm	490 – 515mm	510 – 535mm	500mm	500mm	480mm	480mm	540mm	500mm
Seat Width	480mm	510mm	540mm	490mm	540mm	480mm	480mm	650mm	650mm
Overall Height	980mm	1050mm	1110mm	1120mm	1100mm	1100mm	1100mm	1050mm	1270mm
Overall Width	740mm	780mm	780mm	770mm	780mm	800mm	800mm	930mm	900mm
Backrest Height	580mm	620mm	650mm	690mm	630mm	620mm	620mm	650mm	810mm
SWL	145kg	145kg	155kg	150kg	130kg	150kg	150kg	180kg	300kg and 450kg
Warranty	2 Years	2 Years	2 Years	2 Years	2 Years	2 Years	2 Years	2 Years	2 Years

## UPHOLSTERY & SWATCHES


Numbering before codes refers to colour - see swatches

\* 2 Years Handset /Upholstery and Electronics,  
5 Years Frame

Tested to International Standards, DVA Contracted

1300 133 120 | aidacare.com.au