

COMPARE OUR LIFTER RANGE



MODEL	A150F	200	A205	AURA	A320	200S	A200S	NOVA	STAR
Lifter Type	Folding Floor Lifter	Floor Lifter	Floor Lifter	Floor Lifter	Floor Lifter	Standing Lifter	Standing Lifter	Standing Lifter	Patient Mover
Frame Material	Aluminium	Steel	Aluminium	Aluminium	Aluminium	Steel	Aluminium	Aluminium	Aluminium
Product Code	LSS390600	LSS390620	LSS390610	LSS390615	LSS390630	LSS390650	LSS390640	LSS390652	LSS390580
Safe Working Load – SWL	150kg	200kg	205kg	250kg	320kg	200kg	200kg	250kg	185kg
Product Weight	32.9kg	40kg	38.1kg	38kg	63kg	40kg	33kg	39.6kg	40kg
Base Length	1280mm (1220mm folded)	1320mm	1238mm	1378mm	1500mm	1277mm	1070mm	1100mm	925mm
External Base Width	660 - 1040mm (615mm folded)	665 - 1120mm	651 - 1173mm	682 - 1300mm	765 - 1190mm	656 - 1203mm	590 - 1040mm	682 - 1138mm	520mm
Internal Base Width	580 - 950mm	580 - 1040mm	547 - 1073mm	518 - 1150mm	640 - 1065mm	561 - 1103mm	490 - 940mm	519 - 986mm	740mm
Base Width Adjustment	Manual	Electric	Electric	Electric	Electric	Electric	Electric	Electric	Manual
Base Height / Under Bed Clearance	133mm	110mm	110mm	131mm	155mm	110mm	110mm	136mm	60mm
Minimum Height of Spreader Bar **	428mm	520mm	415 - 515mm	437mm	540 - 640mm	857mm	880mm	692mm	-
Maximum Height of Spreader Bar**	1640mm	1720mm	1700 - 1800mm	1788mm	1840 - 1940mm	1777mm	1415mm	1770mm	-
Under Base Height	72mm	50mm	50mm	49mm	155mm	50mm	50mm	53mm	-
Maximum Overall Height	1873mm (395mm folded)	1935mm	2040mm	2048mm	2155mm	1832mm	1480mm	1800mm	1030mm
Front Castor Diameter	75mm	75mm	75mm	100mm	100mm	75mm	75mm	100mm	100mm
Rear Locking Castor Diameter	100mm	100mm	100mm	125mm	125mm	100mm	75mm	125mm	100mm
Turning Radius	1365mm	1398mm	1332mm	1453mm	1595mm	1361mm	1160mm	1245mm	-
Heaviest Part Weight – Base Assembly	32.9kg	23.0kg	19.6kg	19.6kg	38.9kg	19.6kg	19kg	17.w6kg	-
Handle Height	1170mm	1220mm	1085 - 1185mm	883 - 1305mm	1200 - 1300mm	1045mm	1160mm	810 - 1095mm	1030mm
Warranty – Frame	5 Years	5 Years	5 Years	5 Years	5 Years	5 Years	5 Years	5 Years	2 Years
Warranty – Electronics	3 Years	3 Years	3 Years	3 Years	3 Years	3 Years	3 Years	3 Years	-
Warranty – Battery Handset & Cables	1 Year	1 Year	1 Year	1 Year	1 Year	1 Year	1 Year	1 Year	-

ACCESSORIES

In-line Weigh Scale – SWL 320kg	LSS390680	LSS390680	LSS390680	Not Available	LSS390680	Not Available	Not Available	Not Available	Not Available
Integrated Weigh Scale – SWL 220kg	LSS390675	LSS390675	LSS390675	Available in Q3 2025	LSS390675	Not Available	Not Available	Not Available	Not Available
Manual Pivot Frame – SWL 227kg	Not Available	LSS390670	LSS390670	Not Available	LSS390670	Not Available	Not Available	Not Available	Not Available
Powered Pivot Frame – SWL 227kg	Not Available	LSS390665	LSS390665	Available in Q3 2025	LSS390665	Not Available	Not Available	Not Available	Not Available
4 Point Sling Bar	Not Available	Not Available	Not Available	LSS361000	Standard Inclusion	Not Available	Not Available	Not Available	Not Available
Clip Sling Bar Attachment	Not Available	Not Available	Not Available	Not Available	Not Available	LSS390655	Not Available	LSS390660	Not Available
Leg Protection Strips / Shrouds*	Not Available	GSP390816	GSP390815	GSP242902*	GSP390817	Not Available	GSP390815	GSP242806*	Not Available

**Three mast height settings are possible hence the range of minimum and maximum heights. | To view our comprehensive list of spare parts visit aidacare.com.au