

PATIENT LIFTER RANGE

Mobile patient lifters and slings designed for a range of environments, from community to facility.



COMPARE OUR LIFTER RANGE



MODEL	A150F	200	A205	AURA
Lifter Type	Folding Floor Lifter	Floor Lifter	Floor Lifter	Floor Lifter
Frame Material	Aluminium	Steel	Aluminium	Aluminium
Product Code	LSS390600	LSS390620	LSS390610	LSS390615
Safe Working Load – SWL	150kg	200kg	205kg	250kg
Product Weight	32.9kg	40kg	38.1kg	38kg
Base Length	1280mm (1220mm folded)	1320mm	1238mm	1378mm
External Base Width	660 - 1040mm (615mm folded)	665 - 1120mm	651 - 1173mm	682 - 1300mm
Internal Base Width	580 - 950mm	580 - 1040mm	547 - 1073mm	518 - 1150mm
Base Width Adjustment	Manual	Electric	Electric	Electric
Base Height / Under Bed Clearance	133mm	110mm	110mm	131mm
Minimum Height of Spreader Bar **	428mm	520mm	415 - 515mm	437mm
Maximum Height of Spreader Bar**	1640mm	1720mm	1700 - 1800mm	1788mm
Under Base Height	72mm	50mm	50mm	49mm
Maximum Overall Height	1873mm (395mm folded)	1935mm	2040mm	2048mm
Front Castor Diameter	75mm	75mm	75mm	100mm
Rear Locking Castor Diameter	100mm	100mm	100mm	125mm
Turning Radius	1365mm	1398mm	1332mm	1453mm
Heaviest Part Weight – Base Assembly	32.9kg	23.0kg	19.6kg	19.6kg
Handle Height	1170mm	1220mm	1085 - 1185mm	883 - 1305mm
Warranty – Frame	5 Years	5 Years	5 Years	5 Years
Warranty – Electronics	3 Years	3 Years	3 Years	3 Years
Warranty – Battery Handset & Cables	1 Year	1 Year	1 Year	1 Year

ACCESSORIES

In-line Weigh Scale – SWL 320kg	LSS390680	LSS390680	LSS390680	Not Available
Integrated Weigh Scale – SWL 220kg	LSS390675	LSS390675	LSS390675	LSS390690 (250kg)
Manual Pivot Frame – SWL 227kg	Not Available	LSS390670	LSS390670	Not Available
Powered Pivot Frame – SWL 227kg	Not Available	LSS390665	LSS390665	LSS391250 (250kg)
4 Point Sling Bar	Not Available	Not Available	Not Available	LSS361000
Clip Sling Bar Attachment	Not Available	Not Available	Not Available	Not Available
Leg Protection Strips / Shrouds*	Not Available	GSP390816	GSP390815	GSP242902*

**Three mast height settings are possible hence the range of minimum and maximum heights. | To view our comprehensive list of spare parts visit aidacare.com.au



A320	200S	A200S	NOVA	STAR
Floor Lifter	Standing Lifter	Standing Lifter	Standing Lifter	Patient Mover
Aluminium	Steel	Aluminium	Aluminium	Aluminium
LSS390630	LSS390650	LSS390640	LSS390652	LSS390580
320kg	200kg	200kg	250kg	185kg
63kg	40kg	33kg	39.6kg	40kg
1500mm	1277mm	1070mm	1100mm	925mm
765 - 1190mm	656 - 1203mm	590 - 1040mm	682 - 1138mm	520mm
640 - 1065mm	561 - 1103mm	490 - 940mm	519 - 986mm	740mm
Electric	Electric	Electric	Electric	Manual
155mm	110mm	110mm	136mm	60mm
540 - 640mm	857mm	880mm	692mm	-
1840 - 1940mm	1777mm	1415mm	1770mm	-
155mm	50mm	50mm	53mm	-
2155mm	1832mm	1480mm	1800mm	1030mm
100mm	75mm	75mm	100mm	100mm
125mm	100mm	75mm	125mm	100mm
1595mm	1361mm	1160mm	1245mm	-
38.9kg	19.6kg	19kg	17.w6kg	-
1200 - 1300mm	1045mm	1160mm	810 - 1095mm	1030mm
5 Years	5 Years	5 Years	5 Years	2 Years
3 Years	3 Years	3 Years	3 Years	-
1 Year	1 Year	1 Year	1 Year	-
LSS390680	Not Available	Not Available	Not Available	Not Available
LSS390675	Not Available	Not Available	Not Available	Not Available
LSS390670	Not Available	Not Available	Not Available	Not Available
LSS390665	Not Available	Not Available	Not Available	Not Available
Standard Inclusion	Not Available	Not Available	Not Available	Not Available
Not Available	LSS390655	Not Available	LSS390660	Not Available
GSP390817	Not Available	GSP390815	GSP242806*	Not Available



Weigh Scale Sling Bar
LSS390690

Tailored for the Aura
In use with general-purpose slings featuring loop attachments.

Optimised Lifting Performance
By reducing height loss compared to conventional hoist scale systems.

High Weight Capacity
Up to 250kg.



In-line Weigh Scale
LSS390680

Easy-to-read
LCD display.

Diverse Compatibility
Designed to be fitted on A150F, A205, 200 and A205.

High Weight Capacity
Up to 320kg.



Integrated Weigh Scale
LSS390675

Easy-to-read
LCD display.

Diverse Compatibility
Compatible with A150F, 200 and A205, in use with general-purpose slings featuring loop attachments.

Optimised Lifting Performance
By reducing height loss compared to conventional hoist scale systems.

Infection Control
With non-absorbent polyurethane padding designed for easy cleaning.



Powered Pivot Frame
LSS391250

Tailored for the Aura
In use with slings featuring clip attachments.

Minimises Carer Effort
Assists staff to conduct repetitive transfers in institutional care environments.

Integrated Electronics
Seamlessly integrates with lifter control system for use with one handset.



Pivot Frame - Manual
LSS390670

Diverse Compatibility
Designed to be fitted on 200, A205 and A320.

Immediate Usage
Battery-free design allows immediate usage at all times with minimal maintenance required.



Pivot Frame - Powered
LSS390665

Diverse Compatibility
Designed to be fitted on 200, A205 and A320.

Minimises Carer Effort
Assists staff to conduct repetitive transfers in institutional care environments.

Integrated Electronics
Seamlessly integrates with lifter control system for use with one handset.

