MATISSE QUATTRO



FEATURES & BENEFITS

- Independent Adjustable Headrest, Back, Leg and Lumbar support functions
- Gentle Curve in Seating Surface for an improved seating experience
- Large Side Storage Pockets, ideal for books and newspapers
- Integrated Footrest Extension suits most user heights
- Integrated Lift and Recline functions for easy sit-to-stand
- Backlit Hand Control with integrated USB charging point

- Adjustable Headrest, ideal for watching TV whilst reclined
- **Lithium Ion Battery Backup** ensures chair can be used in the event of a power outage or if moved between rooms during the day
- 4 x Swivel Castors allow the chair to be moved around with ease







MATISSE quattro	SEAT HEIGHT (mm)	SEAT DEPTH (mm)	SEAT WIDTH (mm)	BACKREST HEIGHT (mm)	OVERALL HEIGHT (mm)	OVERALL WIDTH (mm)
	Ī				Ī	
CHP227560	<u> </u>				<u> </u>	ı i
	500	545	570	730	1065	965

^{*} Components: Handset, upholstery and electronics.