



MONTANA X MAXI SINGLE ACTION 1 MOTOR LIFT RECLINE CHAIR

Stylish Bariatric Seating for Your Home

Created for comfort and bariatric support, the Montana X Maxi features a zip-access adjustable backrest for personalised postural adjustments. Its sleek design fits effortlessly into modern homes, blending clinical function with elegant style.



Effortless Control
Simple 2-button handset



**Vertical
Lift**



**Full
Recline**



**Golden Grain
Microsuede**



- **Trusted Bariatric Design:** 220kg high safe working load and generous seating surface
- **Refined, Contemporary Design:** Larger, plush armrests provide extra comfort while maintaining a streamlined look that suits home environment
- **Integrated Footrest Extension:** Reduces fatigue and lower limb swelling by improving healthy circulation
- **Protective Covers for High-wear Areas:** Removable headrest and armrest covers keep your chair clean and long-lasting
- **Commercial-grade Seat Cushioning:** Enhances comfort and durability, providing lasting support
- **Promotes Independent Daily Living:** Simple 2-button handset assists users with cognitive or visual challenges
- **All-day Battery Backup:** Reliable Lithium-ion battery ensures continuous operation during power outages, maintaining your safety and independence



2 YEAR LIMITED WARRANTY*



Code	Colour	Size	Seat Height	Seat Depth	Seat Width	Backrest Height	Overall Height	Overall Depth	Overall Width
CHP227401	Golden Grain	Maxi	510mm	530mm	640mm	670mm	1060mm	940mm	1010mm

*Steel Frame; Motor; Hand Control; Scissor Mechanisms; Steel Lift Mechanisms; Fabric and Upholstery Care