



AIRLINE TRAVEL AND BATTERY INFORMATION

Travelling with Mobility Scooters & Power Wheelchairs

Airlines have specific safety regulations for mobility scooters and power wheelchairs due to the batteries used to power them. These regulations are based on international aviation safety standards and may vary between airlines. Before travelling, it is essential to contact your airline well in advance to confirm their requirements.



APOLLO SCOOTER

12V 12AH (LEAD ACID) **STANDARD** ✓



DELTA SCOOTER

24V 20AH (LITHIUM-ION) **STANDARD** ✗
24V 10AH (LITHIUM-ION) **OPTIONAL** ✓



ELIO SCOOTER

24V 20AH (LITHIUM-ION) **STANDARD** ✗
24V 10AH (LITHIUM-ION) **OPTIONAL** ✓



SUPALITE SCOOTER

24V 12AH (LITHIUM-ION) **STANDARD** ✓



AEROLITE POWERCHAIR

24V 12AH (LITHIUM-ION) **SPARE** ✓