



aspire

SIGNATURE CHAIRS

Quality Adjustable
Lift Recline Chairs
designed with
a gentle glide
movement and
luxurious Warwick
Fabric for comfort,
support and style.



Space Saving Design - Single Motor

Infinite Position - Dual Motor

Available in 3 Sizes

Customisable Seat Depth

High Performance Easy Clean Fabrics

Backlit Handset with USB Charger



aspire | SIGNATURE LIFT RECLINE CHAIR

Aspire Signature Lift Chairs are tailored with premium upholstery and high grade durable foams to provide years of postural support, comfort and safety.

The Signature Space Saver can be positioned closer to walls due to its subtle forward movement, thereby making the best use of your living space.

For stand-up support and reclining comfort the Signature Space Saver will surpass your lifestyle and home convenience expectations.

SPACE SAVER



Space saving recline for close wall positioning

Space Saving Positions - Single Motor



UPRIGHT POSITION LEG RECLINE FULL RECLINE FORWARD LIFT

Designed to enable leg elevation prior to back recline



SINGLE MOTOR

Back and legs move simultaneously

- ✓ DVA Contracted
- ✓ Signature is registered in the ARTG, AUST R ARTG 230186
- ✓ Tested to International Standards



DEPTH ADJUSTABLE SEAT
A choice of two depth settings in 1 inch (25mm) increments for custom comfort and positioning.



ACCESSIBLE STORAGE
Large side pocket for convenient storage. The chair's handset is conveniently fitted in the pocket, allowing users to reel in a dropped handset, discouraging leaning and possible falls.



MOBILITY READY
A set of 4 swivel castors facilitate easy movement for cleaning and transportation. Provided separately as standard (NB. Castors increase seat height).

2 YEAR WARRANTY



145kg

SAFE WORKING LOAD

PREMIUM QUALITY ELECTRONICS

European designed electronics provide a quiet and smooth lift, complete with an integrated battery backup. The large backlit handset buttons with USB charger point is suitable for mobile phones and other portable electronic devices.



CUSTOMISED COMFORT

Adjustable padding to support posture for individual users.



DURABLE DESIGN

High quality foam upholstery and steel frame for increased durability and postural support. Will not compress like polyfil fibres.



EASY CLEAN FABRIC

Warwick Global Suede Fabric offers commercial grade durability and resistance to flame, mould, mildew and pilling.



Aspire Signature Lift Chairs are tailored with premium upholstery and high grade durable foams to provide years of postural support, comfort and safety.

The Signature Dual Action features dual control of the back and legrest functions to enable infinite adjustments to suit your varied postural and support needs.

Along with customisable seat depth and padding the Signature Dual Action is an adaptable furniture choice to surpass your lifestyle and home convenience expectations.

DUAL ACTION



Dual control of leg and recline for infinite positioning



- ✓ DVA Contracted
- ✓ Signature is registered in the ARTG, AUST R ARTG 230186
- ✓ Tested to International Standards



DEPTH ADJUSTABLE SEAT
A choice of two depth settings in 1 inch (25mm) increments for custom comfort and positioning.



ACCESSIBLE STORAGE
Large side pocket for convenient storage. The chair's handset is conveniently fitted in the pocket, allowing users to reel in a dropped handset, discouraging leaning and possible falls.



MOBILITY READY
A set of 4 swivel castors facilitate easy movement for cleaning and transportation. Provided separately as standard (NB. Castors increase seat height).

2 YEAR WARRANTY

155kg

SAFE WORKING LOAD

Infinite Positions - Dual Motor



DUAL MOTOR
Back and legs move/recline independently

PREMIUM QUALITY ELECTRONICS

European designed electronics provide a quiet and smooth lift, complete with an integrated battery backup. The large backlit handset buttons with USB charger point is suitable for mobile phones and other portable electronic devices.



CUSTOMISED COMFORT

Adjustable padding to support posture for individual users.



DURABLE DESIGN

High quality foam upholstery and steel frame for increased durability and postural support. Will not compress like polyfil fibres.




EASY CLEAN FABRIC

Warwick Global Suede Fabric offers commercial grade durability and resistance to flame, mould, mildew and pilling.



TECHNICAL SPECIFICATIONS



Space Saver

Model	SEAT HEIGHT	SEAT DEPTH	SEAT WIDTH	SEAT TO TOP OF BACK	OVERALL HEIGHT	OVERALL DEPTH	OVERALL WIDTH	DISTANCE FROM WALL	USER HEIGHT GUIDE	S.W.L
CHP197700 Small - Mink	470mm	450 - 475mm	450mm	590mm	1010mm	690mm	780mm	140mm	5'2" and under	145kg
CHP197705 Small - Lagoon										
CHP197710 Medium - Mink	480mm	490 - 515mm	480mm	620mm	1050mm	715mm	780mm	150mm	5'3" to 5'8"	145kg
CHP197715 Medium - Lagoon										
CHP197716 Medium - Heritage										
CHP197720 Large - Mink	510mm	510 - 535mm	500mm	670mm	1110mm	740mm	780mm	160mm	5'9" to 6'2"	145kg
CHP197725 Large - Lagoon										
CHP197730 Large - Heritage										

Dual Action

Model	SEAT HEIGHT	SEAT DEPTH	SEAT WIDTH	SEAT TO TOP OF BACK	OVERALL HEIGHT	OVERALL DEPTH	OVERALL WIDTH	DISTANCE FROM WALL	USER HEIGHT GUIDE	S.W.L
CHP197735 Small - Mink	470mm	450 - 475mm	450mm	590mm	1010mm	710mm	780mm	430mm	5'2" and under	155kg
CHP197740 Small - Lagoon										
CHP197745 Medium - Mink	480mm	490 - 515mm	480mm	620mm	1050mm	735mm	780mm	455mm	5'3" to 5'8"	155kg
CHP197750 Medium - Lagoon										
CHP197751 Medium - Heritage										
CHP197755 Large - Mink	510mm	510 - 535mm	500mm	670mm	1110mm	760mm	780mm	480mm	5'9" to 6'2"	155kg
CHP197760 Large - Lagoon										

- * User Height Guide is a reference tool and may not apply to all users
- * Chair measurements may vary +/- 20mm due to upholstery tolerances
- * Seat heights shown are without castors on chair frames
- * Seat width is taken between the thickest sections of armrest padding


AVAILABLE IN 3 SIZES:



FABRIC RANGE



Warwick Global Suede - Mink



Warwick Global Suede - Lagoon



Vinyl - Heritage

GENUINE WARWICK GLOBAL SUEDE FABRICS

- High performance, easy clean fabrics tested to "Heavy Commercial" grade performance standards
- Proven colour fastness and abrasion and stain resistance
- Resistant to piling, mould and mildew; and impervious to moths and insects
- Fire retardant

VINYL HERITAGE FABRIC

- Easy clean
- Water resistant / stain resistant
- Durable and long lasting

CLINICAL CONSIDERATIONS

FALLS PREVENTION

- Raised posture maximises entry or exit access with mobility aids
- Handset cable runs through side wall pocket to remain in reach even if dropped

CHANGING NEEDS MANAGEMENT

- Adjustable seat depth and backrest upholstery with easy to clean fabric
- Lift feature facilitates reduced assistance as mobility improves in line with clinical rehab goals

ANTHROPOMETRIC CONSIDERATIONS

- Three sizes with inbuilt adjustability to ensure correct fit and postural support

OEDEMA MANAGEMENT

- Raised leg positioning to assist with lower limb oedema

VISION & HAND FUNCTION

- Backlit handset with large buttons

ASSISTED SITTING

- Independent backrest function allows gravity assisted upright sitting without raising legs

DAY BED POSITIONING

- Lie flat posture for improved patient comfort whilst maintaining mobility and independence

ACCESSORIES



Armrest Cover Protector



Head Cover Protector

CHP197770	Head/Armrest Cover Set - Aspire Signature Chairs - Mink
CHP197771	Head/Armrest Cover Set - Aspire Signature Chairs - Lagoon
CHP197772	Head/Armrest Cover Set - Aspire Signature Chairs - Heritage

ASSOCIATED PRODUCTS



LIFTERS & SLINGS



TRANSFER AIDS



DAILY LIVING AIDS



THERAPY & EXERCISE



CRUTCHES,
FRAMES & STICKS



SCOOTERS



OVERBED TABLES



BED STICKS & POLES



FOAM MATTRESSES



AIR MATTRESSES



COMMODOES
URINALS & BED PANS



OVER TOILET AIDS



SHOWER STOOLS & CHAIRS



RAISED TOILET SEATS




SHOWER COMMODOES



BATH BOARDS & BENCHES



SEAT WALKERS



WHEELCHAIRS



STD POWER WHEELCHAIRS



PRESSURE CUSHIONS




DAY CHAIRS



HEMOCARE BEDS

Distributed by:



ASBRT11.27 300119

