

TOILET SEAT RAISERS GUIDE

TOILET SEAT RAISERS (NO ARMS)

DESCRIPTION	SIZE	CODE
Toilet Seat Raiser	50mm	BTT149800
	100mm	BTT150200
	150mm	BTT150300
Toilet Seat Raisers with Lid	50mm	BTT149900
	100mm	BTT150100
	150mm	BTT150400

TOILET SEAT RAISERS WITH ARMS


DESCRIPTION	CODE
Freedom Toilet Seat Raiser	BTT164100
Raised Toilet Seat With Swing Back Arms	BTT150450
Raised Toilet Seat With Swing-Back Arms Wide	BTT150460
Etac Cloo Fixed Height Toilet Seat Raiser With Arms	BTT162030
Height Adjustable Toilet Seat Raiser	BTT157950

FITTING GUIDE FOR TOILET SEAT RAISERS

DESCRIPTION	DIMENSION	CODE	COMPATIBLE BOWL WIDTH (A)	COMPATIBLE BOWL LENGTH (B)	ADDITIONAL SEAT HEIGHT FROM BOWL (C)
Aspire Toilet Seat Raiser	50mm	BTT149800	210-370mm	265-400mm	50mm
	100mm	BTT150200			100mm
	150mm	BTT150300			150mm
Aspire Toilet Seat Raisers with Lid	50mm	BTT149900			50mm
	100mm	BTT150100			100mm
	150mm	BTT150400			150mm



DESCRIPTION	CODE	COMPATIBLE BOWL WIDTH (A)	COMPATIBLE BOWL LENGTH (B)	MIN DISTANCE TO BOWL EDGE FROM MOUNTING SCREW (C)	ADDITIONAL SEAT HEIGHT FROM BOWL (D)	DISTANCE BETWEEN ARMRESTS
Freedom Toilet Seat Raiser	BTT164100	240-380mm	250-390mm	110mm	50-130mm	470mm
Raised Toilet Seat with Swing Back Arms	BTT150450	210-340mm	270-360mm	95mm	60-150mm	370mm
Etac Cloo Fixed Height Toilet Seat Raiser With Arms	BTT162030	Min: 205mm	250-400mm	100mm	110mm	400mm
Raised Toilet Seat with Swing-Back Arms – Wide	BTT150460	255-380mm	240-420mm	115mm	80-155mm	435mm
Etac Cloo Fixed Height Toilet Seat Raiser with Arms	BTT162030	205-380mm	250-350mm	70mm	60-140mm	490mm
Height Adjustable Toilet Seat Raiser	BTT157950	260-390mm	260-400mm	115mm	50-130mm	430mm



OTHER RELEVANT INFORMATION

DESCRIPTION	CODE	STYLE	LID	ARMS	SEAT HINGES CONCEALED?
Aspire Toilet Seat Raiser – 50mm	BTT149800	Pressure locked	No	None	No hinges
Aspire Toilet Seat Raiser – 100mm	BTT150200	Pressure locked	No	None	No hinges
Aspire Toilet Seat Raiser – 150mm	BTT150300	Pressure locked	No	None	No hinges
Aspire Toilet Seat Raisers with Lid – 50mm	BTT149900	Pressure locked	Yes	None	No
Aspire Toilet Seat Raisers with Lid – 100mm	BTT150100	Pressure locked	Yes	None	No
Aspire Toilet Seat Raisers with Lid – 150mm	BTT150400	Pressure locked	Yes	None	No
Freedom Toilet Seat Raiser	BTT164100	Bowl mounted*	Yes	Swing back	
Raised Toilet Seat with Swing Back Arms	BTT150450	Bowl mounted*	Yes	Swing back	Yes
Raised Toilet Seat with Swing-Back Arms – Wide	BTT150460	Bowl mounted*	Lid and raisable seat	Swing back	No
Etac Cloo Fixed Height Toilet Seat Raiser with Arms	BTT162030	Bowl mounted*	Yes	Swing back	Yes
Height Adjustable Toilet Seat Raiser	BTT157950	Bowl mounted*	Lid and raisable seat	Swing back	No

*Bowl mounted toilet seat raisers require access to the existing seat hinge mounting point on the toilet.