VOGUE CARBON FIBRE WALKER

Aspire for... Comfort, Stability & Style



FEATURES & BENEFITS

- Luxury look and feel Danish design, profiled tubing, moulded fittings, ergonomic handles and refined colours
- Super lightweight and ergonomic design with small turning circle for excellent manoeuvrability
- **AutoLocate** Smooth height adjustment that automatically locates to pre-set heights, making adjustment simple and accurate
- Locks when folded 225mm folding width minimises storage space and makes transportation easier
- Advanced locking handbrakes for increased user safety when walking and sitting
- Unbreakable TPE super soft wheels offers high durability and a smooth ride for increased comfort
- Removable utility bag simple and stylish to remove for shopping
- Supplied with cane holder attachment
- Deep Seat and Wide Backrest for improved comfort and support

	OVERALL HEIGHT	OVERALL LENGTH	OVERALL WIDTH	FOLDED WIDTH	SEAT DEPTH	SEAT WIDTH	SEAT HEIGHT	WHEEL SIZE	TOTAL WEIGHT	BAG S.W.L	S.W.L
						Ä	5	d	kg	kg	kg
WAF705350 Aspire Vogue Carbon Fibre - Black - Medium	695 - 920mm	680mm	600mm	225mm	240mm	460mm	550mm	8"	5kg	5kg	150kg
WAF705360 Aspire Vogue Carbon Fibre - Candy Red - Tall	770 - 1040mm	680mm	605mm	225mm	240mm	460mm	620mm	8"	6kg	5kg	150kg