


WALKER RANGE

Premium Mobility Solutions



VOGUE CARBON FIBRE SEAT WALKER

One of the lightest seat walkers on the market! Our Vogue Carbon Fibre Seat Walker epitomises function and style with a super lightweight carbon fibre frame for a superior walking experience.

CARBON FIBRE FRAME: Extremely lightweight to maximise portability

COLOURS:



MODELS

WAF705350 WAF705360 Vogue Carbon Fibre - Medium - Black Vogue Carbon Fibre - Tall - Candy Red



VOGUE LIGHTWEIGHT SEAT WALKER

This modern seat walker has a clean and crisp frame design, making it a stand out when going for a walk or browsing the shops!

LIGHTWEIGHT ALUMINIUM BODY: Improves manoeuvrability and portability.

COLOURS:



Red



Champagne



Leopard





Discreet Cable Management



MODELS

WAF705310RE WAF705310CH

WAF705310LE

Vogue Lightweight Aluminium - Red
Vogue Lightweight Aluminium - Champagne

Vogue Lightweight 2 - Leopard

Luxury Look & Feel Danish Designed

FEATURES & BENEFITS SHARED BY THE VOGUE CARBON FIBRE ...



Ergonomic Design
Offers personalised height
adjustment which promotes
correct posture



Compact
Turning circle for improved mobility in tight spaces



Simple
Folding mechanism locks
for easy transportation
and storage



Handbrakes Lock for user safety and stability

... AND LIGHTWEIGHT ALUMINIUM SEAT WALKERS



Removable bag carries up to 5kg of essentials



Backrest
Deep seat and wide backrest for improved comfort and support



Bonus Walking cane holder







Typical



Minimum

	Overall Length	Overall Width	Overall Height	Seat Depth	Seat Width	Seat Height	Folded Width	Wheel Size	Minimum Functional Weight	Typical Functional Weight**	Bag SWL	
						Toto		ď	£	1	7	
Medium Tall	680mm	600mm	695–920mm	240mm	460mm	550mm	225mm	0.11	4.9kg	6.2kg	5kg	
	26.8"	23.7"	27.4"-36.2"	9.4"	18.1"	21.7"	8.9"	8"	10.8lbs	13.7lbs	11lbs	
	680mm	605mm	770–1040mm	240mm	460mm	620mm	225mm	Oll	5.2kg	6.3kg	5kg	
	26.8"	23.8"	30.3"-40.9"	9.4"	18.1"	24.4"	8.9"	8"	11.5lbs	13.7lbs	11lbs	

Overall Length	Overall Width	Overall Height	Seat Depth	Seat Width	Seat Height	Folded Width	Wheel Size	Functional Weight	Functional Weight**	Bag SWL
							3	£	20	£
680mm	600mm	695–920mm	240mm	460mm	550mm	225mm	Oll	4.9kg	6.2kg	5kg
26.8"	23.7"	27.4"-36.2"	9.4"	18.1"	21.7"	8.9"	8"	10.8lbs	13.7lbs	11lbs

^{*} Carbon Fibre Walker - Medium | ** Includes accessories, bag and backstrap.



aspirecare.com.au

VOGUE ADVENTURE SEAT WALKER

Rugged and durable, the Adventure Seat Walker is the perfect all terrain walker. Super sized front wheels can manage any surface, perfect for the adventurous user!

COLOURS:



Black



Wider Seat supports up to 200kg

MODELS

Vogue Adventure - Black





Large Multi-Terrain Wheels



Folds to 245mm



VOGUE FOREARM WALKER

For those who need a little more support, the Forearm Walker is for you. Soft and contoured forearm pads ensure maximum comfort and a smooth ride.

LIGHTWEIGHT ALUMINIUM BODY: Improves manoeuvrability and portability.

COLOURS:



White / Black







WAF705400 Vogue Forearm - White/Black



Padded Forearm Supports



Folds to 225mm



FEATURES & BENEFITS



WAF705450

Ergonomic Design Offers personalised height adjustment which promotes correct posture



Handbrakes Lock for user safety and stability



Simple Folding mechanism locks for easy transportation and storage



Compact Turning circle for improved mobility in tight spaces



Removable bag carries up to 5kg of essentials



Bonus Walking cane holder







FEATURES & BENEFITS

Bonus

Walking cane holder



Ergonomic Design Offers personalised height adjustment which promotes correct posture



Handbrakes Lock for user safety and stability



Simple Folding mechanism locks for easy transportation and storage





Removable bag carries up to 5kg of essentials







Minimum Functional Weight	Typical Functional Weight**	Bag SWL	
_	_	_	

Overall Length	Overall Width	Overall Height	Seat Depth	Seat Width	Seat Height	Folded Width	Wheel Size	Functional Weight	Functional Weight**	Bag SWL
							3	£	£	£
770mm	695mm	900–1035mm	230mm	540mm	530mm	245mm	Front 11"	8.8kg	10kg	5kg
30.3"	27.4"	35.4"-40.7"	9.1"	21.3"	20.9"	9.6"	Rear 9"	19.4lbs	22lbs	11lbs

^{**} Includes accessories, bag and backstrap.



aspirecare.com.au

Overall Overall Overall Wheel Length Width Height Depth Width Height Width 660mm 620mm 1020-1270 240mm 460mm 620mm 225mm 24.4" 40.2"-50" 9.4" 18.1" 24.4" 8.9" 20.1lbs

VOGUE INDOOR WALKER

Ideal for use within the home, the Indoor Walker is perfect for portable dining and in home mobility assistance.

COLOURS:



Red / Black



Champagne / Black

MODELS

Vogue Indoor - Champagne/Black





Folds to 270mm

VOGUE WALKER ACCESSORIES



Oxygen Bottle Holder Bag WAA695600



Night Light - LED Lamp

WAA695630



Single Hand Brake System* – Compatible with Aspire Vogue Carbon Walkers

WAA695640

*Single Hand Brake System are not compatible with the Vogue Indoor Walker

WAF705100RE

WAF705100CH

Vogue Indoor - Red/Black

FEATURES & BENEFITS



Ergonomic Design

Offers personalised height adjustment which promotes correct posture



Handbrakes

Lock for user safety and stability



Simple

Folding mechanism locks for easy transportation and storage



Portable Dining

Safely collect and transport meals or hot drinks



Compact

Turning circle for improved mobility in tight spaces



Storage

Removable bag carries up to 5kg of essentials




Bonus

Walking cane holder



Turn Around. Slide Over Knee. Lock & Enjoy.







Overall Length	Overall Width	Overall Height	Folded Width	Wheel Size	Minimum Functional Weight	Typical Functional Weight**	Bag SWL	Tray SWL
				3				
650mm	540mm	860–970mm	270mm	CII	4.4kg	5.1kg	5kg	5kg
25.6"	21.3"	33.9"-38.2"	10.7"	6"	9.7lbs	11.2lbs	11lbs	11lbs

^{**} Includes accessories, bag and tray.

IV Pole Lamp

WAA695651

ASSOCIATED PRODUCTS





MOBILITY SCOOTERS

Wide range of Aspire Scooters for all user needs.

Shown Aspire HS268 Manual Folding Mobility Scooter SCT671135



COMFIMOTION

COLLECTION

Range of adjustable beds and accessories for both clinical needs and comfort at home.

Shown Aspire ComfiMotion Activ Care Bed in Portebello Grey - Long Single - BEB046740



LIFT RECLINE CHAIRS

Luxury Aspire recliner range has four, dual and single motor functions.

Shown Aspire Raphael Quattro Powerlift Recline Chair Sizes Small – Maxi Available



OVERBED

TABLES

Lightweight mobile table positions easily over the lifestyle community beds.

Shown Over Bed Table with Gunmetal Frame & Sublime Teak Laminate Top BEA006705