

	Scooter Code	Image
SMALL FOLDING	HS268 SCT671135 - White	
	SUPALITE SCT670616 - White	
SMALL PULL-APART	SONIC SCT670665 - Red	
	GT SCT670670 - Red	
MIDI	HS520 SCT671210 - White	
	HS559 SCT671240 - Red	
	HS589 SCT671260 - Blue SCT671270 - Red	
	HS828 SCT671285 - Red	
LARGE	PURSUIT SCT670571 - Red	
	HS898 BRAVO SCT671290 - Black SCT671300 - Blue	
	HS928 SCT671310 - Silver	
MAXI	RANGERIDER SCT670581 - Stone Grey	

Terrain	Battery	Max Range Per Charge*	Safe Working Load
Level & smooth terrain. Short distances.	1 x 12V. 11.6AH (Lithium Ion)	20 km	115 kg
Level & smooth terrain. Short distances.	1 x 24V. 12AH (Lithium Ion)	18.5 km	120 kg
Level & smooth terrain. Short distances.	2 x 12V. 20AH	20 km	136 kg
Able for more varied terrain. Short to medium distances.	2 x 12V. 35AH	30 km	159 kg
Able for more varied terrain. Short to medium distances.	2 x 12V. 50AH	39 km	136 kg
Able for more varied terrain. Short to medium distances.	2 x 12V. 50AH	42.5 km	136 kg
Able for more varied terrain. Short to medium distances.	2 x 12V. 50AH	40 km	135 kg
Can handle rougher terrain well. Long distances.	2 x 12V. 75AH	27 km	160 kg
Can handle rougher terrain well. Long distances.	2 x 12V. 55AH	25 km	159 kg
Can handle rougher terrain well. Long distances.	2 x 12V. 75AH	40 km	181 kg
Can handle rougher terrain well. Long distances.	2 x 12V. 100AH	50 km	226 kg
Can handle rougher terrain well. Long distances.	2 x 24V. 100AH	30 km	159 kg


SCOOTER ACCESSORIES



Oxygen Tank Holder




Cane Holder



Crutch Holder



Rear Shopping Bag



Scooter Flag



Sun Shade Canopy



Cupholder



Dust Cover



Rear Bucket

*Not all accessories are compatible with all models - see website for details

SAFE USE VIDEO

We value the safety of our customers above all else. There are a few safety considerations we'd like you to understand before operating your device.

We strongly recommend all new mobility device customers watch our 7 minute YouTube video, 'Getting to Know your Power Mobility Device' that covers the above 8 points in detail and more.



aspire for..... mobility comfort style

MOBILITY SCOOTER SELECTION GUIDE



HS268

SCT671135 - White



FEATURES

- Manual folding via hand lever
- LED headlights and brakelights
- Intuitive user control interface
- Approved for air travel
- Electro Mechanical Brakes
- Fast charging Lithium Ion battery

SPECIFICATIONS

Length: 950 mm (non-folded) / 780 mm (folded)
Width: 500 mm (non-folded) / 775 mm (folded)
Product Weight: 30.6 kg (incl batteries)
Maximum Speed: 8 km/h
Range: 20 km
Battery Size: 1 x 12V. 11.6AH (Lithium Ion)

115kg
S.W.L

SUPALITE

SCT670616 - White



FEATURES

- Disassembles into 6 lightweight pieces
- Compact, foldable seat with armrests
- Easy-to-grip tiller adjustment knob
- Stable front frame-mounted seat post
- Standard front basket provides storage

SPECIFICATIONS

Length: 1022 mm
Width: 480 mm
Product Weight: 34 kg (incl batteries)
Maximum Speed: 6.5 km/h
Range: 18.5 km
Battery Size: 1 x 24V. 12AH (Lithium Ion)

120kg
S.W.L

SONIC

SCT670665 - Red



FEATURES

- Disassembles into 5 light weight pieces
- Portable and durable scooter
- Puncture-proof tyres
- LED headlights and brakelights
- Detachable battery pack
- Easy access battery compartment,
- Adjustable angle tiller with delta

SPECIFICATIONS

Length: 1040 mm
Width: 495 mm
Product Weight: 52 kg (incl batteries)
Maximum Speed: 6.8 km/h
Range: 20 km
Battery Size: 2 x 12V. 20AH

136kg
S.W.L

GT

SCT670670 - Red



FEATURES

- Advanced suspension for a smooth ride
- Front and rear LED lighting
- Directional signals
- Deluxe tiller with USB charger
- Cup holder

SPECIFICATIONS

Length: 1200 mm
Width: 605 mm
Product Weight: 103.3 kg (incl batteries)
Maximum Speed: 10 km/h
Range: 30 km
Battery Size: 2 x 12V. 35AH

159kg
S.W.L

HS520

SCT671210 - White



FEATURES

- Sleek, Sporty design
- Stylish mag wheels
- LED headlights and brakelights
- Front and rear wheel suspension
- Rearview mirror
- Adjustable swing away armrests
- Electro mechanical brakes

SPECIFICATIONS

Length: 1200 mm
Width: 640 mm
Product Weight: 94 kg (incl batteries)
Maximum Speed: 10 km/h
Range: 39 km
Battery Size: 2 x 12V. 50AH

136kg
S.W.L

HS559

SCT671240 - Red



FEATURES

- Robust design with a large range
- LED headlights and brakelights
- Front and rear wheel suspension
- Rearview mirror
- Adjustable swing away armrests
- Electro mechanical brakes

SPECIFICATIONS

Length: 1200 mm
Width: 625 mm
Product Weight: 101 kg (incl batteries)
Maximum Speed: 10 km/h
Range: 42.5 km
Battery Size: 2 x 12V. 50AH

136kg
S.W.L

HS589

SCT671260 - Blue
SCT671270 - Red



FEATURES

- Classic design with a large range
- LED headlights and brakelights
- Front and rear wheel suspension
- Rearview mirror
- Adjustable swing away armrests
- Electro mechanical brakes

SPECIFICATIONS

Length: 1290 mm
Width: 625 mm
Product Weight: 98 kg (incl batteries)
Maximum Speed: 10 km/h
Range: 40 km
Battery Size: 2 x 12V. 50AH

135kg
S.W.L

HS828

SCT671285 - Red



FEATURES

- Modern design with digital user interface
- Stylish mag wheels
- Headlights, brakelights and indicators
- Front and rear wheel suspension
- LED daytime running headlights
- Rearview mirrors
- Adjustable swing away armrests
- Electro mechanical brakes

SPECIFICATIONS

Length: 1400 mm
Width: 680 mm
Product Weight: 124 kg (incl batteries)
Maximum Speed: 10 km/h
Range: 27 km
Battery Size: 2 x 12V. 75AH

160kg
S.W.L

PURSUIT

SCT670571 - Red



FEATURES

- High performance, sleek and sporty
- Large 13" pneumatic tyres
- LED headlights, brakelights and indicators
- Armrests head rest and lumbar support
- Rugged front and rear suspension

SPECIFICATIONS

Length: 1360 mm
Width: 635 mm
Product Weight: 91 kg (incl batteries)
Maximum Speed: 10 km/h
Range: 25 km
Battery Size: 2 x 12V. 55AH

159kg
S.W.L

HS898 BRAVO

SCT671290 - Black
SCT671300 - Blue



FEATURES

- Heavy Duty Hill Climber
- Stylish mag wheels
- LED headlights, brakelights and indicators
- Front and rear wheel suspension
- Rearview mirrors
- Adjustable swing away armrests
- Electro mechanical brakes

SPECIFICATIONS

Length: 1440 mm
Width: 720 mm
Product Weight: 143 kg (incl batteries)
Maximum Speed: 10 km/h
Range: 40 km
Battery Size: 2 x 12V. 75AH

181kg
S.W.L

HS928

SCT671310 - Silver



FEATURES

- MAXI scooter to cater for larger users
- Stylish mag wheels
- LED headlights and brakelights
- Front and rear wheel suspension
- Rearview mirrors
- Adjustable swing away armrests
- Electro mechanical brakes

SPECIFICATIONS

Length: 1600 mm
Width: 720 mm
Product Weight: 178 kg (incl batteries)
Maximum Speed: 10 km/h
Range: 50 km
Battery Size: 2 x 12V. 100AH

226kg
S.W.L

RANGERIDER

SCT670581 - Stone Grey



FEATURES

- Limited recline, black vinyl, high-back seat with lumbar and pillow-top headrest, flip-up height, and angle adjustable armrests
- CTS Suspension includes adjustable shocks for greater comfort
- Backup sensors with LED light
- Dual hydraulic brakes

SPECIFICATIONS

Length: 1560 mm
Width: 830 mm
Maximum Speed: 10 km/h
Range: 30 km
Battery Size: 2 x 24. 100AH

159kg
S.W.L