



30/08/21

## ASPIRE LIFESTYLE COMMUNITY BEDS – NOW AVAILABLE!

At long last, our brand new Aspire Lifestyle Community Bed range has arrived!

With warm timber finishes paired with a modern dark frame, the Lifestyle bed not only looks sophisticated in the home environment but with five profiling functions it also caters to a broad spectrum of comfort and positioning needs for the ultimate in home care provision.

The robust and easily assembled frame makes for the perfect rental fleet bed, built to withstand the rigors of regular transportation and assembly.

Now available in both Single & King Single sizes too!

### KEY DIFFERENCES & IMPROVEMENTS

- o **Quick setup time:** Tool free with no requirement to remove electronics during assembly and set-up, ensuring quick and easy readiness upon delivery or for transportation and storage. For someone working alone it should take a maximum of 8 minutes to assemble a Lifestyle bed!
- o **Additional sizes:** now available in both single and king single widths to suit most customer needs (Please Note - MAXI (385kg SWL) version is also coming soon!)
- o **Homecare upgrade options:** Timber laminate head and foot board shrouds and full-length timber rails can be fitted for an appealing homely look (pictured below).
- o **Supplied Transport-Ready:** The Lifestyle beds come pre-assembled on transport brackets, so you can simply roll it out of the carton on its wheels, ready to deliver and set-up!

#### ***Aspire Lifestyle Community Bed***

*Shown fitted with Homecare upgrade options:*

- Timber Laminate Bed End Shrouds
- Timber Laminate Full Length Side Rails



New Product Launch

## BED & ACCESSORY CODES

Product Code	Description
<b>Beds</b>	
BEB046030	Bed - Aspire LIFESTYLE Community Bed w/Transport Bracket - Gunmetal Grey - Single - 220kg
BEB046040	Bed - Aspire LIFESTYLE Community Bed w/Transport Bracket - Gunmetal Grey - King Single - 220kg
<b>Shrouds</b>	
BEA029850	Bed End Shrouds - Single - Aspire LIFESTYLE Community Bed
BEA029855	Bed End Shrouds - King Single - Aspire LIFESTYLE Community Bed
<b>Rails</b>	
BEA016586	Bed Side Rails - 3/4 Folding - Aspire LIFESTYLE Community Beds
BEA029863	Bed Side Rails - Full Length - Aspire LIFESTYLE Community Care Bed
BEA029864	Bed Side Rail Extension Kit - Aspire LIFESTYLE Community Bed
<b>Other Accessories</b>	
BEA016592	Bed IV Pole - Aspire LIFESTYLE & Community Care Beds
BEA016596	Bed Self Help Pole - Aspire LIFESTYLE Community Care Bed
BEA029871	Bed Transport Kit - Aspire LIFESTYLE Community Bed

