

07/10/21

## ASPIRE MAXI LIFESTYLE COMMUNITY BED - WITH ELECTRIC FOLDING!

This all new ELECTRICALLY FOLDING bariatric bed solution is now available in stock.

## This product solves several problems:

- Simple and economical solution for up to 350kg bariatric patient care, most of the market is offering non-foldable hospital ward style beds into this category.
- Safe moving and handling of heavy equipment for resellers and healthcare staff.
- Rapid installation/removal of bariatric care solutions which is needed by hospitals and Aged Care providers who provide respite services.
- Compact folded size provides efficient storage as a hire/ward asset.

## **KEY FEATURES**

- **Quick up time:** Single button press electrical folding mechanism requires no lifting and improves workplace health & safety outcomes for delivery staff and maintenance teams.
- **True bariatric size:** Super-King size mattress platform caters to true bariatric clients reducing falls risk and moving and handling complexity.
- **Auto-regression:** in the backrest reduces skin shear forces and removes a significant cause of patient slippage. This reduces the need for care staff to reposition oversized patients.
- **Finishes:** Stylish wooden end shrouds are included as standard and offer a simple solution for cleaning and infection control.
- Mattresses: Compatible with our new Maxi Lifestyle Foam Mattress (BEM046710) and Active Air Bari 10 Super King (PTM639520) to provide static and dynamic sleep surface options.









## BED, ACCESSORY AND MATTRESS CODES

Product Code	Description
Beds	
BEB046020	Bed - Aspire MAXI LIFESTYLE Community Bed - Electric Folding - Gunmetal Grey - 385kg
Rails	
BEA016586	Bed Side Rails - 3/4 Folding - Aspire LIFESTYLE Community Beds
Accessories	
BEA016592	Bed IV Pole - Aspire LIFESTYLE & Community Care Beds
BEA016598	Bed Self Help Pole - Aspire MAXI LIFESTYLE Community Care Bed
Mattresses	
BEM046710	Mattress - Aspire Lifecomfort - Maxi Lifestyle - Super King - 2060x1150x200mm - 350kg
PTM639520	Mattress Replacement - Aspire ACTIVE AIR BARI 10 - Super King - 20 Cells - 2000x1230x270mm - 450kg

