



21/04/2022

ADJUSTABLE STROLLER HANDLE - ASPIRE REHAB RX WHEELCHAIRS

We are pleased to let you know about a new accessory for our Aspire Rehab RX and RX JNR Wheelchairs, the Adjustable Stroller Handle.

The Adjustable Stroller Handle is mounted on the REHAB RX/RX JNR backrest, replacing the standard handle canes. The attendant brakes are removed from the standard handle canes and re-attached to the Adjustable Stroller Handle (mounting instructions are included in the product carton).

The Stroller Handle together with brakes, can be pivoted up and down to suit the carers positioning needs.

PRODUCT FEATURES

- Ergonomic hand grip for carers - especially helpful for tilt and recline movements
- Height and width adjustable to suit all Rehab RX/RX JNR chair sizes
- Padded push handle with brakes
- Pivots up and down via button press positioning
- **Product Code: WCA848590**



New Product Launch