

03/03/2022

ASPIRE COMFIMOTION – NEW ACTIV CARE BED + RANGE UPDATES!

Introducing the newest addition to the Aspire ComfiMotion range, the ACTIV CARE Bed.


The ACTIV CARE bed combines the absolute best of in-home care functionality, positioning and ultra low height safety with luxurious contemporary style for enhanced comfort, care and wellbeing.

ACTIV CARE – NEW AND DISTINGUISHING FEATURES

The following key features or functions distinguish the top of the range ACTIV CARE bed from the CARE bed:

- **Ultra Low Height:** 195mm low mattress platform height significantly reduces the risk of injury in the event of a fall from bed.
- **Massage:** Relaxing vibration modes for personalised neck, shoulder and back relief.
- **Customisable Safe Entry/Exit Height:** A carer or family member can pre-set the safe exit height ensuring it is at the safest level for the particular user.
- **Under Bed Clearance:** Provides clear access for patient lifters and overbed tables.
- **Head Tilt:** Relaxes neck muscles, providing optimal support whilst reading or watching TV.
- **Battery Backup:** Have confidence that the bed will still operate during a power outage (*included as standard, previously an accessory*).
- **Auto-Regression in the Backrest:** Stops the user sliding down the bed during movement while improving comfort and maximising pressure redistribution.
- **USB Charge:** Connect and charge devices conveniently.






OTHER COMFIMOTION RANGE UPDATES

The success of the ComfiMotion range to date has inspired the design and development of the ACTIV CARE bed, but it is not the only story:

- **King Size Headboards - Now Available** – 1 x King headboard for Split Kings instead of 2 x Singles.
- **Split King Upgrade Kits – Now Available** - Facilitates synchronised functioning of Split Kings.
- **Simplified Handsets – Now available** - For clients with reduced cognitive ability.
- **Bed Rails – Now Available** – Now specific to either the CARE or ACTIV CARE beds, please ensure the correct version is selected (see code table below).
- **LUXE Bed: Discontinued** – Massage, Head Tilt, USB, and more are now available on the ACTIV CARE.
- **CARE Bed: Continues** – Fewer functions/features than the ACTIV CARE but very strong value.
- **COMPANION Bed: Continues.**





PRODUCT CODES - AVAILABLE NOW

Product Code	Description	Status
Beds		
BEB046740	Bed - Aspire ComfiMotion - ACTIV CARE - Long Single	New!
BEB046750	Bed - Aspire ComfiMotion - ACTIV CARE - King Single	New!
BEB046710	Bed - Aspire ComfiMotion - LUXE - King Single	Discontinued
BEB046720	Bed - Aspire ComfiMotion - CARE - Long Single	Existing
BEB046730	Bed - Aspire ComfiMotion - CARE - King Single	Existing
Headboards		
BEA830100	Headboard - Aspire ComfiMotion - CHESTERFIELD - Long Single	Existing
BEA830110	Headboard - Aspire ComfiMotion - CHESTERFIELD - King Single	Existing
BEA830116	Headboard - Aspire ComfiMotion - CHESTERFIELD - King	New!
BEA830120	Headboard - Aspire ComfiMotion - ALLURE - Long Single	Existing
BEA830130	Headboard - Aspire ComfiMotion - ALLURE - King Single	Existing
BEA830136	Headboard - Aspire ComfiMotion - ALLURE - King	New!
Rails		
BEA830200	Bed Rail - Crook Handle - Aspire ComfiMotion - LUXE & CARE	Existing
BEA830210	Bed Rail - Aspire ComfiMotion - LUXE & CARE	Existing
BEA830211	Bed Rail - Crook Handle - Aspire ComfiMotion - ACTIV CARE	New!
BEA830212	Bed Rail - Aspire ComfiMotion - ACTIV CARE	New!
Accessories		
BEA830214	Handset - Simplified - Aspire ComfiMotion - ACTIV CARE	New!
BEA830230	Split King Upgrade Kit - Aspire ComfiMotion - ACTIV CARE	New!
BEA830220	Battery Backup & Power Transformer- Aspire ComfiMotion - LUXE & CARE	Existing
MWT451080	Large Accessibility Castor - Aspire ComfiMotion - LUXE & CARE	Existing
Mattresses		
BEM469800	Mattress - Aspire ComfiMotion - POCKET SPRING - Long Single	Existing
BEM469810	Mattress - Aspire ComfiMotion - POCKET SPRING - King Single	Existing
BEM469910	Mattress - Aspire ComfiMotion - MEMORY FOAM - Long Single	Existing
BEM469920	Mattress - Aspire ComfiMotion - MEMORY FOAM - King Single	Existing
BEM469930	Mattress - Aspire ComfiMotion - MEMORY FOAM w/Incontinence Cover - Long Single	Existing
BEM469940	Mattress - Aspire ComfiMotion - MEMORY FOAM w/Incontinence Cover - King Single	Existing
Pillows		
BEP050100	Pillow - Aspire ComfiMotion - BREEZE	Existing
BEP050110	Pillow - Aspire ComfiMotion - PLUSH	Existing

