

NEW PRODUCT LAUNCH



Aspire Genesis HD Seat Walker – Available Soon!

The new Aspire Genesis HD is height-adjustable with a seat height range from 500mm to 590mm as well as adjustable handle height.

As such it replaces all 3 sizes of the previous model Freedom Genesis HD which have not been available for some time.



Aspire Genesis HD Walker - WAF750336

Adjustable Seat Height: Suitable for variable user sizes, optimising seat height for the user assists with sitting and standing easily and safely.

Adjustable Handle Height & Soft-Touch Brakes: Correct handle height promotes walking posture while soft brakes are nice and easy to operate.

Wide Padded Seat & High Weight Capacity: Ensures Comfort and stability for larger users, with the added benefit of a discreet storage pocket beneath the seat.

Large Padded Backrest: Provides generous and comfortable lumbar support when seated.

Oversized Under-Seat Bag: Ample storage space for personal items, accessible by lifting the seat.

Easy Folding Mechanism and Lightweight: Lift to fold for compact storage and easy handling at only 11.3kg total weight.

Solid 8" Wheels with Inline Rear Brakes: Tough treaded wheels for a smooth ride in variable terrain, with inline rear brakes maximise manoeuvring space for the user.

Vibrant Blue Marble Finish