

Aspire Montana HD Lift Recline Chair – Available!


The Montana HD Lift Recline Chair has recently undergone an important revision. The weight capacity of the Montana HD350 (**CHP227405**) model has changed to 300kg, while the Montana HD400 (**CHP227410**) model now supports 450kg.



Commercial grade foams and fabric offer durability and support for prolonged sitting and reclining.



HUGE 160KG LEGREST WEIGHT CAPACITY!



CHP227405- Montana HD - Barley Vinyl - 300kg CHP227410- Montana HD - Barley Vinyl - 450kg

For all custom codes, refer to the custom order form.



PRODUCT UPDATES

PRODUCT FEATURES

- True Bariatric Seating: 300kg & 450kg options available.
- **Perfect Positioning -** The Montana HD allows independent movement of the legrest and backrest to achieve optimal posture for comfort and oedema management
- **Configurable Backrest Options:** Comfort, Lateral Supported and WaterFall designs with rear zippers for posture and gluteal shelf management.
- Durable Construction: Commercial grade foams and fabric for long duration sitting and reclining.
- **Custom Upholstery:** Stocked in Barley vinyl but available in over 70+ colours (made to order) in Leather, Vinyl, Microsuede and Soft Touch fabric. All these colours are in the same range we have on our popular custom Raphael chair offering. These are available via the custom order form.
- **Custom Sizing:** Capability to offer clients a custom-sized lift chair with the ability to specify the seat width, seat depth, seat height and armrest height. This option is available via the custom order form.

CUSTOM ORDER & QUOTE FORM

We have created a custom order form (attached to this release) which works in a very similar way to the custom Raphael Chair order form. Pick and choose your colour and chair size for made-to-order options or choose a standard size which will result in a quicker lead time (modest stock available – *CHP227330*)

The product team ran through this chair and how to arrange orders in the May 2023 product presentation. Please check the Aidacare intranet or website out for more details!





Issued: 22nd September 2023