

Aspire PostureFit Lift Recline Chair – Available Now!

This innovative seating solution is a game-changer for those with unique and changing needs. The PostureFit chair offers a uniquely adjustable and adaptable seating solution. Through its seat width and depth adjustment as well as broad range of backrest and seat cushion options, the PostureFit chair provides unparalleled comfort and support.



CHP198250 – PostureFit – Small/Medium CHP198260 – PostureFit – Tall Medium/Large **For scripting and customisation, refer to the assessment form.**



Issued: 11th December 2023

NEW PRODUCT LAUNCH

PRODUCT FEATURES

- **Tool-free Seating Adjustment:** Therapists or carers can quickly and easily make adjustments without the need for any tools.
 - **Seat Depth:** Backrest can be quickly and easily adjusted to achieve two seat depths with a 50mm range.
 - Seat Width: Armrests are interchangeable for left/right side, swap sides to achieve two seat widths with a 50mm range.
- **Zero Gravity Positioning:** Creates a feeling of weightlessness and redistributes pressure off the lower back for pain management.
- **Oedema Management:** Leg elevation above the heart promotes drainage of lower limb swelling, aids circulation, and reduces fatigue.
- Adjustable Headrest: Electrically adjustable via the handset, the headrest cushion can be tilted forward and backward for optimal support.
- Integrated Leg Extension: Suits most user heights without the need to alter the legrest or height of the chair. The leg extension can also be easily removed for users with specific heel care needs
- Vertical Lift Mechanism: Offers sitting height adjustment for taller users and a stable sit to stand motion that removes the sensation of sliding forward for improved confidence. Available on all models without compromising SWL, unlike competing products.
- Effortless mobility and Safety: Large directional castors are engaged by activating a rear carer foot pedal. Castors are not in contact with the floor when in a neutral seated position, ensuring a strong base of support. This also removed any risk of sideways movement as opposed to other products with front swivelling castors
- **Transport mode:** For safety, all electrical functions are inactivated while castors are engaged to prevent the chair from being activated while mobile.

INTERCHANGEABLE BACKRESTS AND SEAT CUSHIONS



All backrest cushions have adjustable filling and support wedges are also available to meet the unique postural requirements of the user

- Waterfall (standard): For users who require targeted backrest support, the waterfall option allows for all 3 x cushions to be individually positioned to best suit a user's requirements. Can also be used for kyphosis positioning.
- Classic (accessory): Traditional backrest shape with lateral adjustment pockets for gentle contouring and upright sitting support
- Advanced Immersion (accessory): Four adjustable chambers create an immersive and highly supporting backrest for more pronounced upright sitting support and accommodation of kyphosis.



Hybrid Cushion Technology: Unique foam-in-cell air cushion automatically redistributes pressure and maximises user comfort.



After-Market Cushion Compatibility: Pre-cut seat structure accommodates third-party cushions such as ROHO without compromising performance or infection control.



NEW PRODUCT LAUNCH

ASSESSMENT FORM

We have created an assessment form. This is for therapists to fill out setup and accessory options.

CHAIR SIZE	Part Number	Seat Depth (mm)	Seat Height (mm)	Backrest Height (mm)	Seat Width (mm)	SWL (kg)
Small / Medium	CHP198250	🖌 475 D1	490	740	485 W1	100
		525 D2			✓ 535 W2	160
Tall - Medium / Large	CHP198260	525 D1	490	790	535 W1	160
		575 D2			585 W2	

PRESSURE



✓ Pressure reducing foam (Standard) PR1



Small/Medium - 15 x 15in PR2 Tall Medium/Large - 17 x 17in PR2

POSTURE









Waterfall backrest PO1 Kyphosis backrest PO2 (Standard) (Standard) (Standard)

Classic style backrest PO3 Sml/Med - CHP198290 Tall Lrg/ Med - CHP 198291

Advance Immersion backrest P ✓ Medium - CHP198295 Large - CHP198296

ACCESSORIES *Items highlighted in blue



Channel legrest CL CHP198320



Lateral support blocks LS CHP198315



Profiled Headrest PH Sml/Med - CHP198310 Tall Lrg/Med - CHP198311

CUSTOM MADE TO ORDER OPTIONS



ROHO over seat option * 6 x weeks lead-time Doesn't include ROHO, option available as an upcharge

HANDSET CONTROL LOCATION



Left Hand Side HC2

Attendant Position - CHP198325 HC3



NEW PRODUCT LAUNCH

			CHP198250	CHP198260
	Seat Depth	Adjustment Required?	Sml/Med	Tall Med/Lrg
D1	Bracket at rear of backrest	Yes, refer to UM or video	475mm	525mm
D2	Bracket at front of backrest	Yes, refer to UM or video	525mm	575mm
	Seat Width	_	_	
W1	Standard position out of the box	No	485mm	535mm
W2	Adjusted position	Yes, refer to UM or video	535mm	585mm
	Pressure (seat Cushion)	_	_	
PR1	Standard cushion out of the box	No		
PR2	Fitment of after market ushion	Yes, refer to UM or video		
	Posture (Backrest Cushion)			
PO1	Standard backest out of the box	Yes, refer to UM or video		
PO2	Additional accessory	Yes, refer to UM or video		
PO3	Additional accessory	Yes, refer to UM or video		
PO4	Additional accessory	Yes, refer to UM or video		
	Accessories		_	
CL	Channel Legrest	Yes, refer to UM or video		
LS	Lateral Support Blocks	Yes, refer to UM or video		
PH	Profiled Headrest	Yes, refer to UM or video		
	Handset Control Location	_	_	
HC1	Right Hand Side, standard out of the box	No		
HC2	Left Hand Side	Yes, refer to UM or video		
HC3	Attendant (Rear Side)	Yes, refer to UM or video		

ACCESSORIES



Profiled Headrest CHP198310 - Medium CHP198311 - Large



Lateral Support Blocks CHP198315



Channel Legrest CHP198320

