

# NEW PRODUCT LAUNCH



## Aspire PostureFit Lift Recline Chair – Available Now!

This innovative seating solution is a game-changer for those with unique and changing needs.

The PostureFit chair offers a uniquely adjustable and adaptable seating solution. Through its seat width and depth adjustment as well as broad range of backrest and seat cushion options, the PostureFit chair provides unparalleled comfort and support.



CHP198250 – PostureFit – Small/Medium

CHP198260 – PostureFit – Tall Medium/Large

**[For scripting and customisation, refer to the assessment form.](#)**

# NEW PRODUCT LAUNCH

## PRODUCT FEATURES

- **Tool-free Seating Adjustment:** Therapists or carers can quickly and easily make adjustments without the need for any tools.
  - **Seat Depth:** Backrest can be quickly and easily adjusted to achieve two seat depths with a 50mm range.
  - **Seat Width:** Armrests are interchangeable for left/right side, swap sides to achieve two seat widths with a 50mm range.
- **Zero Gravity Positioning:** Creates a feeling of weightlessness and redistributes pressure off the lower back for pain management.
- **Oedema Management:** Leg elevation above the heart promotes drainage of lower limb swelling, aids circulation, and reduces fatigue.
- **Adjustable Headrest:** Electrically adjustable via the handset, the headrest cushion can be tilted forward and backward for optimal support.
- **Integrated Leg Extension:** Suits most user heights without the need to alter the legrest or height of the chair. The leg extension can also be easily removed for users with specific heel care needs
- **Vertical Lift Mechanism:** Offers sitting height adjustment for taller users and a stable sit to stand motion that removes the sensation of sliding forward for improved confidence. Available on all models without compromising SWL, unlike competing products.
- **Effortless mobility and Safety:** Large directional castors are engaged by activating a rear carer foot pedal. Castors are not in contact with the floor when in a neutral seated position, ensuring a strong base of support. This also removed any risk of sideways movement as opposed to other products with front swivelling castors
- **Transport mode:** For safety, all electrical functions are inactivated while castors are engaged to prevent the chair from being activated while mobile.

## INTERCHANGEABLE BACKRESTS AND SEAT CUSHIONS



All backrest cushions have adjustable filling and support wedges are also available to meet the unique postural requirements of the user.

1. **Waterfall (standard):** For users who require targeted backrest support, the waterfall option allows for all 3 x cushions to be individually positioned to best suit a user's requirements. Can also be used for kyphosis positioning.
2. **Classic (accessory):** Traditional backrest shape with lateral adjustment pockets for gentle contouring and upright sitting support
3. **Advanced Immersion (accessory):** Four adjustable chambers create an immersive and highly supporting backrest for more pronounced upright sitting support and accommodation of kyphosis.



**Hybrid Cushion Technology:** Unique foam-in-cell air cushion automatically redistributes pressure and maximises user comfort.



**After-Market Cushion Compatibility:** Pre-cut seat structure accommodates third-party cushions such as ROHO without compromising performance or infection control.

# NEW PRODUCT LAUNCH

## ASSESSMENT FORM

We have created an assessment form (attached to this release). Provide this to your therapists for setup options and all optional accessories.

- Each section of this form has a linking code for the chair setup
- These alphanumeric codes allow for chair set-up preferences to be shown as notes in GP (for example overleaf). It is up to the person entering the order to add this into the comments section of the order, the warehouse will know what position to set the chair up in at the customer's site
- Chairs are to be shipped out in boxes and setup onsite, not in warehouses
- For setup instructions, please refer to the [VIDEO](#)

CHAIR SIZE	Part Number	Seat Depth (mm)	Seat Height (mm)	Backrest Height (mm)	Seat Width (mm)	SWL (kg)
Small / Medium	CHP198250	<input checked="" type="checkbox"/> 475 D1	490	740	<input type="checkbox"/> 485 W1	160
		<input type="checkbox"/> 525 D2			<input checked="" type="checkbox"/> 535 W2	
Tall - Medium / Large	CHP198260	<input type="checkbox"/> 525 D1	490	790	<input type="checkbox"/> 535 W1	160
		<input type="checkbox"/> 575 D2			<input type="checkbox"/> 585 W2	

### PRESSURE



Pressure reducing foam (Standard) PR1

Small/Medium - 15 x 15in PR2  
 Tall Medium/Large - 17 x 17in PR2

### POSTURE



Waterfall backrest PO1 (Standard)



Kyphosis backrest PO2 (Standard)



Classic style backrest PO3  
 Sm/Med - CHP198290  
 Tall Lrg/ Med - CHP 198291



Advance Immersion backrest PO4  
 Medium - CHP198295  
 Large - CHP198296

### ACCESSORIES \*Items highlighted in blue



Channel legrest CL  
 CHP198320



Lateral support blocks LS  
 CHP198315



Profiled Headrest PH  
 Sm/Med - CHP198310  
 Tall Lrg/Med - CHP198311

### CUSTOM MADE TO ORDER OPTIONS



ROHO over seat option \* 6 x weeks lead-time  
 Doesn't include ROHO, option available as an upcharge

### HANDSET CONTROL LOCATION

Right Hand Side HC1

Left Hand Side HC2

Attendant Position - CHP198325 HC3

# NEW PRODUCT LAUNCH

## LINKING CODE GLOSSARY

- All staff are required to watch the setup video on the chair to understand what each linking code refers to
- Linking codes are only done for ease of information delivery from script form to GP

			CHP198250	CHP198260
Seat Depth		Adjustment Required?	Sml/Med	Tall Med/Lrg
D1	Bracket at rear of backrest	Yes, refer to UM or video	475mm	525mm
D2	Bracket at front of backrest	Yes, refer to UM or video	525mm	575mm
Seat Width				
W1	Standard position out of the box	No	485mm	535mm
W2	Adjusted position	Yes, refer to UM or video	535mm	585mm
Pressure (seat Cushion)				
PR1	Standard cushion out of the box	No		
PR2	Fitment of after market ushion	Yes, refer to UM or video		
Posture (Backrest Cushion)				
PO1	Standard backrest out of the box	Yes, refer to UM or video		
PO2	Additional accessory	Yes, refer to UM or video		
PO3	Additional accessory	Yes, refer to UM or video		
PO4	Additional accessory	Yes, refer to UM or video		
Accessories				
CL	Channel Legrest	Yes, refer to UM or video		
LS	Lateral Support Blocks	Yes, refer to UM or video		
PH	Profiled Headrest	Yes, refer to UM or video		
Handset Control Location				
HC1	Right Hand Side, standard out of the box	No		
HC2	Left Hand Side	Yes, refer to UM or video		
HC3	Attendant (Rear Side)	Yes, refer to UM or video		

## ACCESSORIES



**Profiled Headrest**

CHP198310 - Medium  
CHP198311 - Large

**Lateral Support Blocks**

CHP198315

**Channel Legrest**

CHP198320