

Aspire Lift Chairs - Handest on Left or Right

We are pleased to announce that there is now provision on all Aspire Lift Chairs to position the handset on the left or right-hand side, currently, all chairs are supplied with the handset on the right-hand side (from the position of a user sitting within the chair).

The ability to change sides is intended to assist clients who may have limited cognition, movement or dexterity on either side of their body.

This changeover can be achieved with a small modification under the chair. An extension cable (details are given below) will need to be added for some chairs to move the handset from the right-hand side to the left-hand side of the chair.

Contact your local Aidacare store for more information.



