



22/10/21

## ASPIRE STAR PATIENT MOVER – NOW AVAILABLE!

The Aspire Star Patient Mover (the Star Mover) allows patients to be transferred from beds and chairs with comfort and safety. The compact and versatile design is suitable for Homecare, Institutional and Aged Care environments.

The Star Mover replaces the Aspire Swift Mover but with many design improvements

**Product Code: LSS390580, replacing previous Swift Mover (LSS390590).**

### KEY FEATURES

- **Lower footplate** for easier patient entry and exit.
- **Wider space between seat pads** improves access for a greater range of users.
- **Narrower overall width** improves accessibility in tight spaces.
- **Larger 100mm** twin tante castors reduce rolling resistance for improved manoeuvrability.
- **Height adjustable knee support** improves comfort and is adaptable to a larger range of user heights.
- Low maintenance moving parts for ease of servicing.
- **Knock-Down design** for improved freight cost management.
- Tapered closed leg position to access between chair legs.



New Product Launch