

FLX BED TRANSPORT TROLLEY RELEASE!

The New EasyMove Transport Trolley for the FLX beds (And FL250) is now in stock! - BEA030055

The trolley is near identical to the FL250 version but now accommodates both FLX and FL250 beds while also allowing the bed to be locked and transported at intermediate angles which means side rails or side trims do not have to be removed before transporting.

Note: Due to the angle changes and the need to make the trolley compatible with both FLX and FL250 beds, the transport trolley comes as a pair with a LEFT and RIGHT sided trolley.



