

Hybrid Comfort Bariatric Mattress – Now Available!

We are happy to announce the addition of a Bariatric Hybrid Mattress in King Single size, as presented in the August national product presentation.

This new model ensures that individuals with higher weight requirements can also experience the same benefits and comfort provided by our Hybrid Comfort Mattress. The Hybrid Bariatric Mattress has a safe working load of 380kg, ensuring the highest level of safety and support for bariatric care.



Hybrid Comfort Bariatric Mattress – PTM638007

Hybrid Comfort Bariatric Pump – PTM638032

PRODUCT FEATURES

- **High SWL:** Up to 380kg SWL, ensuring true bariatric care.
- **Immersive Visco-elastic Memory Foam Surface:** Thick immersive top layer redistributes pressure and optimizes user comfort
- **Intelligent Hose Management:** When not in use the air-hose can be stored internally, reducing infection control and OH&S risks
- Integrated Shear Reduction: Memory foam layer moves independently of the cell surface reducing shear when bed is adjusted
- **Optimised Heel Protection:** Tapered heel zone redistributes pressure to lower leg, away from the vulnerable heel region. The castellated surface collapses in a wave-like movement reducing shear

