

MALVESTIO ARIA ICU BED - Now Available!

The Aria ICU bed comes with Integrated Drive Wheel, Weigh Scale, Patient Alarms and Lateral Turning as standard.

Modest stock quantities of the Aria Bed will be located in distribution centres nationally with access for the Acute Care Team. As ICU beds are a large capital purchase, we will always hold a small amount of stock for demo/trial/one-off purchases with larger quantities available on a four-month lead time.

BEB050160 - Malvestio Aria ICU Bed











NEW PRODUCT LAUNCH

Key Features

- Lateral Turning: Programmable lateral rotation up to 18°controlled via the touch screen.
- **Incentro System:** Combined back and pelvic translation nullifies horizontal forces, prevents slipping, and keeps the body at a 30° angle, optimal for ventilation
- **DeltaDrive:** Electric centre 5th wheel with integrated controls in the headboard provides increased maneuverability for hospital staff.
- Weigh Scale: An integrated system allows for precision weight measurements.
- Integrated Patient Alarms:
 - o **Restless Patient:** Detect different levels of patient movement, alerting the operator if the patient becomes restless.
 - o **Predictive Exit:** By detecting the patient's movements, it is possible to receive an alert when the patient is about to leave the bed unattended.
 - o **Patient Exit:** The weight monitoring system can be programmed to alert if the patient leaves the
- **Biocote:** metal components are epoxy powder coated and finished with an antibacterial additive for improved infection control.
- Safe Exit Lights: 2 x LED lights signal whether the bed is at a safe exit height, promoting safe egress for the patient (orange is unsafe, and green is safe).
- Four Section Safety Sides with Electric Adjustment: Integrated bed movement controls in side panels allow for quick and easy activation by the patient or carer.





