



11/11/21

## ASPIRE PRESSURE REDUCING COMMUNE SEATS – UPDATE

We have made some upgrades to our already successful Aspire Pressure Reducing Commune Seats.

These seats have been made in QLD at Geebung over the past few years but we have moved production to our upholstery partners in Sydney, freeing up the Geebung team to focus on Complex Rehab and for stock flow of the commune seats to be streamlined.

Most sizes will be held in stock however lower volume and custom sizes will be made to order only (via the Custom Made order form).

### KEY DIFFERENCES & IMPROVEMENTS

- **Carrflex multi stretch fabric and Australian made foam** reduces shear and promotes pressure redistribution (already used on previous seats).
- **Completely sealed solution** for maximum infection control management. Now made with an internal routed board that allows fabric to be stretched around the seat, stapled and then sealed with an anti-mould silicone, totally preventing any fluid ingress.
- **Custom Made** solutions are still available. Please refer to the updated (now interactive) 'Custom Made' order form [attached and also available on the intranet and Aidacare website, see page 2 for details]. Please note this is now to be emailed to [customfurniture@aidacare.com.au](mailto:customfurniture@aidacare.com.au) with a copy of the sales order, to be confirmed and processed.



### CUSTOM MADE SHOWER COMMUNE SEATS



**CLIENT** Name: \_\_\_\_\_ Date: \_\_\_\_\_  
 Phone: \_\_\_\_\_ Email: \_\_\_\_\_  
 Address: \_\_\_\_\_ Fax: \_\_\_\_\_

**THERAPIST** Name: \_\_\_\_\_ Email: \_\_\_\_\_  
 Phone: \_\_\_\_\_ Fax: \_\_\_\_\_

**AIDACARE REPRESENTATIVE** Name: \_\_\_\_\_ Email: \_\_\_\_\_

Please complete all relevant details below in 1 to 6 and send completed script form to [customfurniture@aidacare.com.au](mailto:customfurniture@aidacare.com.au)

1 Please choose one:  
 Order - Please Provide Order Number: \_\_\_\_\_ or  Quote Request

2 Which brand of shower commode will the seat be attached to? (Required for attachment type):  
 \_\_\_\_\_

3 Choose your Fabric

PU Coated 4 way Stretch Fabric (Standard seat thickness 70mm with 30mm memory foam top layer) 

Medical Grade White Vinyl (Standard seat thickness 50mm) 

4 Choose your Seat Configuration (MUST CHOOSE ONE)

Open Front  ↑  Open Right  ↑

Closed Front  ↑  Open Left  ↑



Product Update





Code	Description	Notes
BTC089130	Aspire PR Commode Seat - 460mm - Open Front	Stocked
BTC089140	Aspire PR Commode Seat - 460mm - Closed Front	Stocked
BTC089180	Aspire PR Commode Seat - 500mm - Closed Front	Stocked
BTC089170	Aspire PR Commode Seat - 500mm - Open Front	Stocked
BTC089190	Aspire PR Commode Seat - 550mm - Open Front	Stocked
BTC089200	Aspire PR Commode Seat - 550mm - Closed Front	Stocked
BTC089210	Aspire PR Commode Seat - 600mm - Open Front	Stocked
BTC089220	Aspire PR Commode Seat - 600mm - Closed Front	Stocked
BTC089250	Aspire PR Commode Seat - Customised (as per script)	Made to Order
BTC089122	Aspire PR Commode Seat - 320mm - Open Front	Made to Order
BTC089240	Aspire PR Commode Seat - 650mm - Closed Front	Made to Order
BTC089124	Aspire PR Commode Seat - 320mm - Closed Front	Made to Order
BTC089126	Aspire PR Commode Seat - 400mm - Open Front	Made to Order
BTC089150	Aspire PR Commode Seat - 460mm - Open Left	Made to Order
BTC089128	Aspire PR Commode Seat - 400mm - Closed Front	Made to Order
BTC089160	Aspire PR Commode Seat - 460mm - Open Right	Made to Order
BTC089230	Aspire PR Commode Seat - 650mm - Open Front	Made to Order

## CUSTOM FORM LOCATIONS

- Aidacare website (uploaded but not live until 12/11/21)
  - [Aspire Custom Commode Seat - Website](#)

