

**Product Awareness – November 2019**

**Model Update – Aspire Vogue Lightweight 2 Walkers**

Out with the old and in with the new! Our Aspire Vogue Lightweight Walkers are to be replaced by an updated model which has now arrived in stock!

The NEW Aspire Vogue Lightweight 2 has a sleek new frame design and other improved features such as:

- Wider backrest (now 190mm vs 80mm);
- Deeper seat (now 240mm vs 185mm);

Both improving the support and comfort provided by these surfaces.

- Reduced product weight from 6.7kg to 6.2kg to assist handling and manoeuvring.

PLEASE NOTE – the improved wider backrest and deeper seat will also be supplied on the Vogue Carbon Fibre and Vogue Forearm models from JAN 2020, which otherwise remain unchanged in their designs.

All other features of the prior model are retained, including folding clip lock, removable utility bag, walking stick holder, super soft wheels, locking hand brakes, easy handle height adjustment and 150kg SWL.

The White colour is out and replaced with Champagne due to strong customer feedback.

**OLD MODEL**

WAF705300RE – Vogue Lightweight (Red)

WAF705300WH – Vogue Lightweight (White)



**NEW MODEL**

WAF705310RE – Vogue Lightweight 2 (Red)

WAF705310CH – Vogue Lightweight 2 (Champagne)

