

12/10/20

ASPIRE ACTIVE AIR ACUTE 8 MATTRESS – NOW WITH TURNASSIST!

Our Active Air Acute 8 mattress systems now offer industry leading TURNASSIST (TM) technology as standard! Easy to operate valves located on the side of the mattress (foot end) engage lateral turning through integrated air cells within the mattress.

Please help to clear old stock (PTM639400 & PTM639410) prior to selling new stock.

FEATURES

- o Lateral turning through integrated air cells, allowing the mattress to turn a patient without the need to add an after-market turning device
- o This patented innovation is unique and saves customers thousands of \$\$ normally required for automated turning therapy
- o Drastically reduces the moving and handling stresses on Care Staff by rotating the patient up to 30 degrees. This technology is ideal for:
 - o Alleviating asymmetric or laterally located pressure areas
 - o Simplified in-bed bathing
 - o Rapid changing of linen
 - o Perfect positioning of manual handling slings and slide sheets



Aspire Active Air 8 Acute w/turnassist
PTM639400T – Single
PTM639410T – King Single

