



07/10/20

CANE HOLDER FOR ASPIRE SEAT WALKERS - UPDATE

Problem:

Existing cane holders are too short, causing the cane to fowl/collide with the hand wheels & frames on some Aspire Seat Walkers.

Walkers Affected:

1. Aspire Deluxe Seat Walker
2. Aspire Classic 8 Seat Walker
3. Aspire Classic 6 Seat Walker
4. Aspire Mini Seat Walker
5. Aspire Push Down Seat Walker

Solution:

The cane holder has been extended, allowing the cane to be mounted without colliding with the frame.

TECHNICAL CHANGE

Left: Original Production**Right: Improved Production**

STOCK PLANNING

Codes have been updated as below:

WAA686924 – Suits Deluxe, Classic 6/8, Mini & Push Down Walkers

WAA686925 – Suits Superlite & XL Walkers

