



08/04/21

## LIFECOMFORT MATTRESS COVERS UPGRADES!

We are pleased to announce some important upgrades to our LIFECOMFORT Foam Mattress covers. All LIFECOMFORT® Foam Mattresses are now fitted with superior hospital-grade polyurethane coated covers designed to better withstand rigorous use and cleaning in the Home Care, Aged Care and Acute Care environments.

**IMPROVED FEATURES** 

- o ALL LIFECOMFORT® Foam Mattresses will now offer superior covers that include:
  - Hospital-grade polyurethane covers suitable for bleach concentrations up to 10,000ppm
  - Welded waterfall covers for improved infection control
  - o Carry handles for safe manual handling of mattresses
  - Four-sided zippers for easy removal, cleaning and refitting of covers
  - o Asset management labels for facility equipment management
  - Rotational turning graphics to assist mattress rotation policies



LIFECOMFORT® Foam Mattresses
\*Refer to brochure for codes & sizes



