

SCOOTER SELECTION GUIDE

	Scooter Code	lmage
	BOOT SCOOTER (SCT670613)	
JIMPACI	GO-GO® LX + CTS SCOOTER (SCT670612)	
	SONIC SCOOTER (SCT670665)	Comp.
5	QUEST FOLDING COMPACT SCOOTER (SCT670630)	
	BOOT SCOOTER SUPALITE (SCT670616)	
MEDIUM	GT SCOOTER (SCT670670)	⊘ DVA
	PATHRIDER® 10 SCOOTER (SCT670615)	⊘ DVA
	PATHRIDER® 10 DX SCOOTER (SCT670620)	⊘ DVA
	VICTORY 10 LX SCOOTER (SCT670640)	⊘ DVA
	GOPHER DX SPORT SCOOTER (SCT670638)	
	PATHRIDER® 130XL PURSUIT SCOOTER (SCT670571)	⊘ DVA
	PATHRIDER® 140 XL SCOOTER (SCT670576)	⊘ DVA
AKGE	CELEBRITY® DX XL SCOOTER (SCT670540)	⊘ DVA
	SPORTRIDER SCOOTER (SCT670585)	
	RANGERIDER SCOOTER (SCT670581)	

Terrain	Battery	Max Range Per	Safe Working
		Charge*	Load
Level & smooth terrain. Short distances.	2x 12V 12AH	11.5 km	118 kg
Level & smooth terrain. Short distances.	2x 12V 18AH	17 km	136 kg
Level & smooth terrain. Short distances.	2 x 12V 20AH	20 km	136 kg
Able for more varied terrain. Short to medium distances.	1x 24V 11.5AH Li-Polymer	20 km	115 kg
Level & smooth terrain. Short distances.	1x 24V 12AH	14.4 km	120 kg
Able for more varied terrain. Short to medium distances.	2 x 12V 35AH	30 km	159 kg
Able for more varied terrain. Short to medium distances.	2x 12V 31AH (40AH optional)	40 km	159 kg
Able for more varied terrain. Short to medium distances.	2x 12V 31AH (40AH optional)	33 km	159 kg
Can handle rougher terrain well. Long distances.	2x 12V 40AH	26 km	181 kg
Able for more varied terrain. Short to medium distances.	2x 12V 50AH	33 km	160 kg
Can handle rougher terrain well. Long distances	2x 12V 55AH	25 km	159 kg
Can handle rougher terrain well. Long distances	2x 12V 70AH	35 km	181 kg
Can handle rougher terrain well. Long distances	2x 12V 40AH	40 km	159 kg
Can handle rougher terrain well. Long distances	2x 12V 70AH	50 km	181 kg
Can handle rougher terrain well. Long distances	2x 100AH	30 km	15 9 kh

Features	
- Up to 6km/h - Disassembles for easy transport	

- Up to 6.4km/h - Disassembles for easy transport - Comfort Trac front & rear suspension (CTS)

> - Up to 6.4km/h - Disassembles in easy steps

- Up to 8.8km/h - Great outdoor stability, ease of use and maximum value

> - Up to 6.5km/h - Disassembles for easy transport - Easy-to-grip tiller adjustment knob

> > - Up to 10km/h - Front & rear suspension

- Up to 9km/h - Comfort Trac front & rear suspension (CTS)

- Up to 8.5km/h - Delivers high performance and feather-touch disassembly - Can upgrade to 40AH batteries (x2) for improved performance

- Up to 9km/h - Front & rear suspension - Can upgrade to 70AH batteries (x2) for improved performance

> - Up to 10km/h Front & rear suspension - High visibility, front and rear lighting package

- Up to 10km/h Front & rear suspension - Can upgrade to 100AH batteries (x2) for improved performance

- Up to 10km/h - Front & rear suspension - Can upgrade to 100AH batteries (x2) for improved performance

- Up to 10km/h - Front & rear suspension - Can upgrade to 100AH batteries (x2) for improved performance

- Up to 10km/h Front & rear suspension - Can upgrade to 100AH batteries (x2) for improved performance

- Up to 10km/h - Front & rear suspension

- CTS Suspension includes adjustable shocks for greater comfort - Backup sensors with LED light

*Range Per Charge varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition & tyre condition. Battery upgrades (where possible) improve Range Per Charge.

SCOOTER MOBILITY

At the heart of Pride's philosophy is the understanding that mobility products empower lives. Through industry-leading, consumer-inspired designs, Pride offers an innovative, distinctive and performance driven comprehensive line of mobility products to help our customers lead the most active, independent lifestyles.

Aidacare works in close association with Pride Australia to promote and support Pride mobility products

- Pride Scooter products are all registered in the ARTG
- Please always read and understand your owners handbook before operating your new mobility scooter
- Please take care to observe safety and the speed limit requirements of your state when driving your mobility scooter
- Specification and other details on this information sheet may be subject to change without notice









5 LED headlight, tail light and indicators indicators on some models

6 Low profile, puncture proof tyres or pneumatic tyres



SAFE USE VIDEO

At Aidacare, we value the safety of our customers above all else. There are a few safety considerations we'd like you to understand before operatina vour device.

We strongly recommend all new mobility device customers watch our 7 minute YouTube video, 'Getting to Know your Power Mobility Device' that covers the above 8 points in detail and more





MOBILITY SCOOTER SELECTION GUIDE











FEATURES

- Feather touch disassembly into 5 pieces for easy transportation & storage
- Comfortable front-mounted post seat with arm rests for stability Easy to remove rechargeable
- Soft and easy to adjust tiller
- Modular design for easy servicing

SPECIFICATIONS

Colour: Length: 1003 mm Width: 495 mm Red Product Weiaht: 36.6 kg Maximum Speed: 6.4 km/h Range: 11.5 km HP: 1.6 Battery Size: 2 x 12V (12 AH)

GO-GO® LX + CTS **SCOOTER** SCT670612



aspire

SONIC SCOOTER

FEATURES

- CTS Suspension (Comfort-Trac Suspension) brings a new level of performance to travel mobility
- Lightweight mobility scooter • Disassembles into 5 pieces for easy
- transportation and storage • Easy to remove rechargeable battery
- Soft and easy to adjust tiller

SPECIFICATIONS

Length: 1040 mm Colour Width: 520 mm Product Weight: 39.5 kg Maximum Speed: 6.4 km/h Range: 17 km HP: 1.6 Battery Size: 2 x 12V (18AH)





- Portable and durable scooter
- Puncture-proof tyres
- Detachable battery pack • Front and rear LED light
- Easy access battery compartment • Adjustable angle tiller with delta handle

SPECIFICATIONS

Length: 1040 mm Colour Width: 495 mm Product Weight: 45 kg Maximum Speed: 6.8 km/h Range: 20 km HP: 0.241 Battery Size: 2 x 12V 20AH

QUEST FOLDING COMPACT SCOOTER



aspire

GT SCOOTER

SCT670670

PATHRIDER® 10

SCOOTER

SCT670615

FEATURES

MOBILITY SCOOTER SELECTION GUIDE

- Sleek, stylish and extremely easy to operate
- Ready for transport and folded in just 3 easy steps
- Can be stowed upright or laid down Great performance both indoors and out
- Articulated front end
- Front and rear LED headlights

SPECIFICATIONS

Length: 930 mm Colour: As Shown Width: 485 mm Product Weight: 23.5 kg Maximum Speed: 6 km/h Range: 20 km HP: 0.362 Battery Size: 1 x 24V (11.5 AH) Li-Polymer



• Offers full advanced suspension for a smoother, more comfortable ride over varied terrain

FEATURES

- Full advanced suspension
- Front and rear LED lighting illuminates your path
- Deluxe tiller with USB charger and cup holder

SPECIFICATIONS

Lenath: 1200 mm Colour Width: 605 mm Product Weight: 82 kg Maximum Speed: 10 km/h Range: 30 km HP: 0.536 Battery Size: 2 x 12V (35AH)

FEATURES



- · High performance, sleek and sporty mobility scooter
- Low profile, puncture proof tyres • Back lit battery gauge
- Comfortable front-mounted post seat with armrests for stability
- Eight-position seat rotation • Front basket as standard
- Wrap around easy drive tiller

SPECIFICATIONS

Length: 1190 mm Colour Width: 565 mm • Blue Product Weight: 63.5 kg Maximum Speed: 8.5 km/h Range: 21 / 25 km HP: 1.8 Battery Size: 2 x 12V (31AH) (40AH optional)

SCOOTER SCT670640

VICTORY 10 LX



aspire

PATHRIDER® 130XL

PURSUIT SCOOTER

SCT670571

SPECIFICATIONS

Colour: Width: 565 mm Battery Size: 2 x 12V (40AH)

FEATURES



- rests, head rest and lumbar support • Rugged front and rear suspension

1300 133 120

www.aspirecare.com.au

• Long lasting high intensity LED kerb light, headlight, tail light and indicators

Width: 635 mm Product Weight: 73 kg Maximum Speed: 12.4 km/h **Range:** 25 km **HP:** 3.1 Battery Size: 2 x 12V (55AH)

PATHRIDER® 140 XL **SCOOTER**





CELEBRITY® DX XL

SCOOTER

SCT670540

- **FEATURES**
- Large 13" pneumatic tyres • Luxurious suspension seat with arm
- Front and rear suspension and sealed
- Automatic full lighting system


Length: 1428 mm Width: 692 mm

Maximum Speed: 10 km/h • Silver **Range:** 35 km **HP:** 4.6

- arm rests for stability
- State of the art rear suspension



Length: 1300 mm Width: 640 mm Product Weight: 91 kg Maximum Speed: 10 km/h Range: 40 km HP: 2.83 Battery Size: 2 x 12V (40AH)




- speedometer/odometer

aspire

- Powerful front/rear lighting kit

 Black Battery Size: 2 x 12V (70AH)





- hydraulic brakes
- Back lit battery gauge
- Wrap around easy drive tiller



 Red Product Weight: 118 kg • Yellow

Battery Size: 2 x 12V (70AH)




- Deluxe front-mounted high back swivel seat with headrest and
- Lumbar support and manual recline
- · Unique easily adjusted steering mount with modern LED lighting





FEATURES



- Full suspension system
- High back, adjustable seat
- Large rear storage pod included as standard



Colour



aspire

BOOT SCOOTER

SUPALITE

SCT670616

FEATURES Compact, foldable seat with armrests

- Easy-to-grip tiller adjustment knob
- Stable, front frame-mounted seat post A standard front basket provides
- Easy to transport at home, hotels, or



Length: 1022 mm Width: 490 mm Product Weight: 8.5 kg

Maximum Speed: 6.5 km/h Range: 14.4 km HP: 8.9 Battery Size: 24V (12AH)

FEATURES

- Speed reduction technology automatically reduces speed when turning corners Independant front suspension, mono
- shock rear suspension • Easy-drive tiller with wraparound
- · High visibility, front and rear lighting
- package



Length: 1370 mm Colour: Width: 670 mm Product Weight: 102.5 kg Maximum Speed: 10 km/h

Range: 33 km HP: 20.62 Battery Size: 12V (50AH)

FEATURES

seat with lumber and pillow-top headrest, flip-up height, and angle adjustable armrests

shocks for greater comfort

 Backup sensors with LED light Dual hydraulic brakes



Length: 1560 mm •Stone Grev







PATHRIDER® 10 DX

SCOOTER

SCT670620

SPECIFICATIONS

FEATURES

High performance sleek and

Low profile, puncture proof tyres

seat with armrests for stability

Comfortable front-mounted post

• Back-lit battery gauge on monitor

• Rugged front and rear suspension

sporty mobility scooter

Length: 1193 mm Colour: Width: 565 mm Product Weiaht: 85 kg Maximum Speed: 9 km/h Range: 21 / 25 km HP: 1.8 Battery Size: 2 x 12V (31AH) (40AH optional)

FEATURES CTS Suspension (Comfort-Trac front

simple frame separation with only one

- and rear suspension) • Sleek, modern design Feather-touch disassembly permits
- Full LED lighting package that includes curb and rear lights and directional signals



Product Weight: 71 kg Maximum Speed: 10 km/h **Range:** 26 km **HP:** 0.570





Automatic full lighting system















Length: 1650 mm Width: 780 mm Product Weight: 120.5 kg Maximum Speed: 12 km/h Range: 50 km HP: 4.1



















Width: 830 mm Maximum Speed: 10 km/h Range: 30 km HP: 18.6 Battery Size: (2) (100AH)

