

ASPIRE REHAB RS

MASS QLD ORDER FORM SUMMARY PAGE

Effective:

ORDERING INFORMATION: Phone 1300 133 120

| DEALER INFORMATION | | |
|--------------------------------------------------------------------|-----------------------------------------------|--|
| Date Account Code Sales Rep Mobile Signature | City / State | |
| CUSTOMER INFORMATION | | |
| Customer Name MASS Plaque # MASS PO# Signature Rep Name & Ph # | Address City / State Postcode | |
| PRESCRIBER INFORMATION | | |
| Prescriber NameMobile | Address City / State | |
| Signature | Postcode ———————————————————————————————————— | |

PRODUCT INFORMATION

The REHAB RS is an ideal solution for Aged Care and Community users who requires increased postural support and bodyweight redistribution. A reclined position is easily achieved, making the REHAB RS also a great option for those who fatigue easily.

| Max. User Weight | S.W.L - 180kg | |
|----------------------------------------------------------------------------|--------------------------------------|--|
| Chair Weight | 40.7 - 46.7kg | |
| Seat Height - adjustable (floor to top of fabric) | Front: 410 - 440mm Rear: 400 - 475mm | |
| | 410 - 460mm (16-18") | |
| Seat Width | 460 - 510mm (18-20") | |
| | 510 - 560mm (20-22") | |
| Seat Depth - adjustable | 440 - 530mm, 470 - 560mm | |
| Backrest Height - adjustable (top of seat to top of backrest handle) | 460 - 710mm | |
| Overall Width | 690mm 740mm 790mm | |
| Armrest Height (top of seat to top of armrest pad) | 310 - 460mm | |
| Approvals | AUST R ARTG 291393 | |
| | Crash tested to AS/NZS3696.19 | |





ASPIRE REHAB RS

MASS QLD ORDER FORM

Effective:

ORDERING INFORMATION: Phone 1300 133 120

| ACCOUNT INFORMATION | |
|---------------------|------------|
| Invoice to | Deliver to |
| Account # | Name |
| Address | Address |
| | |
| Email | Email |
| Phone | Phone |
| Order Date | Marked for |
| PO # | Order # |
| Number Required | |

Removable, Height & Depth Depth, Width & Height Colour: Silver Armrest: Seat: Adjustable Adiustable Quick Release, Solid 24" Tyres, Height, Tension & Angle Rigid Steel Frame Wheels: **Backrest:** Frame:

Height & Depth Adjustable Adjustable Swing-away, Removable,

PU Breathable Multi-Stretch Legrest: Padded & Elevating Castors: 6", Height Adjustable **Upholstery:** Fabric Covers

Articulated

Multilayer Castallated Back Attendant Operated Drum **Footrest:** Height & Angle Adjustable Brakes & Push-to-lock **Cushion:** and Seat Cushion Covered in **Brakes:** Breathable PU Fabric

Handbrakes

FOR OPTIONAL EXTRAS, PLEASE CHOOSE ONE OF THE PACKAGES BELOW:

All customisations will be quoted as labour in addition to the base price, customisations can be carried out by Aidacare or by an authorised Aidacare Network Partner.

Seat Width 410 - 460mm 460 - 510mm 510 - 560mm

Code **Elevating Legrests** Price Qty Code **Drum Brakes** Price Qty WCA848100 Left - 410mm WCA848510 Rear Wheels w/Drum Brakes Kit - 16in WCA848110 Right - 410mm WCA848520 Rear Wheels w/Drum Brakes Kit - 20in WCA848120 Left - 460/510mm Code Depth Extensions Price Qty WCA848130 Right - 460/510mm WCA848500 WCA848530 Base Frame Extension Kit - 50mm - Black Rear Wheels w/Drum Brakes Kit - 24in WCA848270 Anti-Tip Bars WCA848540 Base Frame Extension Kit - 100mm - Black **Articulating Elevating Legrests** Code Qty Code Other Price Qty Price WCA848160 Left - 410mm WCA848550 Therapist-Friendly Adjustment Kit Right - 410mm WCA848170 WCA848560 Meal Tray - 320mm to 410mm WCA848180 Left - 460/510mm WCA848570 Meal Tray - 460mm to 510mm WCA848190 Right - 460/510mm WCA848280 IV Pole Code **Stump Supports** Price Qty WCA848300 Safety Lap Belt - Metal Seatbelt Buckle WCA848220 Left - Silver - 410/460mm WCA848580 Calf Strap WCA848230 Right - Silver - 410/460mm WCA848240 Left - Silver - 510mm WCA848250 Right - Silver - 510mm Code Oxygen Bottle Holder Price Qty WAA695600 Oxygen Bottle Bag MASS Plaque No. Code **Chest Supports** Price Qty MASS Funded WCA848450 Left Client Contribution WCA848460 Right

TOTAL PRICE