

ASPIRE REHAB RX JNR

MASS QLD ORDER FORM SUMMARY PAGE

Effective:

ORDERING INFORMATION: Phone 1300 133 120

DEALER INFORMATION		
Date Account Code Sales Rep Mobile Signature	City / State Postcode Fax	
CUSTOMER INFORMATION		
Customer Name MASS Plaque # MASS PO# Signature Rep Name & Ph #	Address City / State Postcode	
PRESCRIBER INFORMATION		
Prescriber NameMobile	Address City / State	
Signature	Postcode Phone Fax Email	

PRODUCT INFORMATION

The Aspire REHAB RX JNR is designed for advanced postural support and pressure care for children via easy repositioning and advanced adjustability.

Max. User Weight	S.W.L - 180kg
Chair Weight	39.5 kg
Seat Height - adjustable (floor to top of fabric)	Front: 410 - 440mm Rear: 400 - 475mm
Seat Width	320 - 370mm (13-15")
Seat Depth - adjustable	440 - 490mm
Backrest Height - adjustable (top of seat to top of backrest handle)	460 - 710mm
Overall Width	600mm
Armrest Height (top of seat to top of armrest pad)	310 - 460mm
Approvals	AUST R ARTG 291393 Crash tested to AS/NZS3696.19





ASPIRE REHAB RX JNR

MASS QLD ORDER FORM

Effective:

ORDERING INFORMATION: Phone 1300 133 120

ACCOUNT INFO	RMATION	
Invoice to		Deliver to
Account #		Name
Address		Address
Email		Email
Phone		Phone
Order Date		Marked for
PO #		Order #
Number Required		

EVOKE STANDARD MASS PACKAGE

Depth, Width & Height Removable, Height & Depth Frame: Rigid Steel Frame Armrest: Seat: Adjustable Adjustable Swing-away, Removable, Quick Release, Solid 20" Tyres, Height, Tension & Angle Padded & Elevating Wheels: **Backrest:** Legrest: Height & Depth Adjustable Adjustable Articulated PU Breathable Multi-Stretch **Upholstery:** Footrest: Height & Angle Adjustable **Castors:** 6", Height Adjustable Fabric Covers Multilayer Castallated Back Height, Depth, Angle Attendant Operated Drum & Width Adjustable, Brakes & Push-to-lock **Headrest: Cushion:** and Seat Cushion Covered in **Brakes:** Removable Handbrakes Breathable PU Fabric

FOR OPTIONAL EXTRAS, PLEASE CHOOSE ONE OF THE PACKAGES BELOW:

ACCESSORIES - PLEASE SPECIFY QUANTITY

Code	Standard Included Feature		
WCA848160 WCA848170 WCA848520 WCA848270	Elevating Legrest - Articulating - Left - 320mm Elevating Legrest - Articulating - Right - 320mm Rear Wheels w/Drum Brakes Kit - 20in Anti-Tip Bars		
Code	Elevating Legrests	Price	Qty
WCA848100	Left - 320mm		
WCA848110	Right - 320mm		
Code	Stump Supports	Price	Qty
WCA848220	Left - Silver - 320mm		
WCA848230	Right - Silver - 320mm		

Code	Drum Brakes	Price	Qty
WCA848510	Rear Wheels w/Drum Brakes Kit - 16in		
Code	Other	Price	Qty
WAA695600	Oxygen Bottle Bag		
WCA848550	Tool-Free Adjustment Kit		
WCA848560	Meal Tray - 320mm to 410mm		
WCA848280	IV Pole		
WCA848300	Safety Lap Belt - Metal Seatbelt Buckle		
WCA848580	Calf Strap		

MASS Plaque No.	
MASS Funded	
Client Contribution	
TOTAL PRICE	