

ASPIRE REHAB RX

MASS QLD ORDER FORM SUMMARY PAGE

Effective:

ORDERING INFORMATION: Phone 1300 133 120

DEALER INFORMATION		
Date Account Code Sales Rep Mobile Signature	Name Address City / State Postcode Fax Email	
CUSTOMER INFORMATION		
Customer Name MASS Plaque # MASS PO# Signature Rep Name & Ph #	Address City / State Postcode	
PRESCRIBER INFORMATION Prescriber Name Mobile Signature	Organisation Address City / State Postcode Phone Fax	
	Email	

The REHAB RX is geared towards the high needs, rehabilitation environment. Highly configurable and with pressure care cushioning surfaces, the REHAB RX is designed for the user with more complex needs.

Max. User Weight	S.W.L - 180kg
Chair Weight	41 - 47.2 kg
Seat Height - adjustable (floor to top of fabric)	Front: 410 - 440mm Rear: 400 - 475mm
Seat Width	410 - 460mm (16-18") 460 - 510mm (18-20") 510 - 560mm (20-22")
Seat Depth - adjustable	440 - 530mm, 470 - 560mm
Backrest Height - adjustable (top of seat to top of backrest handle)	460 - 710mm
Overall Width	690mm 740mm 790mm
Armrest Height (top of seat to top of armrest pad)	310 - 460mm
Approvals	AUST R ARTG 291393 Crash tested to AS/NZS3696.19





ASPIRE REHAB RX

MASS QLD ORDER FORM

Effective:

ORDERING INFORMATION: Phone 1300 133 120

ACCOUNT INFORMATION	
Invoice to	Deliver to
Account #	Name
Address	Address
Email	Email
Phone	Phone
Order Date	Marked for
PO#	Order#
Number Required	

EVOKE STANDARD MASS PACKAGE

Frame: Rigid Steel Frame Armrest: Removable, Height & Depth Adjustable Seat: Depth, Width & Height Adjustable

Swing-away, Removable,

Legrest: Padded & Elevating Wheels: Quick Release, Solid 24" Tyres, Padded & Elevating Wheels: Height, Tension & Angle

Articulated Articulated Height & Depth Adjustable Adjustable

Footrest: Height & Angle Adjustable Castors: 6", Height Adjustable Upholstery: PU Breathable Multi-Stretch

gnt & Angle Adjustable Castors: 6, Height Adjustable Upholstery: Fabric Covers

Height, Depth, Angle Attendant Operated Drum Multilayer Castallated Back

Headrest: & Width Adjustable, Brakes: Brakes & Push-to-lock Cushion: and Seat Cushion Covered in

Removable Handbrakes Push-to-lock Cushion: and Seat Cushion Covered in Removable

TOTAL PRICE

FOR OPTIONAL EXTRAS, PLEASE CHOOSE ONE OF THE PACKAGES BELOW:

CHAIR CUSTOMISATION without cushion - PLEASE CIRCLE BELOW

All customisations will be quoted as labour in addition to the base price, customisations can be carried out by Aidacare or by an authorised Aidacare Network Partner.

Seat Width 410 - 460mm 460 - 510mm 510 - 560mm

ACCESSORIES - PLEASE SPECIFY QUANTITY

Standard Included Feature

Code

WCA848160	Elevating Articulating Legrest - Left - 410m	m		
WCA848170	Elevating Articulating Legrest - Right - 410mm			
WCA848180	Elevating Articulating Legrest - Left - 460/510mm			
WCA848190	Elevating Articulating Legrest - Right - 460/510mm			
WCA848500	Rear Wheels w/Drum Brakes Kit - 24in			
WCA848270	Anti-Tip Bars			
Code	Elevating Legrests	Price	Qty	
WCA848100	Left - 410mm			
WCA848110	Right - 410mm			
WCA848120	Left - 460/510mm			
WCA848130	Right - 460/510mm			
Code	Stump Supports	Price	Qty	
WCA848220	Left - Silver - 410/460mm			
WCA848230	Right - Silver - 410/460mm			
WCA848240	Left - Silver - 510mm			
WCA848250	Right - Silver - 510mm			
Code	Chest Supports	Price	Qty	
WCA848450	Left			
WCA848460	Right			

Code	Drum Brakes	Price	Qty
WCA848510	Rear Wheels w/Drum Brakes Kit - 16in		
WCA848520	Rear Wheels w/Drum Brakes Kit - 20in		
Code	Depth Extensions	Price	Qty
WCA848530	Base Frame Extension Kit - 50mm - Black	FIICE	Gity
WCA848540	Base Frame Extension Kit - 100mm - Black		
Code	Other	Price	Qty
WCA848550	Tool-free Adjustment Kit		
WCA848560	Meal Tray - 320mm to 410mm		
WCA848570	Meal Tray - 460mm to 510mm		
WCA848280	IV Pole		
WCA848300	Safety Lap Belt - Metal Seatbelt Buckle		
WAA695600	Oxygen Bottle Bag		
WCA848580	Calf Strap		
WCA848490	Camber Kit		
MASS Plaque No.			
MASS Funded			
Client Contribution			