NSW HEALTH C819 EQUIPMENT CATALOGUE

FEBRUARY 2019

VERSION: 2 UPDATED: MARCH 2019



Quality Equipment Enterprise Solutions Home & Community





Aidacare & NSW Health

Aidacare are proud to be recognised as a

NSW HealthShare / EnableNSW contracted service provider for a large range of products included in the recently announced Beds and Mattresses contract. Our commitment to all NSW HealthShare stakeholders including prescribing health care professionals and most importantly, the NSW HealthShare client, is to provide the very best service and product quality available.

For the healthcare prescriber, Aidacare offer industry leading support to assist in the provision of equipment for your clients.

- All contracted products are available for trial, anywhere in New South Wales, to enable assessment of suitability for your clients specific needs. Our equipment consultants can meet with you at any location, including client's private homes, to demonstrate and trial items.
- Dedicated NSW HealthShare display stock located in each of the Aidacare showrooms is available to enable prescribers to meet with clients on-site to trial and assess all required items in a professional environment.
- Many Aidacare showrooms even offer prescribers use of a private office and meeting room to utilise while meeting with clients or just to catch up on admin and reports.
- Specialist Complex Rehab & Seating staff to assist and support with prescription of customised and specialised equipment.
- Clinical education and support via the Aidacare Professional Training program and workshops.
- Service and repair of all contracted items with our workshop and mobile technicians.

General customer service and support including the booking of trails and any general advice and information required.

Our electronic catalogue is available for download at www.aidacare.com.au/globalassets/resources

For information on any of the extensive range of products available from Aidacare, please browse the complete online Aidacare catalogue at www.aidacare.com.au

For additional information, please contact your local Aidacare office directly on 1300 133 120 or email nswhealth@aidacare.com.au

EnableNSW provides equipment and services to people in NSW with chronic health conditions or disability to assist them with mobility, communication and self-care. We also provide assistance for people who need to travel significant distances for specialist medical care.

Aids and Equipment Program - The Aids and Equipment Program is a NSW Government initiative to assist eligible residents of NSW, who have a life-long or long-term disability, to live and participate within their community by providing appropriate equipment, aids and appliances.

Information on eligibility for the program, the range of equipment available, the lodgement of applications, mutual responsibilities and expectations, and steps for complaints and appeals is outlined in the Policy Directive - "Assistive Technology for Communication, Mobility, Respiratory Function & Self-Care". https://www1.health.nsw.gov.au/ PDS/pages/doc.aspx?dn=PD2011_027

People with a disability interested in seeking further information about the program should discuss their query with a healthcare provider, or Enable NSW on 1800 ENABLE (1800 362 253)

As part of HealthShare NSW, EnableNSW is responsible for the administration of NSW Health disability support and other assistance programs including:

- Aids and Equipment Program
- Continence Assistance
- Home Respiratory Program (HRP)
- Prosthetic Limb Service (formerly known as ALS)
- Specialised Equipment Essential for Discharge (SEED)

For more information, visit the EnableNSW website: http://www.enable.health.nsw.gov.au/

State Wide Distribution Network

AIDACARE

SYDNEY SOUTH 4 Noonan Rd, Ingleburn NSW 2565 Tel: (02) 8706 2300 Fax: (02) 9618 5111

SYDNEY NORTH 2/4 Stanton Rd, Seven Hills NSW 2147 Tel: (02) 8706 2330 Fax: (02) 9624 4256

WOLLONGONG 140 Church St Wollongong NSW 2500 Tel: (02) 8706 2359 Fax: (02) 9618 5111 **ALBURY** 347 Urana Rd, Lavington NSW 2641 Tel: (02) 6040 3444

Fax: (02) 6040 3222 BALLINA 18 Endeavour Close, Ballina NSW 2478

Tel: (02) 6686 6644 Fax: (02) 6686 9383 COFFS HARBOUR 2/177 Orlando St,

2/177 Orlando St, Coffs Harbour NSW 2450 Tel: (02) 6651 5333 Fax: (02) 6652 1272

PORT MACQUARIE

31 Jindalee Rd, Port Macquarie NSW 2444 Tel: (02) 6581 0018 Fax: (02) 6581 0038



TAMWORTH 3/248 Peel St,

3/248 Peel St, Tamworth NSW 2340 Tel: (02) 6761 2018 Fax: (02) 6761 2090



CALL 1300 133 120

Call this number to be automatically diverted to your local Aidacare office (for the price of a local call).

Product Index

Hospital Beds	Pg 04	(7)
Aspire Elite Ward Bed Malvestio Botero Bariatric Ward Bed - 400kg		
Home Care Beds	Pg 06	(j

Aidacare FL250 Floorline Bed	06
Aspire Homecare Bed	07

Mattress Overlays & Replacements

Aspire Lifecomfort Active Air 5 Mattress Overlay	
Aspire Lifecomfort Active Air 8 Mattress Replacement	09
Aspire Lifecomfort Active Air Acute 8 Mattress Replacement	09
Aspire Lifecomfort Active Air Bari 10 Bariatric Mattress Replacement	10
Aspire Lifecomfort Hybrid Acute Mattress Replacement	10
Aspire Lifecomfort Foam Mattress	
Aspire Lifecomfort Deluxe Foam Mattress	
Aspire Lifecomfort Pressure Reducing Foam Mattress	
Aspire Lifecomfort Acute Care Foam Mattress	
Aspire Lifecomfort Bariatric Foam Mattress	
Aspire Lifecomfort Premier Care Foam Mattress	
Aspire Lifecomfort Premier Care Acute Foam Mattress	

Patient Trolleys

Aspire Advanced Electric Emergency Stretcher

Maternity & Birthing

Paediatric Cot - Manual..... Aspire Gynaecology Delivery Bed





Pg 08

Pg 13

Pg 14

Hospital Beds



Aspire Elite Ward Bed



- Designed and Manufactured in Europe with both patient and Caregiver in mind
- Provides clinical benefits for patients and easy and intuitive to use
- Split safety sides with integrated controls as standard
- Angle indicators provide easy to read information for crucial head elevation measurement requirements
- Easy to service and maintain
- Strong, durable and long lasting
- Integrated Bed Extension
- ABS Bed Surfaces



- Electric platform, height and Trendelenburg adjustment
- Locking of functions by nurse control panel
- Mechanical CPR by the handles on both sides of back rest
- Electrical CPR in 1 push
- Fowler and Vascular Positions
- Manual knee mechanism by gradual adjustment
- Easily removable fixed head / foot boards and side rails made from polyoropylene material
- Standard Inclusions: IV Pole, Catheter Hooks & Integrated Controls in Side Panels

:22



Split Safety Sides with Integrated Control -Included

Knee Bend Angle

BEB050130
2000 - 2200mm
1060mm
385 - 775mm
2000mm
860mm
65°
35°
-14° to 14°
150mm
250kg



Included Accessory



IV Pole (BEA826105) **Optional Accessories**





Side Rails

(BEA826110)



(GSP504505)

Hospital Beds

Malvestio Botero Bariatric Ward Bed - 400kg

MALVESTIO



The specific bariatric design focus offers width and length adjustable deck to manage changing needs in both ward and respite environments.

- Integrated tool-free bed width extension of 1020-1220mm
- Integrated bed length extension of 2000-2180mm
- Low height of 380mm
- Antibacterial powder coating for improved infection control management
- Backrest translation (autoregression) reduces compression of the thoracic and abdominal areas as the backrest is raised and alleviates shear and pressure for improved skin integrity
- Compatible with all leading bed movers

Code	BEB050200
Overall Length	2260 - 2440mm
Overall Width	1020 - 1220mm
Overall Height	380 - 850mm
Patient Surface Length	2000 - 2180mm
Patient Surface Width	1020 - 1220mm
Backrest Angle	0° to 70°
Knee Bend Angle	0° to 30°
Trendelenburg	-17° to 17°
Castor Size	150mm
SWL	400kg



Malvestio Botero Bariatric Ward Bed Mattress - 400kg (BEM050200)

- Attendant control panel with intuitive pre-programmed single touch functions for speed and simplicity
- 150mm Tente Castors
- Central locking castors with integrated directional lock
- Full Trendelenburg functionality
- Integrated folding side rails
- · Removable width adjustable head and foot panels
- Standard Inclusions: ACP, Battery Backup, Width Adjustment & Antibacterial Frame





Aidacare FL250 Floorline Bed



- Low Height With it's impressive low height of 110mm, injury to the resident can be significantly reduced or avoided in the event of a fall from bed, therefore reducing the need for side rails
- ExitSafe By pushing the ExitSafe button the FL250 moves to a flat and ergonomic position at the optimum height for resident ingress/exit
- SmartDescent Automatically slows descent speed over the final 300mm of lowering to the floor level, allowing possible impingements to be identified and moved
- BrakeFirst Brakes automatically lock when bed reaches it's lowest height. Brakes are located within the footprint of the bed minimising any trip hazard for he caregiver and resident alike
- AutoRegression Helps to reduce sacral pressure, shear forces and abdominal compression during bed profiling





- Chair Positioning For improved pressure management and respiration. Also a more comfortable dining and sitting position
- Integrated underbed lighting, extendable patient surface and full +/- Trendelenburg functionality

Code	BEB046100 (Single)	BEB046150 (King Single)
Patient Surface	2000mm x 900mm	2000mm x 1060mm
Overall Length	2088mm to 2288mm	
Length Extension	100mm/200mm	
Overall Width	900mm 1060mm	
Patient Surface Height Range	110mm to 740mm	
SWL	250kg	
Backrest Range	0 - 60°	

Knee Bend Range Underbed Clearance **Mattress Platform** Product Weight

(Single)	(King Single)	
2000mm x 900mm	2000mm x 1060mm	
2088mm te	o 2288mm	
100mm/	/200mm	
900mm	1060mm	
110mm to 740mm		
250kg		
0 - 60°		
0 - 35°		
Variable. Minimum 30mm		
4 Sections		
105kg 115kg		

8 xTwin 75mm



Optional Accessories



IV Pole (BEB046530)



Split Folding Side Rails (Pair) (BEB046525)



Self Help Pole (BEB046540)



Castor - Total Locking (BEB046650)



(BEB046520)



Tahoe Bed End Set Sublime Teak Single (BEB046200) King Single (BEB046220)





Mattress **Overlays** & Replacements











- Digital Pump The Easy-to-read LCD screen provides clear and precise information on cycle time, operating mode and fault codes to assist with troubleshooting
- AutoLock Control panel automatically locks after 30 seconds, peventing accidental mode changes
- AutoDetect Based on the user's weight and size, appropriate pressure output is automatically set, ensuring optimal care and protection



- Continuous Low Pressure (CLP) When activated, this mode increases user immersion and envelopment
- 24hr Transport Mode In the event of a power failure, this feature provides a security buffer allowing patients to remain on an inflated surface until power is restored
- Workplace Health & Safety The power cord is retained within the overlay eliminating it's potential as a trip hazard
- Pressure Care Outcomes Intelligent cell design maximises the duration of pressure reduction in the alternation cycle





Code	PTM639200 (Single)	PTM639210 (King Single)
Overall Length	2000)mm
Overall Width	850mm	1050mm
Overall Height	130mm	
Number of Cells	16 (incl 3x static head cells)	
S.W.L	175kg	
Pump	Digital	
Warranty	2 years pump / 1 year mattress	









Rapid Deflation CPR Valve

Digital Pump

Aspire Lifecomfort Active Air 8 Mattress Replacement

- Digital Pump The Easy-to-read LCD screen provides clear and precise information on cycle time, operating mode and fault codes to assist with troubleshooting
- AutoLock Control panel automatically locks after 30 seconds, peventing accidental mode changes
- AutoDetect Based on the user's weight and size, appropriate pressure output is automatically set, ensuring optimal care and protection
- AutoProtect When the head if the bed is raised to 45° or more, air pressure automatically increases to prevent bottoming out
- Continuous Low Pressure (CLP)

 When activated, this mode increases user immersion and envelopment



- 24hr Transport Mode In the event of a power failure, this feature provides a security buffer allowing patients to remain on an inflated surface until power is restored
- Workplace Health & Safety The power cord is retained within the overlay eliminating it's potential as a trip hazard
- Pressure Care Outcomes Intelligent cell design maximises the duration of pressure reduction in the alternation cycle



Code	PTM639300 (Single)	PTM639310 (King Single)
Overall Length	2000mm	
Overall Width	900mm	1050mm
Overall Height	200mm	
Number of Cells	20 (incl 3x static head cells)	
S.W.L	200kg	
Pump	Digital	
Warranty	2 years pump / 1 year mattress	









Rapid Deflation CPR Valve

Digital Pump

Aspire Lifecomfort Active Air Acute 8 Mattress Replacement

 Dual Compressor Digital Pump

 The Easy-to-read LCD screen provides clear and precise information on cycle time, operating mode and fault codes to assist with troubleshooting

- One in three alterntaing air cycle provides acute level therapy for maximum pressure redistribution and wound management
- AutoDetect Based on the user's weight and size, appropriate pressure output is automatically set, ensuring optimal care and protection
- AutoProtect When the head if the bed is raised to 45° or more, air pressure automatically increases to prevent bottoming out
- Continuous Low Pressure (CLP)

 When activated, this mode increases user immersion and envelopment





- 24hr Transport Mode In the event of a power failure, this feature provides a security buffer allowing patients to remain on an inflated surface until power is restored
- Workplace Health & Safety The power cord is retained within the overlay eliminating it's potential as a trip hazard
- Pressure Care Outcomes Intelligent cell design maximises the duration of pressure reduction in the alternation cycle



Code	PTM639400 (Single)	PTM639410 (King Single)
Overall Length	2000mm	
Overall Width	900mm	1050mm
Overall Height	200mm	
Number of Cells	24 (incl 3x static head cells)	
S.W.L	250kg	
Pump	Digital	
Warranty	2 years pump / 1 year mattress	













Rapid Deflation **CPR** Valve

Digital Pump

..........

Aspire Lifecomfort Active Air Bari 10 Bariatric Mattress Replacement - 450kg

- Quick Inflation Blower Pump High powered mattress pump ensures rapid inflation of mattress for quick use and ensures constant high powered air flow which is extremely important for bariatric care.
- Turning Cycle With the ability to turn a patient, the BARI 10 ensures for safe patient transfer and lateral movement, improving manual handling and wound therapy outcomes.
- AutoDetect Based on the user's weight and size, appropriate pressure output is automatically set, ensuring optimal care and protection
- AutoProtect When the head if the bed is raised to 45° or more, air pressure automatically increases to prevent bottoming out

Continuous Low Pressure (CLP) -When activated, this mode increases user immersion and envelopment

2) aspire

- 1 in 2, 3 or 4 alternating air cycles provides acute level therapy for maximum pressure redistribution and wound management.
- Workplace Health & Safety The power cord is retained within the overlay eliminating it's potential as a trip hazard
- Pressure Care Outcomes Intelligent cell design maximises the duration of pressure reduction in the alternation cycle



Code	PTM639500 (Single)	PTM639510 (King Single)	PTM639520 (Super King Single)
Overall Length	2000mm		
Overall Width	900mm 1050mm 1230mm		1230mm
Overall Height	270mm		
Number of Cells	20 (incl 3x static head cells)		
S.W.L	450kg		
Pump	Digital		
Warranty	2 years pump / 1 year mattress		





Additional Digitial Pump (PTM638030) for an active support surface

Aspire Lifecomfort Hybrid Acute Mattress Replacement

- Self Adjusting Pressure Redistribution -• In non-powered mode, nine interconnected foam filled cells allow air to move freely between them via the user's body movements
- Active Support Surface Solution -By connecting the pump, the mattress converts to an active support surface
- Active Pressure Redistribution In powered function, the Alternating and Static modes become available
- Memory Foam Top Layer Integrated memory foam overlay improves immersion, envelopment and comfort
- Intelligent Hose Management When not in use the air-hose can be stored internally, reducing infection control and OH&S issues
- Integrated Shear Reduction The memory foam overlay moves independently of the underlying cells, reducing shear as the bed profiles



- Optimised Heel Protection A 12° tapered heel zone redistributes pressure up the lower leg, away from the vulnerable heel region. The castellated surface collapses in a wave-like movement reducing shear
- CPR Compatible In non-powered mode, the foam-core cells allow CPR to take place. When in powered mode, simply disconnect the pump to facilitate resuscitation



















Aspire Lifecomfort Foam Mattress

- Pressure care mattress
- Two layer foam support design for comfort and support
- Full size memory foam overlay
- Two way waterproof stretch cover
- Vinyl bottom cover
- Three sided zip for easy cover removal

Aspire Lifecomfort

Pressure care mattress

comfort and support

Vinyl bottom cover

Deluxe Foam Mattress

• Three layer foam support design for

 Full size memory foam overlay Two way waterproof stretch cover

Three sided zip for easy cover

 Stitched zipper cover with waterproof flap for excellent hygiene

	_	ZUKG	
Code	BEM046410 (Single)	BEM046460 (King Single)	
Overall Length	2000mm	2030mm	
Overall Width	900mm	1050mm	
Overall Height	130mm		
S.W.L	120kg		
Warrantv	5 Years Foam, 2 Years Cover		

🗊 aspire comfort



- removal Stitched zipper cover with waterproof flap for excellent hygiene
- PUPPS 2 (2004) compliant

Aspire Lifecomfort Pressure Reducing Foam Mattress

- Pressure care mattress
- Four layer foam support design for maximum comfort and support
- Full size memory foam overlay
- Two way waterproof stretch cover
- Vinyl bottom cover
- Three sided zip for easy cover removal
- Stitched zipper cover with waterproof flap for excellent hygiene
- PUPPS 2 (2004) compliant

Aspire Lifecomfort Acute Care Foam Mattress

- Strong U-Shaped Foam Outer Structure Promotes Ease of Patient Transfers
- Welded High Temp Washable Fabric for Optimal Infection Control
- Vapour Permeable Four Way Stretch Fabric Improves Skin Microclimate
- Castellated Patient Surface for Improved Air Flow & Efficient Pressure Relief

Code	BEM046615 (Single)	BEM046616 (King Single)
Overall Length	1980mm	2030mm
Overall Width	880mm	1050mm
Overall Height	150mm	150mm
S.W.L	250kg	
Warranty	5 Years Foam, 2 Years Cover	



11





Overall Height	150mm
S.W.L	150kg
Warranty	5 Years Foam, 2 Years Cover
	-

Code

Code

Warranty







5 Years Foam, 2 Years Cover









Aspire Lifecomfort Bariatric Foam Mattress



Four layer foam, designed for larger users with a maximum user weight of 300kg. Full size memory foam overlay and side walls for added stability and ease of transfers. European fabric cover with fully welded seams, water resistant zippers and reinforced waterfall flaps for optimum infection control.

- Pressure reducing bariatric mattress
- Extra weight capacity up to 300 kg
- Four layer foam support 200 mm high for bariatric use with maximum comfort and support
- Full size memory foam overlay
- Foam side walls for extra stability on the edge of the mattress
- Two way waterproof stretch cover with welded seams
- PVC bottom cover
- Four sided zip for easy cover removal

- Stitched zipper cover with waterproof flap for excellent hygiene
- Independently tested and certified for durability and weight rating



27 aspire

ecomfort

BEM046700 (King Single)
2030mm
1050mm
200mm
300kg
5 Years Foam, 2 Years Cover





Aspire Lifecomfort Premier Care Foam Mattress

- Pressure care mattress suitable for high risk patients
- Four layer foam support for maximum comfort and support
- Full size memory foam top layer
- Integrated cotton stockinette promotes foam longevity
- Reinforced side walls for patient stability and ease of transfers
- Inbuilt handles for safe manual handling
- Two way waterproof stretch top cover
- Heavy duty PVC base cover
- Three sided zipper for easy cover removal



• PUPPS 2 (2004) Compliant

	s.w.L 180kg	
Code	BEM046725 (Single)	BEM046730 (King Single)
Overall Length	2000mm	2030mm
Overall Width	900mm	1050mm
Overall Height	150mm	150mm
S.W.L	180kg	
Warranty	5 Years Foam, 2 Years Cover	





Aspire Lifecomfort Premier Care Acute Foam Mattress

- Pressure care mattress suitable for high risk patients
- Four layer foam support for maximum comfort and support
- Full size memory foam top layer with MEMair[™] Breeze open cell foam for improved microclimate
- Integrated cotton stockinette promotes foam longevity
- Reinforced side walls for patient stability and ease of transfers
- Inbuilt handles for safe manual handling
- Two way waterproof stretch top cover
- Heavy duty PVC base cover

- Three sided zipper for easy cover
- Stitched zipper cover with waterproof flap for improved infection control
- PUPPS 2 (2004) Compliant

removal



🧊 aspire

Code	BEM046735 (Single)	BEM046740 (King Single)
Overall Length	2000mm	2030mm
Overall Width	900mm	1050mm
Overall Height	150mm	150mm
S.W.L	250kg	
Warranty	5 Years Foam, 2 Years Cover	



Aspire Advanced Electric Emergency Stretcher

🦪 aspire







- Electric operated lift, tilt, back raise, Trendelenburg and reverse Trendelenburg functions
- Heavy Duty 225kg Safe Working Load
- Full length sliding X-ray deck ensures complete spinal coverage
- Comfortable 100mm foam mattress with multi-stretch cover
- Safety straps for staff and patient security
- CPR functions integrated on either side for rapid response time
- Collapsible vertical side rails fold away below the patient surface
- Fold-away carer push handles to maximise workspace

- 200mm castors minimise rolling resistance
- Fifth wheel improves steering and reduces manual handling burden in transportation
- Central and Directional Locking functions
- Molded base cover facilitates cleaning and supports infection control processes
- Large corner bumpers minimise wear and tear on equipment and facilities
- Intuitive design and function layout from Red Dot design winning team
- Standard Inclusions: Mattress, IV Pole, O2 Holder, Uring Bag Hooks & Accessory Slots



13

Code	TSE674230
Overall Length	2140mm
Overall Width	840mm
Overall Height	600 - 910mm
Patient Surface Length	1920mm
Patient Surface Width	600mm
Backrest Angle	0° to 70°
Trendelenburg	-12° to 12°
Castor Size	200mm
SWL	225kg

Maternity & Birthing

Paediatric Cot - Manual





Code	WAT830550
Overall Length	1430mm
Overall Width	875mm
Overall Height	1560mm
Patient Surface Length	1215-700mm
Patient Surface Width	700mm
Trendelenburg	12°
Castor Size	150mm
SWL	40kg

- Stainless steel & powder coated aluminium construction
- Soft internal edge design
- Smooth glide side and end panels
- 12° Trendelenburg
- Lightweight aluminium side and end panels
- Buffers
- Hand or knee operated panel latch
- Obstruction free base design for access on all sides and ends
- Castor Options:
 - 100mm castors
 - 125mm castor
 - 125mm turn wheel castors
 - Central locking





Aspire Gynaecology Delivery Bed

🦪 aspire



Designed to deliver safe, fast, effective care in the busiest of ward environments

- Retractable leg section, easily and quickly slides under the bed and out of the way for delivery.
- Single action gas assisted side rails minimise transfer gaps for patient transfer
- Integrated gynae cut-out promotes accessibility and effective positioning
- · Ergonomically designed patient handles
- Removable padded leg supports maximise patient postural support and can be safely stored when not in use.
- Electrically adjustable for patient comfort and clinically correct positioning



BACK RAISE

HEIGHT TRENDELENBURG ADJUSTABLE

Code	BEB050310
Overall Length	2175mm
Overall Width	1060mm
Overall Height	620 - 890mm
Patient Surface Length	2160mm
Patient Surface Width	955mm
Backrest Angle	0° to 63°
Trendelenburg	-11º to 11º
Castor Size	125mm
S.W.L	175kg

Optional Accessories



Aspire Delivery Bed Mattress (BEM050310)



· Comfortable 2-section mattress includes multi-stretch cover for

• Intuitive Linak hand control and electronics for reliability and simplicity

• Integrated Trendelenburg capabilities support clinical outcomes and

• Tente 125mm Central locking castors with directional functionality

• Standard Inclusions: Mattress, IV Pole, Wall Buffers, Urine Bag Holder,

patient comfort and improved infection control

Stainless steel bowl with easy access positioning

of servicing

facilitate manual handling

Belt Hooks & Accessory Slots

Placenta Bowl (GSP501220)





#