

## **OSD ASPIRE EVOKE**

## **SWEP SCRIPT FORM**

Effective:

ORDERING INFORMATION: Phone 1300 133 120

DEALER INFORMATION			
Date Account Code Sales Rep Mobile Signature	Name Address City / State Postcode Fax Email		
CUSTOMER INFORMATION			
Customer Name  MASS Plaque #  MASS PO#  Signature  Rep Name & Ph #	Conact Address City / State Postcode Phone Fax Email		
PRESCRIBER INFORMATION			
Prescriber Name	Organisation  Address  City / State  Postcode		
Signature	Phone Fax Email		

#### PRODUCT INFORMATION

The OSD Aspire EVOKE is a highly adjustable lightweight aluminum wheelchair designed to provide a personalised fit and feel for every user.

Extensive seat and wheel adjustability offers maximum independence for the user without compromising clinical outcomes.

With a comprehensive range of accessories available, the EVOKE provides maximum performance and is suitable for indoor and outdoor use.

Max. User Weight	S.W.L 130kg				
Weight	15kg (chair weight) 8.5kg (transport weight)				
Seat Height - adjustable (floor to top of fabric)	FRONT: 50cm 47cm 45cm REAR: 45cm 43cm 41cm				
Seat Width	35cm (14") 40cm (16") 45cm (18") 50cm (20")				
Seat Depth - adjustable	40cm (16") 43cm (17") 45cm (18") 48cm (19") 50cm (20") 53cm (21")				
Backrest Height - adjustable (top of seat to top of backrest handle)	40cm - 45cm				
Overall Width	Seat width plus 20cm				
Armrest Height (top of seat to top of armrest pad)	22cm fixed type, variable height and depth on adjustable option				
	TGA# ARTG224574				
Approvals	AS/NZS tested and approved				











# **OSD ASPIRE EVOKE**

### **SWEP SCRIPT FORM**

Effective:

ORDERING INFOR	RMATION: Phone 13	00 133 120				
ACCOUNT INFO	RMATION					
Invoice to Account # Address			Nan	ver to ne dress		
Email Phone Order Date PO # Number Required						
EVOKE STANDA	ARD MASS PACKA	GE				
Colour Silv	/er	Seat	Standard pa	dded nylon	Footplate	Angle adjustable footplat
Backrest Ter Rear Wheels 24'	nsion adjustable ,	Back Upholstery  Armrest	Tension adju	istable back	Castor Brake	8" Push to lock
		Legrest	70° swing in	/ out	Accessories	Tool kit
FOR OPTIONAL E	XTRAS, PLEASE C	HOOSE ONE OF THE P.	ACKAGES BE	LOW:		
CHAIR MEASURE	EMENTS without c	ushion - PLEASE SELE	СТ			
	s will be quoted as I are Network Partne	abour in addition to the	e base price, c	ustomisations can	be carried out by	y Aidacare or by an
Seat Width	35 cm	40 cm	45 cm	50 cm		
Seat Depth	원 8 40 cm	43 cm	45 cm	¥ 48 ст	50 cm	53 cm
Seat Height	50 cm	47 cm	45 cm	45 cm	43 cm	41 cm
ACCESSORIES -	PLEASE SPECIFY	QUANTITY				
Code	Accessory				Pri	ce Qty
WCA858370	Elevating Legre	est - Black Left			\$104	
WCA858390	Elevating Legrest - Black Right \$104		)4			
WCA858440	Stump Suppor	t - Black Left			\$141	
WCA858450	Stump Support	t - Black Right			\$141	
MSP421020	Armrests (Height & Depth Adj.) - Pair \$99			9		
WCA858460	Anti Tip Bars \$52			2		
WCA858463	Transport Tie D	Downs			\$52	2
MWT469865	24" Pneumatic	Rear Wheels (Pair)			\$11.	5
					NOTES	
Client Contribution TOTAL PRICE	າ					
TOTAL FRICE						