MASS EQUIPMENT CATALOGUE

FEB 2020 VERSION: 2 FEBRUARY 2020



Quality Equipment Enterprise Solutions Home & Community



Introduction to MASS



Aidacare & MASS

Aidacare, in conjunction with our supply partners are proud to be recognised suppliers of a large number of products listed on the current MASS SOA (Standing Offer Arrangement).

The SOA is a list of preferred products which MASS is offered for purchase at contracted prices. Each product has a MASS-specific order form, which must be completed and submitted with the application.

Aidacare offer a broad range of services to MASS prescribers and clients including:

- · Access to an extensive range of SOA and other demonstration stock for trial throughout Queensland
- Specialist sales staff to assist with trials and support prescription of customised and specialised equipment
- Clinical specialists to provide additional advice
 and support for healthcare professionals
- Clinical education and support
- Service & repair field technicians
- · General customer service and support including the booking of equipment trials

Our electronic catalogue is available for download at www.aidacare.com.au/contracts/mass.aspx.

For information on non-SOA products please browse the complete online Aidacare Catalogue at www.aidacare.com.au.

For additional information please contact the Queensland Customer Service Team on 1300 133 120 or email mass@aidacare.com.au.



Contents



Bathroom & Toilet		05
Aspire Shower Commodes05 Verto Shower Commodes	Aspire Shower Commodes Accessories	09
Bedroom		10
Bedside Commodes10 Lifecomfort Foam Mattresses11	ROHO Mattress Section Lifecomfort Active Air Mattresses	
Cushions		16
ROHO Pressure Care Cushions16	Promedicare Pressure Care Cushions	
Manual Handling - Lifters & Slings		19
Lifting Devices	Orbi-Turn Patient Turner Aspire Slings Oxford Slings	24
Walking Frames & Aids - Mobility		27 ()
Seat Walkers / Rollators27		
Wheelchairs & Seating		30
Aspire Manual Wheelchairs	Aspire Wheelchair Accessories	
Resources		33








MASS Order & Approval Process

Aidacare offers a service to Prescribers that is based on "Partnering to provide the best solution". We listen, respond and work with therapists to provide the highest possible customer service solution for the client.



Shower Commodes



Aspire Bariatric Shower Commode

- Super Heavy Duty Bariatric design with 790mm seat and 350kg SWL
- Anti-bacterial frame coating for improved infection control
- Removable arms to assist with positioning and allow lateral slide transfers
- Height adjustable, swing away footrests for optimal positioning and ease of transfers
- Padded seat and arms for maximum comfort



- 4x lockable castors maximise safety by providing carer access to a brake at all angles around the commode
- Comes standard with pan and carrier

Code	Description	Overall Height	Overall Width	Seat Width	Backrest Height	Armrest Height	Anterior Seat Height	Posterior Seat Height	SWL
BTC061200	Bariatric Shower Commode	990 mm	900 mm	790 mm	430 mm	170 mm	570 mm	570 mm	350 kg



Aspire 460mm Swing Away Footrest Commode



- 8° backrest angle ensures patient comfort
- Improved over toilet clearance
- Integrated access trough improves hygiene whilst maintaining toilet bowl clearance height of 465mm
- 4x total locking stainless steel castors provide brake access from all angles
- Safety arms prevent falls and swing back out of the way for transfers
- Open armrests style accommodates a wider sitting surface whilst maintaining a smaller base footprint for improved accessibility
- Angle adjustable backrest from 8° to 18° for improved seating comfort
- Ultra low seat height of 590mm which can be reduced to 490mm when customised
- Pan carrier and open seat included



Self propelled model available

Code	Description	Overall Height	Overall Width	Seat Width	Backrest Height	Armrest Height	Anterior Seat Height	Posterior Seat Height	S.W.L
BTC065900	Swing Away - 460mm - Att Prop - S/Steel	895 - 995 mm	535 mm	460 mm	305 mm	225 mm	590 mm	560 mm	200 kg
BTC066100	Swing Away - 460mm - Self Prop - S/Steel	895 - 995 mm	665 mm	460 mm	305 mm	225 mm	590 mm	560 mm	200 kg





Aspire 460mm Weight Bearing Commode

- 8° backrest angle ensures patient comfort
- Improved over toilet clearance
- Integrated access trough improves hygiene whilst maintaining toilet bowl clearance height of 465mm
- 4x total locking stainless steel castors
- provide brake access from all anglesSafety arms prevent falls and swing back out of the way for transfers
- Open armrests style accommodates a wider sitting surface whilst maintaining a smaller base footprint for improved accessibility
- Easy-slide weight bearing platform footrest slides completely away for transfers
- Pan carrier and open seat included





Self propelled model available

Code	Description	Overall Height	Overall Width	Seat Width	Backrest Height	Armrest Height	Anterior Seat Height	Posterior Seat Height	S.W.L
BTC066000	WB Platform - 460mm - Att Prop - S/Steel	995 mm	535 mm	460 mm	305 mm	225 mm	590 mm	560 mm	200 kg
BTC066200	WB Platform - 460mm - Self Prop - S/Steel	995 mm	665 mm	460 mm	305 mm	225 mm	590 mm	560 mm	200 kg



Aspire 530mm Swing Away Footrest Commode

- 8° backrest angle ensures patient comfort
- Improved over toilet clearance
- Integrated access trough improves hygiene whilst maintaining toilet bowl clearance height of 465mm
- 4x total locking stainless steel castors provide brake access from all angles
- Safety arms prevent falls and swing back out of the way for transfers
- Open armrests style accommodates a wider sitting surface whilst maintaining a smaller base footprint for improved accessibility
- Angle adjustable backrest from 8° to 18° for improved seating comfort
- Ultra low seat height of 590mm which can be reduced to 490mm when customised
- Pan carrier and open seat included



🦪 aspire

Self propelled model available

Code	Description	Overall Height	Overall Width	Seat Width	Backrest Height	Armrest Height	Anterior Seat Height	Posterior Seat Height	S.W.L
BTC066065	Swing Away - 530mm - Self Propelled - S/Steel	995mm	735mm	530mm	305mm	225mm	590mm	560mm	200kg



Aspire 530mm Weight Bearing Commode

- 8° backrest angle ensures patient comfort
- Integrated access trough improves hygiene whilst maintaining toilet bowl clearance height of 465mm
- Safety arms prevent falls and swing back out of the way for transfers
- Easy-slide weight bearing platform footrest slides completely away for transfers
- 4x total locking stainless steel castors provide brake access from all angles
- Open armrests style accommodates a wider sitting surface whilst maintaining a smaller base footprint for improved accessibility
- Improved over toilet clearance
- Pan carrier and open seat included





Self propelled model available

Code	Description	Overall Height	Overall Width	Seat Width	Backrest Height	Armrest Height	Anterior Seat Height	Posterior Seat Height	S.W.L
BTC066050	WB Platform - 530mm - Att Propelled - S/Steel	995mm	605mm	530mm	305mm	225mm	590mm	560mm	200kg
BTC066055	WB Platform - 530mm - Self Propelled - S/Steel	995mm	735mm	530mm	305mm	225mm	590mm	560mm	200kg



Aspire 600mm Swing Away Footrest Commode

- 8° backrest angle ensures patient comfort
- Safety arms prevent falls and swing back out of the way for transfers
- Open armrests style accommodates a wider sitting surface whilst maintaining a smaller base footprint for improved accessibility
- 4x total locking stainless steel castors provide brake access from all angles
- Self Propelled Model Available
- Ultra low seat height of 590mm which can be reduced to 490mm when customised

- Angle adjustable backrest from 8° to 18° for improved seating comfort
- Integrated access trough improves hygiene whilst maintaining toilet bowl clearance height of 465mm
- Pan carrier and open seat included

Code	Description	Overall Height	Overall Width	Seat Width			Anterior Seat Height	Posterior Seat Height	S.W.L
BTC061300	Swing Away - 600mm - Att Propelled - S/Steel	995mm	675mm	600mm	305mm	225mm	590mm	560mm	200kg
BTC061305	Swing Away - 600mm - Self Propelled - S/Steel	995mm	675mm	600mm	305mm	225mm	590mm	560mm	200kg



Aspire 600mm Weight Bearing Commode - 400kg

- Heavy duty reinforced wide frame can support bariatric users up to 400kg
- 4x total locking stainless steel castors provide brake access from all angles
- Height adjustable reinforced armrests offer drop flat positioning to assist with lateral slide transfers

🦪 aspire

🧊 aspire

- Easy-slide weight bearing platform footrest slides completely away for transfers
- Pan carrier and open seat included

Code	Description	Overall Height	Overall Width		Backrest Height		Anterior Seat Height	Posterior Seat Height	S.W.L
BTC061400	WB Platform - 600mm - Att Propelled - S/Steel - 400kg	995mm	725mm	600mm	300mm	290mm	590mm	560mm	400kg



Aspire Tilt & Recline Shower Commode

• 30° of backrest angle adjustment in 10°

maximum user comfort

brake access from all angles

showering process

bathing needs

•

increments for effortless repositioning and

• -5° to 35° of gas-assisted tilt-in-space functionality

minimises manual handling efforts during the

• 4x total locking stainless steel castors provide

Open back straps provide greater access for

Backrest may be configured as full length or 3/4

🦪 aspire

length to accommodate an advanced headrest

- Open armrests style accommodates a wider sitting surface whilst maintaining a smaller base footprint for improved accessibility
- Height adjustable, swing away footrests with calf pads for optimal positioning, support and ease of transfers
- Tool-free height adjustment to 4 settings allow the commode to be easily adapted to varying needs
- Seat Backrest Armrest Posterior Overall S.W.L Overall Height Code Description Seat Height Width Width Height Seat Heigh Height BTC066030 460mm - Attendant Propelled 1315 - 1485mm 460mm 560mm 900mm 275mm 460 - 630mm 460 - 620mm 250ka BTC066040 530mm - Attendant Propelled 1315 - 1485mm 630mm 530mm 900mm 275mm 460 - 630mm 460 - 620mm 250kg





Verto 475mm Height Adjustable Shower Commode - Attendant Propelled

- Adjustable seat height in increments of 25mm
- Depth and height adjustable polyurethane backrest which is removable for cleaning purposes
- Height adjustable swing-away legrests
- Robust, lightweight stainless steel construction
- Height adjustable armrests
- 5 year warranty
- *Seat Sold Separately



Self Propelled Model

Code	Description	Seat Width	Seat Height	Overall Width	Backrest Height	S.W.L
BTC069490	Swing Away - 475mm - Att Prop	475 mm	495 - 620 mm	545 mm	485 mm	180 kg
BTC069480	Swing Away - 475mm - Self Prop	475 mm	495 - 620 mm	660 mm	485 mm	180 kg



Verto Height Adjustable Shower Commode - Attendant Propelled

- Adjustable seat height in increments of 25mm
- Depth and height adjustable polyurethane backrest which is removable for cleaning purposes
- Weight bearing footplate
- Robust, lightweight stainless steel construction
- Height adjustable swing-away legrests
- Height adjustable armrests
- 5 year warranty
- *Seat Sold Separately

Code	Description	Seat Width	Seat Height	Overall Width	Backrest Height	S.W.L
BTC069495	WB Platform - 475 mm - Att Prop	475 mm	500 - 625 mm	550 mm	485 mm	180 kg





Aspire Shower Commode Accessories



Code	Description
BTC066500	Shower Commode Seat - Aspire - Closed 410mm
BTC066700	Shower Commode Seat - Aspire - Closed 460mm
BTC066600	Shower Commode Seat - Aspire - Open 460mm
BTC066750	Shower Commode Seat - Aspire - Open 530mm
BTC066800	Shower Commode Seat - Aspire - Open 600mm
BTC089124	Shower Commode Seat - Aspire - Hand Made Pressure Reducing - 320mm - Closed Front
BTC089128	Shower Commode Seat - Aspire - Hand Made Pressure Reducing - 400mm - Closed Front
BTC089140	Shower Commode Seat - Aspire - Hand Made Pressure Reducing - 460mm - Closed Front
BTC089180	Shower Commode Seat - Aspire - Hand Made Pressure Reducing - 500mm - Closed Front
BTC089200	Shower Commode Seat - Aspire - Hand Made Pressure Reducing - 550mm - Closed Front
BTC089240	Shower Commode Seat - Aspire - Hand Made Pressure Reducing - 650mm - Closed Front
BTC089220	Shower Commode Seat - Aspire - Hand Made Pressure Reducing - 600mm - Closed Front
BTC089122	Shower Commode Seat - Aspire - Hand Made Pressure Reducing - 320mm - Open Front
BTC089126	Shower Commode Seat - Aspire - Hand Made Pressure Reducing - 400mm - Open Front
BTC089130	Shower Commode Seat - Aspire - Hand Made Pressure Reducing - 460mm - Open Front
BTC089170	Shower Commode Seat - Aspire - Hand Made Pressure Reducing - 500mm - Open Front
BTC089190	Shower Commode Seat - Aspire - Hand Made Pressure Reducing - 550mm - Open Front
BTC089210	Shower Commode Seat - Aspire - Hand Made Pressure Reducing - 600mm - Open Front
BTC089230	Shower Commode Seat - Aspire - Hand Made Pressure Reducing - 650mm - Open Front
BTC089150	Shower Commode Seat - Aspire - Hand Made Pressure Reducing - 460mm - Open Left
BTC089160	Shower Commode Seat - Aspire - Hand Made Pressure Reducing - 460mm - Open Right
BTC067000	Shower Commode Armrest - Aspire - Padded Offset
BTC066400	Shower Commode - Bowl & Pan Carrier - Aspire (410mm wide) - Fits 460/530/600mm
BTC067110	Shower Commode Elevating Legrest - Aspire - Left
BTC067115	Shower Commode Elevating Legrest – Aspire - Right
BTC089100	Hand Made-Padded Seat Cde Closed front sizes 1- 4 MSC
BTC089110	Hand Made Padded Shower Cmmd Open front MSC
BTC089115	Hand Made Padded Shwr Bariatric Cde St Opn/cls MSC
BTC089120	Hand Made Prssur Rliving Shwr Cmmd St Opn/clsd frn



Bedside Commodes



Raised armrests support the user to sit down and stand up steadily








Lifecomfort Standard Foam Mattress

- Two layer foam support design for comfort and support
- Full size memory foam overlay
- Two way waterproof stretch cover





Serial Tracked



- Three sided zip for easy cover removal
- Welded zipper cover with waterproof flap for the best hygiene



Infection Control

Vinyl bottom cover

Cover

Code	Description	Height	Width	Length	S.W.L
BEM046410	Standard Foam Mattress - AC Series Single	130 mm	900 mm	2000 mm	120 kg
BEM046460	Standard Foam Mattress - King Single	130 mm	1050 mm	2030 mm	120 kg



Lifecomfort Foam Mattress with Safety Sides

- Sidewall support to help minimise the risk of falls from bed
- Two layer foam support design for comfort and support
- Full size memory foam overlay
- Two way waterproof stretch cover





- Vinyl bottom cover
- Three sided zip for easy cover removal
- Welded zipper cover with waterproof flap for the best hygiene

		Foam Levels	(High Quality Standard Foam)	SerialTracked	I Infectio	n Control	Foam	Cover
Code	Description			Height	Width	Length	Safety Sides	S.W.I
BEM046475	Foam Mattress with Safety Sides			130 mm	900 mm	1940 mm	220 mm	120 k
-	Medium Risk High Risk Very High	ıh Risk		130 mm	900 mm	1940 mm	220 mm	120




kg

Foam Mattresses Bedroom



Lifecomfort Deluxe Foam Mattress

- Three layer foam design for enhanced comfort and support
- Full size memory foam overlay
- Two way waterproof stretch cover



50 mm (Hospital Grade) Serial Tracked Foam Levels





Three sided zip for easy cover removal

Welded zipper cover with waterproof flap

Infection Control

Vinyl bottom cover

for the best hygiene

D	
Year	Year
/arranty	Warranty
Foam	Cover

27 asoire

lifecomfort

Code	Description	Height	Width	Length	S.W.L
BEM046500	Deluxe Foam Mattress - Short	150 mm	900 mm	1450 mm	150 kg
BEM046525	Deluxe Foam Mattress - Small	150 mm	820 mm	1920 mm	150 kg
BEM046510	Deluxe Foam Mattress - AC Series Single	150 mm	900 mm	2000 mm	150 kg
BEM046550	Deluxe Foam Mattress - Long	150 mm	850 mm	2000 mm	150 kg
BEM046455	Deluxe Foam Mattress - King Single	150 mm	1050 mm	2030 mm	150 kg



Lifecomfort Pressure Reducing **Foam Mattress**

- Pressure reducing mattress
- Four layer foam support design for maximum comfort and support
- Full size memory foam overlay
- Foam side walls for extra stability on the edge of the mattress





30 mm (Superior Quality Foam) Foam Levels

Serial Tracked



- Two way waterproof stretch cover with welded seams
- PVC bottom cover
- Three sided zip for easy cover removal
- Welded zipper cover with waterproof flap for the best hygiene





Code	Description	Height	Width	Length	S.W.L
BEM046625	Pressure Reducing Foam Mattress - Small	150 mm	820 mm	1920 mm	150 kg
BEM046600	Pressure Reducing Foam Mattress - Regular	150 mm	900 mm	1940 mm	150 kg
BEM046610	Pressure Reducing Foam Mattress - AC Series Single	150 mm	900 mm	2000 mm	150 kg
BEM046650	Pressure Reducing Foam Mattress - Long	150 mm	850 mm	2000 mm	150 kg
BEM046660	Pressure Reducing Foam Mattress - King Single	150 mm	1050 mm	2030 mm	150 kg



Lifecomfort Acute Care Foam Mattress

- Acute pressure care mattress
- Castellated foam for effective pressure redistriubution
- Enhanced microclimate with air channels between castellations
- Durable base facilitates patient transfers on and off the mattress
- Stitched zipper cover with waterproof flap for excellent hygiene



- High-stretch, vapour permeable, waterproof cover
- Heavy duty anthracite bottom cover for durability
- Innovative cover graphics including cleaning instructions, product specifications and 180° turning process



Serial Tracked

BEM046615 Acute Care Mattress - Single 150 mm 880 mm 1980 mm 2	Code	Description	Height	Width	Length	S.W.L
	BEM046615	Acute Care Mattress - Single	150 mm	880 mm	1980 mm	250 kg
BEM046616 Acute Care Mattress - King Single 150 mm 2030 mm 2	BEM046616	Acute Care Mattress - King Single	150 mm	1050 mm	2030 mm	250 kg

High Risk

V

Very High Risk

Low Risk (M) Medium Risk (H)

Foam Mattresses



Lifecomfort Bariatric Foam Mattress

- Pressure reducing bariatric mattress
- Extra weight capacity up to 300kg
- Four layer foam support 200mm high for bariatric use
- with maximum comfort and supportFull size memory foam overlay
- Foam side walls for extra stability on the edge of the mattress

25 mm (Memory Foam)



- Two way waterproof stretch cover with welded seams
- PVC bottom cover
- Four sided zip for easy cover removal
- Welded zipper cover with waterproof flap for the best hygiene
- Independently tested and certified for durability and weight rating

Code	Description	Height	Width	Length	Safety Sides	S.W.L
BEM046700	Bariatric Foam Mattress	200 mm	1050 mm	2030 mm	300 kg	300kg



Note: pictured both as the standard unit as well with the optional ROHO mattress section in place

📘 Low Risk 🚺 Medium Risk 🚹 High Risk

Lifecomfort Pressure Reducing with ROHO Insert Cut Out

- Bed mattress with ultimate pressure reduction solution
- Uses standard size ROHO mattress section (PTM636320) for additional pressure relief in the pelvic section of the mattress
- Support, comfort and pressure reduction all in one
- Four layer foam support design with side walls for extra stability
- Full size memory foam overlay
- Two way welded waterproof stretch cover



- Three sided zip for easy cover removal
- Welded zipper cover with waterproof flap for the best hygiene





Cut Out

ROHO Section (Sold Separately)

Code	Description	Height	Width	Length	S.W.L	Warranty
BEM046720	Pressure Reducing with ROHO® Insert Cut Out	150 mm	900 mm	1940 mm	150 kg	5 Years Foam / 2 Years Cover



ROHO® Mattress Sections

- Dry Flotation® mattress system allows user to lie or sit comfortably for long periods
- Interconnected and flexible air cells provide protection and relief for user at risk of skin breakdown or already with early stage ulcers
- Adjustable air pressure to suit the user
- · Flexibility provides effective friction and shear relief
- Drainage holes allow fluids to escape

Note: Bedsize divided into 4 sections where ROHO Mattress Sections can be applied. 1, 2, 3 or 4 sections can be used in combination with foam sections where risk of pressure is not as critical

Code	Description	Height	Width	Length	S.W.L	Warranty
PTM636335	DRY FLOATATION® Mattress Overlay (4 sections)	80 mm	860 mm	1960 mm	Unlimited	2 Years
PTM636320	DRY FLOATATION® Mattress Section (1 section)	80 mm	860 mm	490 mm	Unlimited	2 Years
PTM636360	Levelling Pad - Long	80 mm	864 mm	2032 mm	100 kg	2 Years
PTM636220	Levelling Pad - Short	80 mm	863 mm	508 mm	100 kg	2 Years

Very High Risk









The Lifecomfort Hybrid Mattress offers foam-in-cell technology - enabling flexible operation as either a high comfort foam static support surface or powered-up alternating pressure support surface.

- Intelligent Hose Management when not in use the air-hose can be stored internally, reducing infection control and OH&S issues.
- **Optimised Heel Protection -** a 12° tapered heel zone redistributes pressure up the lower leg, away from the vulnerable heel region. The castellated surface collapses in a wave-like movement reducing shear.
- Integrated Shear Reduction the memory foam overlay moves independently of the underlying cells, reducing shear as the bed profiles.
- **Memory Foam Top Layer** integrated memory foam overlay improves immersion, envelopment and comfort.

Code	Description	Height	Width	Length	No. Of Cells	Static Head Cells	Cell Cycle	Cycle Time	S.W.L	Warranty	
PTM638010	Hybrid - Single	170 mm	895 mm	1980 mm	9	3	1 in 2	11 mins	180 kg	5Years Foam /	
PTM638015	Hybrid - King Single	170 mm	1050 mm	1980 mm	9	3	1 in 2	11 mins	180 kg	2 Years Cover	
PTM638030	Digital Pump	230 mm	320 mm	130 mm						2 Years	



Lifecomfort Hybrid Acute Static & Alternating Mattress

Lifecomfort Hybrid

Static & Alternating Mattress



The Lifecomfort Hybrid Acute (230kg) Mattress Systems offer foam-in-cell technology - enabling flexible operation as either a high comfort foam static support surface or powered-up alternating pressure support surface.

- Intelligent Hose Management when not in use the air-hose can be stored internally, reducing infection control and OH&S issues.
- **Optimised Heel Protection -** a 12° tapered heel zone redistributes pressure up the lower leg, away from the vulnerable heel region. The castellated surface collapses in a wave-like movement reducing shear.
- **Integrated Shear Reduction** the memory foam overlay moves independently of the underlying cells, reducing shear as the bed profiles.
- **Memory Foam Top Layer** integrated memory foam overlay improves immersion, envelopment and comfort.

Code	Description	Height	Width	Length	No. Of Cells	Static Head Cells	Cell Cycle	Cycle Time	S.W.L	Warranty	
PTM638020	Hybrid Acute - Single	170 mm	895 mm	1980 mm	9	3	1 in 2	11 mins	230 kg	5 Years Foam /	
PTM638025	Hybrid Acute - King Single	170 mm	1050 mm	1980 mm	9	3	1 in 2	11 mins	230 kg	2 Years Cover	
PTM638030	Digital Pump	230 mm	320 mm	130 mm						2 Years	



Lifecomfort Active Air 4 Mattress Overlay



The Active Air 4 is a lightweight, easy to use Overlay that is placed on top of an existing foam mattress. Featuring COMFORT CONTROL and offering a suite of other features, the Active Air 4 caters to users considered to be "At Risk" of Pressure Injury. An ideal solution for Home Care or Facility environments.

- Breathable multi-stretch cover
- Integrated cable management
- Simple and easy to operate
 - semi-digital pump

- Patented cable lock protection
- Welded seams and zip protection
- 4" Overlay + Memory Foam

Code	Description	Height	Width	Length	No. Of Cells	Static Head Cells	Cell Cycle	CycleTime	S.W.L	Warranty
PTM639100	Active Air 4" Mattress Overlay	100 mm	850 mm	2000 mm	19	3	1 in 2	10 / 15 / 20 mins	135 kg	2yr/1yr



Air Mattresses Bedroom





Lifecomfort Active Air 5 Mattress Overlay



🗊 aspire

lifecomfort

The Active Air 5 provides an effective solution for higher weight users considered to be "At Risk" of Pressure Injury. Featuring AUTO DETECT and INTELLIGENT STATIC MODE which provide lower interface pressure and enhanced comfort, the Active Air 5 is suitable for use in a wide range of care settings from Home Care and Long Term care.

- Breathable multi-stretch cover
- Integrated cable management
- Simple and easy to operate digital pump
- Patented cable lock protection
- Welded seams and zip protection
- 5" Overlay + Foam Mattress

Code	Description	Height	Width	Length	No. Of Cells	Static Head Cells	Cell Cycle	CycleTime	S.W.L	Warranty
PTM639200	Active Air 5" Overlay - Single	130 mm	850 mm	2000 mm	16	3	1 in 2	10/15/20 mins	175 kg	2yr/1yr
PTM639210	Active Air 5" Overlay - King Single	130 mm	1050 mm	2000 mm	16	3	1 in 2	10/15/20 mins	175 kg	2yr/1yr



Lifecomfort Active Air 8 Mattress Replacement

The Active Air 8 is a dynamic mattress replacement suitable for users at "High Risk" of Pressure Injury. Featuring AUTODETECT, and AUTOPROTECT, together with CELLIN-CELL technology, the Active Air 8 provides the highest levels of care into Long Term care environments. The Active Air 8's innovative features ensure lowest interface pressures with an easy-to-use digital pump.

- Breathable multi-stretch cover
- Integrated cable management
- Simple and easy to operate digital pump
- Patented cable lock protection
- Welded seams and zip protection
- 8" Mattress replacement

Code	Description	Height	Width	Length	No. Of Cells		Cell Cycle		S.W.L	Warranty
PTM639300	Active Air 8" Mattress Replacement - Single	200 mm	900 mm	2000 mm	20	3	1 in 2	10/15/20/25/30 mins	200 kg	2yr/1yr
PTM639310	Active Air 8" Mattress Replacement - King Single	200 mm	1050 mm	2000 mm	20	3	1 in 2	10/15/20/25/30 mins	200 kg	2yr/1yr





Pressure Care Cushions



Smart Check by ROHO® is a sensing device used in the setup and monitoring of a ROHO HIGH PROFILE® or MID PROFILE™ Sensor Ready Cushion



High Profile Single Valve No Cover

Note: This range can fit chairs from $13.25'' \times 13.25''$ to $22'' \times 22''$. Please contact Aidacare for enquiries relating to sizes outside of the displayed range. Sizes listed are seat dimensions.

All models include

- 2 Way Stretch Cover
- Hand Pump & Repair Kit
- Operating Instructions
- 2 Years Warranty

ROHO[®] - Single Compartment

- Appropriate for clients with very high risk of ischaemic ulcers and for management and recovery from existing pressure injury conditions
- Interconnected cells help disperse pressure and reduce shear risk promoting pressure management

High Profile (HP)

- 4" (10cm) High Cells
- Higher Pressure Disbursement

Low Profile (LP)

- 2" (5cm) High Cells
- Greater Stability



- and durable
- Maximum User Weight: Unlimited

	e		DEPTH (Ce	ells / Size")	
WIDTH (Cells / Size")	Profile	12 / 22"	11 / 20.25"	10 / 18.5"	9 / 16.75"
10 / 00//	HP	-	PTC544560	PTC544500	PTC544420
12 / 22″		PTC544600	PTC544590	PTC544550	PTC544430
11 / 00 05 //	HP	-	PTC544300	PTC544200	PTC544400
11 / 20.25″		-	PTC544350	PTC544250	PTC544410
10 / 10 5"	HP	-	PTC543900	PTC543700	PTC544000
10 / 18.5″		-	PTC543910	PTC543800	PTC544100
9 / 16.75″	HP	-	PTC546100	PTC545900	PTC546600
		-	PTC546200	PTC546000	PTC546700

ROHO® - Dual Compartment

- Appropriate for clients with very high risk of ischaemic ulcers and for management and recovery from existing pressure injury conditions
- Two separate compartments of the cushion allow positioning adjustments to address a pelvic obliquity (side-to-side) or sacral sitting

High Profile (HP)

- 4" (10cm) High Cells
- Higher Pressure Disbursement

Low Profile (LP)

- 2" (5cm) High Cells
- Greater Stability



- The Dual Compartment offers slightly less postural control and flexibility than the Quadtro Select range but retains the pressure and shear risk management standard for ROHO
- Maximum User Weight: Unlimited

All models include

- 2 Way Stretch Cover
- Hand Pump & Repair Kit
- Operating Instructions
- 2 Years Warranty

	ie	DEPTH (Cells / Size")					
WIDTH (Cells / Size")	Profil	11 / 20.25"	10 / 18.5"	9 / 16.75"			
10 / 18.5″	HP	PTC547290	PTC547300	PTC547500			
10 / 10.5	LP	PTC547295	PTC547400	PTC547505			
9 / 16.75″	HP	-	PTC548000	PTC548400			
3710.75	LP	-	PTC548100	PTC548500			



Note: This range can fit chairs from 13.25" x 13.25" to 20.25" x 20.25." Please contact Aidacare for enquiries relating to sizes outside of the displayed range.

Sizes listed are seat dimensions. The actual cushion is 1.5" longer to allow for the cushion to sit between the back canes.

Low Risk M Medium Risk High Risk

V

Very High Risk

ROHO Pressure Care Cushions Cushions





Note: This range can fit chairs from 14" x 14" to 20 x 20". Please contact Aidacare for enquiries relating to sizes outside of the displayed range.

Sizes listed are seat dimensions. The actual cushion is 1.5" longer to allow for the cushion to sit between the back canes.

ROHO® - Air Lite

- Appropriate for everyday comfort and pressure injury prevention
- · Designed for clients with tactile sensation who are at low risk of developing skin breakdown
- · Contoured foam base with an internal sealed air cell provides comfort and support for basic positioning and pressure management needs
- Maximum User Weight: 113kg

All models include

2 Way Stretch Cover

- Hand Pump & Repair Kit
- Operating Instructions
- 1 Year Warranty

			DEPTH	(Size")		
WIDTH (Size")	20″	18″	17″	16″	15″	14″
20″	-	PTC551450	-	-	-	-
19″	-	-	PTC551440			
18″	PTC551430	PTC551400	-	PTC551380	-	-
17″	-	PTC551330	PTC551325	-	-	-
16″	-	PTC551300	PTC551320	PTC551310	-	-
15″	-	-	-	-	PTC551295	-
14″	-	-	-	PTC551290	-	PTC551285

ROHO® - Hybrid Elite

- Appropriate for clients with moderate to high risk of skin breakdown who require positioning support
- Single valve type
- The Hybrid Elite combines the skin protection and shape matching ability of ROHO Dry Floatation with a contoured foam base
- The design allows for improved pelvic alignment and stability while maintaining



optimal immersion and dispersion of pressure away from the pelvic area

- The ROHO Dry Floatation system and 2-way stretch cover help reduce friction and shear whilst maximising blood flow
- The stability provided by the foam base assists "friction & shear" with safer transfers
- Maximum User Weight: 225kg

	DEPTH (Size")								
WIDTH (Size")	20″	18″	17″	16″	15″	14″			
24″	PTC594949	PTC594948	-	-	-	-			
22″	PTC594947	PTC594946	-	-	-	-			
20″	PTC594951	PTC594950	-	-	-	-			
18″	PTC594945	PTC594940	PTC594939	PTC594935	-	-			
17″	-	PTC594931	PTC594930	PTC594929	-	-			
16″	PTC594921	PTC594920	PTC594911	PTC594910	PTC594909	-			
15″	-	-	PTC594901	PTC594900	PTC594899	-			
14″	-	-	-	PTC594860	PTC594850	PTC594845			



Note: This range can fit chairs from 13" x 13" to 22" x 22". Please contact Aidacare for enquiries relating to sizes outside of the displayed range

Medium Risk

Sizes listed are seat dimensions

All models include


- 2 Way Stretch Cover
- Hand Pump & Repair Kit
- Operating Instructions
- 3 Years Warranty

Low Risk

High Risk v Very High Risk

ROHO® Quadtro Select

- Appropriate for clients at very high risk of ischaemic ulcers and for management and recovery from existing pressure injury conditions
- ROHO Dry Floatation promotes circulation, and pressure relief whilst flexible neoprene cells reduce shear friction
- Quick-and-easy air compartment adjustment maximises function and utility



• Maximum User Weight: Unlimited

High Profile (HP)

- 4" (10cm) High Cells
- Higher Pressure Disbursement

• 2" (5cm) High Cells
 Greater Stability

Low Profile (LP)

	e			DEPTH (Ce	ells / Size″)		
WIDTH (Cells / Size")	Profile	12 / 22″	11 / 20″	10 / 18.5″	9 / 16.75″	8 / 15″	7 / 13″
12 / 22″	HP	PTC589900	PTC589800	PTC589600	PTC589550	-	-
12/22	LP	-	-	PTC589700	-	-	-
11 / 20″	HP	PTC589200	PTC589000	PTC588800	PTC589400	-	-
11 / 20	LP	PTC589210	PTC589100	PTC588900	PTC589500	-	-
10 / 18.5″	HP	PTC588200	PTC588000	PTC587700	PTC588500	PTC588300	-
10 / 18.5	LP	PTC588210	PTC588100	PTC587900	PTC588700	PTC588400	PTC588280
9 / 16.75″	HP	-	PTC592400	PTC592200	PTC593000	PTC592800	PTC592700
9/10.75	LP	PTC588000	PTC592500	PTC592300	PTC593200	PTC592900	PTC592710



Note: This range can fit chairs from 14" x 14" to 20" x 20". Please contact Aidacare for enquiries relating to sizes outside of the displayed range

 Operating Instructions 2 Years Warranty

Sizes listed are seat dimensions. The actual cushion is 1.5" longer to allow for the cushion to sit between the back canes.





ΔΔ





Promedicare Inserto Domino

- Foam cushion designed for those who need higher positioning support
- Ideal for individuals with leg discrepancies, wind swept, pelvic obliquity, pelvic rotations, severe posterior/anterior pelvis tilt
- 17 modular and configurable inserts as standard
- Flexible open celled foam with low spring back action
- Comes with washable cover
- 2 Years Warranty (Cushion), 6 Months for (Cover)

Cushion Covers

Code	Description	Code	Description
PTC608280	10 x 12in	PTC608287	16 x 16in
PTC608281	12 x 12in	PTC608288	16 x 18in
PTC608282	12 x 15in	PTC608289	18 x 18in
PTC608283	14 x 15in	PTC608290	18 x 20in
PTC608284	14 x 16.5in	PTC608291	19 x 18in
PTC608285	16.5 x 18in	PTC608292	19 x 20in
PTC608286	16.5 x 20in		

		DEPTH (Size")					
WIDTH (Size")	12"	15"	16"	18"	19.5"	20″	
10"	PTC608015	-	-	-	-		
12"	PTC608020	PTC608025	-	-	-		
14"	-	PTC608030	-	-	-		
16"	-	-	PTC608040	PTC608045	-		
16.5"	-	PTC608035	-	-	PTC608036	PTC608037	
18"	-	-	-	PTC608050	PTC608055		
19"	-	-	-	PTC608056	-	PTC608057	



Promedicare Inserto Tanto

- Foam cushion designed for those who need medium/high postural support
- Ideal for individuals with mild leg length discrepancies, wind-swept, pelvic obliquity, pelvic
 - rotations, posterior and anterior pelvis tilt
- 13 modular and configurable inserts as standard
- Flexible open celled foam with low spring back action
- Comes with washable cover
- 2 Years Warranty (Cushion), 6 Months Warranty (Cover)

Cushion Covers

Shion C	Jovers		
Code	Description	Code	Description
PTC608320	14 x 15in	PTC608325	16.5 x 18in
PTC608321	14 x 16.5in	PTC608326	18 x 18in
PTC608322	16 x 16in	PTC608327	18 x 20in
PTC608323	16 x 18in	PTC608328	19in x 18in
PTC608324	16.5 x 20in		

For Customisable Cushions Please Refer To Our

COMPLEX REHAB CUSHION RANGE

On Page 34 Of This Catalogue

Complex Rehab Wheelchair Range on Page 32 SULT WITH OUA SOULT WITH OUA MARKEN SOULT MARKEN SOULTA

LOOK FOR THE COMPLEX REHAB ICON ON THE TOP RIGHT CORNERS OF PAGE

📘 Low Risk 🚺 Medium Risk Η High Risk 🚺 Very High Risk





P200 Portable Ceiling Hoist

SWL 150KG

Naidacare | A150F

- Lightweight portable hoist
- Suitable for use with fixed ceiling tracks for room to room transfers or mobile gantry systems
 Quick Release System (QRS) allows for the P200 to
- be safely and easily attached or removed from the gantry or track system with an easy reacher
- Safety locking latches to secure slings in place

PRISM MEDICAL

- Perfect for use in community or nursing environments
- Pneumatic handset is ideal for use in wet environments
- Emergency lowering allows for the safe lowering of the in emergency situations

Code	Description	Height Range	Weight	S.W.L
LIC350120	Ceiling Hoist - Prism - P200 (P440) Portable Manual Traverse 200kg & Carry Bar	Up to 2070 mm	4.8 kg	200 kg



Transactive Xtra Ceiling Hoist

- Fixed ceiling hoist allows for a large array of configurations
- Digital display screen provides essential information for the user and caregiver
- Return to Charge (RTC) option automatically returns the hoist to it's charger when required
- Manual wind down and raising allows the hoist to continue to be used during a power outage
- Smooth lifting/running tape with large lifting range

PRISM MEDICAL

- Pneumatic handset is ideal for use in wet environments
- Low profile hoist unit is unobtrusive and provides greater lifting height
- Emergency Lowering allows for the safe lowering of the in emergency situations
- Safety locking latches to secure slings in place
- Available in manual or powered traverse options

Code	Description	Total Height	S.W.L
LIC350150	ManualTraverse & Carry bar - 160 kg	Up to 2080 mm	160 kg
LIC350170	ManualTraverse & Carry bar - 200 kg	Up to 2080 mm	200 kg
LIC350200	ManualTraverse & Carry bar - 270 kg	Up to 2080 mm	270 kg
LIC350180	Powered Traverse & Carry bar - 200 kg	Up to 2080 mm	200 kg

		RollOn Standing c	and Positic	ning Aid			🛛 Liko
(· ·	Sit to stand mobile aid supp mobilising themselves	orting patients ir		afety lock ensures rear supports remain in an id osition		
0		le security and a a lightweight de:	• Lo	-	base for optimal ac tors ensure stable	, ,	
Code	Description		Total Length	Total Width	Total Height	Material	SWL
LSS393360	Liko RollOn		890 mm	620 - 900 mm	1410 mm	Aluminium	160 kg







Aspire A150F Aluminium Folding Lifter

- Folds quickly and conveniently for storage or transport
- Slimline handset technology from linak fits neatly in the hand with an intuitive control panel
- Battery and service indicator lights assist with service scheduling
- Off set leg design provides additional accessibility around chairs and furniture
- Simple foot operation to open or close legs for lifting stability or accessibility
- Quick release yoke attachment: one touch release and attach mechanism, designed for the busy caregiver. Padded yoke for added patient safety

- Strong powder coated aluminium frame with industry leading linak electronics and warranty provision
- Critical features include emergency stop button, emergency lowering and anti-entrapment technology, with both patient and caregiver safety at heart
- State of the art linak jumbo care control system. Tracks usage and service requirements for safety and compliance
- 2 X 100mm locking castors rear locking castors for safety and maneuverability
- Tested and compliant to AS/NZS ISO 10535:2011

Code	Display Name	Total Length	Total Width	Total Height	Lift Range	Leg Gait Range	Material	SWL
LSS390600	A150F Folding Patient Lifter	1220 - 1280 mm	660 - 1040 mm	1873 mm	428 - 1640 mm	580 - 950 mm	Aluminium	150 kg



Aspire A205 Aluminium Lifter

- Transportation or Storage: Is quick and convenient with tool free disassembly
- Ergonomic Handset: The latest slimline technology from Linak fits neatly in the hand with an intuitive control panel. Battery and service indicator lights assist with service scheduling
- Adjustable Height: Three adjustable height settings provides the caregiver with the flexibility to adjust the height range to suit the specific lift and patient scenario
- Powered Leg Function: Opens and closes the legs for lifting stability or accessibility
- Quick Release Yoke Attachment: One touch release and attach mechanism, designed for the busy caregiver. Padded yoke for added patient safety



🦪 aspire

- Build Quality: Strong powder coated aluminium frame with industry leading Linak electronics and warranty provision
- Safety: Critical features include emergency stop button, emergency lowering and antientrapment technology, with both patient and caregiver safety at heart
- Intelligent Diagnostics: State of the art Linak Jumbo Care control system. Tracks usage and service requirements for safety and compliance
- Quality Locking Castors: 2 x 100mm rear locking castors for safety and manoeuvrability
- Item Weight: 38.1 kg
- Tested and compliant to AS/NZS ISO 10535:2011

Code	Description	Total Length	Total Width	Total Height	Lift Range	Leg Gait Range	Material	SWL
LSS390610	A205 Patient Lifter	1135 mm	620 - 1095 mm	2045 mm	395 - 1805 mm	547 - 1015 mm	Aluminium	205 kg



Aspire 200 Steel Patient Lifter

- Ergonomic Handles Provide multiple grip points to suit carers of different heights and strengths
- Ergonomic Handset The latest slimline technology from Linak fits neatly in the hand with an intuitive control panel. Battery and service indicator lights assist with service scheduling
- Powered Leg Function Opens and closes the legs for lifting stability or accessibility
- Quick Release Yoke Attachment One touch release and attach mechanism, designed for the busy caregiver. Padded yoke for added patient safety
- Sturdy steel frame with industry leading Linak electronics and warranty provision



- Safety Critical features include emergency stop button, emergency lowering and anti-entrapment technology, with both patient and caregiver safety at heart
- Intelligent Diagnostics State of the art Linak Jumbo Care control system. Tracks usage and service requirements for safety and compliance
- Quality Locking Castors 2 x 100mm rear locking castors for safety and manoeuvrability
- Item Weight: 43.9 kg
- Tested and compliant to AS/NZS ISO 10535:2011

Code	Description	Total Length	Total Width	Total Height	Lift Range	Leg Gait Range	Material	SWL
LSS390620	200 Patient Lifter	1280 mm	665 - 980 mm	1940 mm	535 - 1725 mm	580 - 900 mm	Steel	200 kg





Aspire 200S Steel Standing Lifter

- Bullhorn Lifting Handles for patient stability and 6 sling attachment points
- Ergonomic Handles Provide multiple grip points to suit carers of different heights and strengths
- Ergonomic Handset The latest slimline technology from Linak fits neatly in the hand with an intuitive control panel. Battery and service indicator lights assist with service scheduling
- Powered Leg Function Opens and closes the legs for lifting stability or accessibility
- Safety Critical features include emergency stop button, emergency lowering and anti-entrapment technology, with both patient and caregiver safety at heart

• Sturdy steel frame with industry leading Linak electronics and warranty provision

🦪 aspire

🦪 aspire

🦪 aspire

- Adjustable Padded Legrests Provide comfort and support during transfers & can be adjusted to suit the patient
- Intelligent Diagnostics State of the art Linak Jumbo Care control system. Tracks usage and service requirements for safety and compliance
- Quality Locking Castors 2 x 100mm rear locking castors for safety and manoeuvrability
- Item Weight: 44.6 kg
- Tested and compliant to AS/NZS ISO 10535:2011

Code	Description	Total Length	Total Width	Total Height	Lift Range	Leg Gait Range	Material	SWL
LSS390650	200S Standing Patient Lifter	1150 mm	656 - 1016 mm	1790 mm	850 - 1740 mm	575 - 936 mm	Steel	200 kg

Aspire A200S Aluminium Standing Lifter

- Adjustable Depth Yoke: Yoke is depth adjustable for optimal patient fit
- Ergonomic Handles: Provide multiple grip points to suit carers of different heights and strengths
- Ergonomic Handset: The latest slimline technology from Linak fits neatly in the hand with an intuitive control panel. Battery and service indicator lights assist with service scheduling
- Powered Leg Function: Opens and closes the legs for lifting stability or accessibility
- Build Quality: Strong powder coated aluminium frame with industry leading Linak electronics and warranty provision
- Safety: Critical features include emergency stop

button, emergency lowering and antientrapment technology, with both patient and caregiver safety at heart


- Adjustable Padded Legrest: Provide comfort and support during transfers & can be adjusted to suit the user
- Intelligent Diagnostics: State of the art Linak Jumbo Care control system. Tracks usage and service requirements for safety and compliance
- Quality Locking Castors: 2 x 75mm rear locking castors for safety and manoeuvrability
- Item Weight: 30.4 kg
- Tested and compliant to AS/NZS ISO 10535:2011

Code	Description	Total Length	Total Width	Total Height	Lift Range	Leg Gait Range	Material	SWL
LSS390640	A200S Standing Patient Lifter	985 mm	590 - 915 mm	1455 mm	890 - 1395 mm	490 - 815 mm	Aluminium	200 kg



Aspire A320 Aluminium Bariatric Lifter

- Ergonomic Handles Provide multiple grip points to suit carers of different heights and strengths
- Ergonomic Handset The latest slimline technology from Linak fits neatly in the hand with an intuitive control panel. Battery and service indicator lights assist with service scheduling
- Powered Leg Function Opens and closes the legs for lifting stability or accessibility
- Quick Release Yoke Attachment One touch release and attach mechanism, designed for the busy caregiver. Padded yoke for added patient safety
- Build Quality Strong powder coated aluminium frame with industry leading Linak electronics and warranty provision



- Intelligent Diagnostics State of the art Linak Jumbo Care control system. Tracks usage and service requirements for safety and compliance
- Quality Locking Castors 2 x 125mm rear locking castors for safety and manoeuvrability
- Item Weight: 63 kg
- Tested and compliant to AS/NZS ISO 10535:2011

Code	Description	Total Length	Total Width	Total Height	Lift Range	Leg Gait Range	Material	SWL
LSS390630	A320 Bariatric Patient Lifter	1505 mm	765 - 1190 mm	2165 mm	575 - 1960 mm	640 - 1065 mm	Aluminium	320 kg









Oxford Advance Folding Homecare Lifter

- Mobile patient lifter for transfer to a seated or lying position
- Compact design can be folded easily for storage or transportation
- Epoxy powder coated frame
- · Single speed handset control • Manual foot operated leg gating and
- tapered leg design for close access
- Emergency lowering
- Over-sized push handles for excellent manoeuvrability
- Suits a large range of patient slings
- Supplied with battery and charger



Code	Description	Length	Width	Height	Leg Gait	Lift Range	Product Weight	S.W.L	Warranty
LSS377300	Advance Folding Lifter	1300 mm	550 mm	1360 - 1860 mm	570 - 1000 mm	350 - 1860 mm	31.7 kg	155 kg	2 Years



Oxford Journey Folding Standing Lifter

- Mobile compact patient lifter for transfer to a seated, lying or standing
- Compact and lightweight allowing carers to manoeuvre in tight spaces
- Folds away for easy storage or
- Unique adjustable 'cow-horn' shape supports a wide range of patient
- Adjustable knee pad
- Epoxy coated frame
- Electric hand control operation
- Emergency operation
- Easily manoeuvrable with 125mm rear locking castors and 75 mm front twin swivel castors
- Electric leg gating

Code	Description	Length	Width	Height	Leg Gait	Lift Range	Prod. Weight	S.W.L	Warranty
LSS377305	Journey Standing Lifter	1000 mm	685 mm	1130 - 1600 mm	590 - 875 mm	720 - 1530 mm	40.2 kg	155 kg	2 Years





- seated or lying position
- High lifting range accomodates transfers to higher surfaces
- Large load capacity of 227kg
- Robust epoxy powder coated frame construction makes it ideal for facilities
- Electric hand control operation
- Emergency lowering and raising

- Over-sized push handles for excellent manoeuvrability
- 100mm rear locking castors and 100mm front twin swivel castors
- Electric leg gating
- Suits a large range of patient slings
- Supplied with battery and charger
- Pictured with optional 4 point pivot frame (LSS377396)

Code	Description	Length	Width	Height	Leg Gait	Lift Range	Prod. Weight	S.W.L	Warranty
LSS377325	Presence Lifter	1460 mm	670 mm	1460 - 2061 mm	550 - 1020 mm	440 - 1910 mm	40 kg	227 kg	2Years





Oxford Stature Lifter

- The Stature lifter recognises the true needs of the modern care environment, able to cope with even the most demanding and technical of patient handling situations
- Safe working load of 227kgs and great lifting ranges available (min 390mm, max 1550mm)
- The Stature can lift patients from the floor as well as accommodate transfers to higher surfaces, while allowing the choice of different sling systems
- Careful attention to product design ensures correct ergonomic usage, protecting both the patient and caregiver from injury
- Manual or powered cradles are available to assist with repositioning of patients
- The powered base allows the caregiver to open and close the legs without bending or twisting
- Ergonomic foot push-pad to assist with initiating movement

Code	Description	Length	Width	Height	Leg Gait	Lift Range	Product Weight	S.W.L	Warranty
LSS377373	Stature Lifter	1390 mm	1360 mm	1945 - 2130 mm	640 - 1230 mm	390 - 1550 mm	73.5 kg	227 kg	2 Years





Orbi-turn Standing Patient Turner

Designed for caregivers to facilitate standing and seat-to-seat transfers while minimising the risk of injury to themselves and the patient. Allows a client with moderate upper and lower limb strength to participate in functional transfers and ensure good working posture for therapists. The turntable effortlessly rotates through 360° for repositioning in even the most confined spaces. Nonslip footprints provide secure foot placement and show clients where to position feet.

- Overall height adjustment of 813 1320mm accommodates variety of patients
- Height-adjustable and angled kneepad results in normal biomechanical hip and knee movement
- Height-adjustable hand bar suits a range of clients
- Thigh pad gives the patient a secure, supportive contact point for increased stability and comfort during transfers.
- Easy-to-operate foot brake prevents rotation and overcompensating on balance during sit-to-stand movement.
- Two rear casters provide easy transport from client to client.
- Front-mounted sprung caster assists in positioning of the Orbi-Turn once the client's feet are in place.
- Folds flat for easy portability and storage.
- Steel frame ensures durability and strength.
- Weighs 15kg.

Code	Description	Height	Product Weight	S.W.L
TSE672935	Orbi-turn Standing Patient Turner	813 - 1320 mm	15 kg	190 kg





Aspire Slings

7	- hand and
oire	Jack
gs	Jort

All slings have reinforc colour coded for eact making them easy to t	n size available	SMALL MEDIUM	LARGE X-LARGE	XX- LARGE XXX- LARGE
Size	Edging	Girth	Height	S.W.L
Small	•	65 cm - 95 cm	120 cm - 160 cm	200 kg
Medium	•	77cm - 110 cm	140 cm - 175 cm	200 kg
Large	•	98 cm - 135 cm	158 cm - 185 cm	300 kg
X-Large		118 cm - 157 cm	165 cm - 193 cm	300 kg
XX-Large	•	155 cm - 184 cm	170 cm - 193 cm	300 kg
XXX-Large		180 cm - 210 cm	170 cm - 193 cm	300 kg



Aspire Deluxe General Purpose Sling

Versatile sling offering body contoured fabric

• For FIM score 1 or 2



🦪 aspire

shaped for	shaped for maximum patient comfort and support						
	Code	Description	Size	Edge	S.W.L	Wty	
	LSS390700	Polyester - Loop	Small		200 kg	1 Year	
	LSS390705	Polyester - Loop	Medium	•	200 kg	1 Year	
Polyester	LSS390710	Polyester - Loop	Large		300 kg	1 Year	
	LSS390715	Polyester - Loop	X-Large		300 kg	1 Year	



Aspire Deluxe General Purpose Sling with Head Support

- This sling has a comfortable padded head support with rigidisers for posture management
- Body contoured fabric
 Ear EIM agora 1 or 2
- For FIM score 1 or 2

	Code	Description	Size	Edge	S.W.L	Wty
	LSS390740	Polyester - Loop	Small		200 kg	1 Year
	LSS390745	Polyester - Loop	Medium		200 kg	1 Year
Polyester	LSS390750	Polyester - Loop	Large		300 kg	1 Year
,	LSS390755	Polyester - Loop	X-Large		300 kg	1 Year





Aspire Deluxe Standing Sling

- Intuitively designed, this sling offers generous underarm padding and a reinforced wide-cut back for comfort.
- Quilted and padded under arm to enhance patient comfort

Polyes

comfort

F

- · Padded back support with stays for stability and rigidity
- Security strap features velcro and buckle fastening for easy positioning and added safety.

🦪 aspire

🦪 aspire

• For FIM score 3 or 4

	Code	Description	Size	Edge	S.W.L	Wty
	LSS390720	Padded Polyester - Loop	Small		200 kg	1 Year
	LSS390725	Padded Polyester - Loop	Medium		200 kg	1 Year
ster	LSS390730	Padded Polyester - Loop	Large		300 kg	1 Year
	LSS390735	Padded Polyester - Loop	X-Large		300 kg	1 Year



Aspire Deluxe Hygiene Access Sling

quilted leg straps and a reinforced wide-cut back.

- Security strap features velcro and buckle
- fastening for easy positioning and added safety
- For FIM score 2 or 3

	Code	Description	Size	Edge	S.W.L	Wty
	LSS390760	Padded Polyester - Loop	Small		200 kg	1 Year
	LSS390765	Padded Polyester - Loop	Medium	•	200 kg	1 Year
Polyester	LSS390770	Padded Polyester - Loop	Large		300 kg	1 Year
	LSS390775	Padded Polyester - Loop	X-Large		300 kg	1 Year



Aspire Deluxe Hygiene Access Sling with Head Support

- This sling offers an adjustable head support for additional patient safety and comfort whilst maintaining all the benefits of the Deluxe Hygiene Access Sling
- Padded poly fabric and generous under arm padding for increased patient comfort



- Intuitive shape allows easier access for undressing and washing
- · Security belt with velcro and buckle
- For FIM score 2 or 3

	Code	Description	Size	Edge	S.W.L	Wty
	LSS390780	Padded Polyester - Loop	Small		200 kg	1 Year
	LSS390785	Padded Polyester - Loop	Medium	•	200 kg	1 Year
Polyester	LSS390790	Padded Polyester - Loop	Large		300 kg	1 Year
	LSS390795	Padded Polyester - Loop	X-Large		300 kg	1 Year



Aspire Deluxe Standing Sling

- This full body sling provides complete comfort and support for high needs patients.

 - Supersoft for long term use

- 🧊 aspire
- Parasilk leg strap lining for easy positioning • For FIM score 1 or 2

	Code	Description	Size	Edge	S.W.L	Wty
	LSS390800	Super Soft - Loop	Small		200 kg	1 Year
	LSS390805	Super Soft - Loop	Medium	•	200 kg	1 Year
Super Soft	LSS390810	Super Soft - Loop	Large		300 kg	1 Year
	LSS390815	Super Soft - Loop	X-Large		300 kg	1 Year





Oxford Slings





Standing Sling

- The standing sling is suitable for those patients who have a greater degree of weight bearing ability
- It allows excellent access for toileting and is easy to fit, providing quick and effective transfers
- Adjustable waist strap and non-slip back pad ensures the sling does not ride up during transfer
- Can be unsafe however, when used with the wrong patient, so a careful risk assessment is advised
- Available in three sizes small, medium and large
- FIM Scores: 3 and 4

Code	Description	S.W.L
LSS371952	Standing Sling - Small	200 kg
LSS371951	Standing Sling - Medium	200 kg
LSS371950	Standing Sling - Large	200 kg



Quickfit Deluxe & Quickfit Deluxe with Head Support

- This sling is a development of the Quickfit design and is suitable for 95% of patients
- Can be used for amputees following assessment, and is available in sizes S to L
- FIM Scores: 1 and 2 and any other person from the floor
- Both Quickfit Deluxe and Quickfit Deluxe with Head Support are available in polyester
- Quickfit Deluxe is also available in mesh.

Code	Description	S.W.L
LSS371698	Quickfit Deluxe Poly - Small	227 kg
LSS371632	Quickfit Deluxe Poly - Medium	227 kg
LSS371699	Quickfit Deluxe Poly - Large	227 kg
LSS371697	Quickfit Deluxe Mesh - Small	227 kg
LSS371700	Quickfit Deluxe Mesh - Medium	227 kg
LSS371702	Quickfit Deluxe Mesh - Large	227 kg
LSS371635	Quickfit Deluxe Poly w/Head Support - Small	227 kg
LSS390490	Quickfit Deluxe Poly w/Head Support - Medium	227 kg
LSS371642	Quickfit Deluxe Poly w/Head Support - Large	227 kg



Comfort Sling (for Pivot Frame)

- This is a more specialised sling
- Allows correct positioning to be made via the 4-point cradle and uses the Securi3 sling connection system ensuring no inadvertent detachment of the sling from the cradle
- Incorporates a removable comfort pad for head support and fits to give full protection for patients who go into extension or have involuntary movements or behavioural problems
- FIM Scores: 1 and 2 and any other person from the floor

Code	Description	S.W.L
LSS371639	Comfort Sling Poly - Extra Small	227 kg
LSS371643	Comfort Sling Poly - Small	227 kg
LSS371644	Comfort Sling Poly - Medium	227 kg
LSS371645	Comfort Sling Poly - Large	227 kg
LSS371640	Comfort Sling Poly - Extra Large	227 kg

Seat Walkers / Rollators



Aspire Vogue Indoor Walker / Rollator

- Luxury look and feel Danish design, profiled tubing, moulded fittings, ergonomic handles and refined colours.
- Now supplied with bonus cane holder (ref. image)
- Lightweight and ergonomic design with small turning circle for excellent manoeuvrability.
- No lifting & spillage meals can be collected and brought to the desired location without the need for lifting or carrying. The tray is also removable.
- Sitting down for a meal when the user is sitting, the walker can be turned around and positioned comfortably in front of them, to enjoy their meal.



- Unique folding design the slimline design reduces the footprint when folded but allows the folded walker to stand upright on its own.
- Brakes the ergonomic 'Soft Touch' lockable design allows the walker to be secured in place with minimal hand pressure.
- Removable storage bag utensils and other items can be stowed securely away, ideal for transporting personal belongings around the house.

Code	Description	Total Height	Unit Weight	SWL
WAF705100CH	Aspire Vogue Indoor Walker – Champagne/Black	860 - 970 mm	5.2 kg	110 kg
WAF705100RE	Aspire Vogue Indoor Walker - Red/Black	860 - 970 mm	5.2 kg	110 kg



Aspire Vogue Lightweight 2 Seat Walker

🌒 aspire

- Luxury look and feel Danish design, profiled tubing, moulded fittings, ergonomic handles and refined colours.
- Now supplied with bonus cane holder (ref. image)
- Lightweight and ergonomic design with small turning circle for excellent manoeuvrability.
- "AutoLocate" Smooth height adjustmentthat automatically locates to pre-set heights, making adjustment simple and accurate.
- Locks when folded makes transportation and storage easier.
- Advanced locking handbrakes for increased user safety when walking and sitting.
- Unbreakable TPE super soft wheels offers high durability and a smooth ride for increased comfort.
- Removable utility bag simple and stylish to remove for shopping. Can hold up to 5 kg.
- New wider curved backrest

Code	Description	Total Length	Total Width	Total Height	Seat Width	Unit Weight	SWL
WAF705310RE	Vogue Lightweight 2 - Red	655 mm	620 mm	785 - 920 mm	465 mm	6.1 kg	150 kg
WAF705310CH	Vogue Lightweight 2 - Champagne	655 mm	620 mm	785 - 920 mm	465 mm	6.1 kg	150 kg





Aspire Vogue Forearm Seat Walker / Rollator

- Smart Posture System equipped with forearm supports that are height, depth and angle adjustable, the user can achieve optimal centre of gravity for improved balance and weight distribution.
- Luxury look and feel Danish design, profiled tubing, moulded fittings, ergonomic handles and refined colours.
- Lightweight and ergonomic design with small turning circle for excellent manoeuvrability.





27 aspire

- "AutoLocate" Smooth height adjustmentthat automatically locates to pre-set heights, making adjustment simple and accurate.
- Locks when folded makes transportation and storage easier.
- Advanced locking handbrakes for increased user safety when walking and sitting.
- Unbreakable TPE super soft wheels offers high durability and a smooth ride for increased comfort.

Code	Description	Seat Width	Total Height	SeatTo Floor Height	Unit Weight	SWL
WAF705400	Aspire Vogue Forearm Walker	460 mm	1020 - 1270 mm	620 mm	9 kg	150 kg



Aspire Push Down Seat Walker / Rollator

- Robust aluminium frame and easy wipe-clean surfaces.
- Weight activated brakes ideal for arthritic conditions. May also may ideal for those with cognitive issues who have difficulty applying brakes (risk assessment advised)
- Anatomically designed, latex-free handles allow for comfortable grip, especially for users with arthritis

• Ultra-soft hand brakes with palm ball - lock

with arthritis or wrist/hand joint pain

effortlessly and offer increased comfort for users

Anatomically designed, latex-free handles - allow

for comfortable grip, especially for users with

surfaces.

arthritis.

- Adjustable handle height for customised comfort and posture.
- Easy folding mechanism to facilitate storage and transportation.
- Padded seat that flips up with a discreet pocket for smaller items and a storage bag supplied.
- Solid, puncture-proof 8" wheels for increased manoeuvrability and peace of mind.

Code	Description	Seat Width	Total Height	Seat To Floor Height	Unit Weight	SWL
WAF705600RE	Aspire Push Down Seat Walker - Red	355 mm	851 - 962 mm	580 mm	7 kg	130 kg
WAF705600BL	Aspire Push Down Seat Walker - Blue	355 mm	851 - 962 mm	580 mm	7 kg	130 kg



Aspire Classic 8" Seat Walker / Rollator

- Robust aluminium frame and easy wipe-clean • Adjustable handle height - for customised comfort and posture.
 - Easy folding mechanism to facilitate storage and transportation.
 - Padded seat that flips up with a discreet pocket for smaller items and a storage bag supplied.
 - Solid, puncture-proof 8" wheels for increased manoeuvrability, particularly outdoors.

Code	Description	Seat Width	Total Height	Seat To Floor Height	Unit Weight	SWL
WAF705700RE	Aspire Classic 8 Seat Walker - Red	355 mm	817 - 972 mm	580 mm	7 kg	130 kg
WAF705700BL	Aspire Classic 8 Seat Walker - Blue	355 mm	817 - 972 mm	580 mm	7 kg	130 kg



Aspire Classic 6" Seat Walker / Rollator

- Robust aluminium frame and easy wipe-clean surfaces
- Ultra-soft hand brakes with palm ball lock effortlessly and offer increased comfort for users with arthritis or wrist/hand joint pain
- Anatomically designed, latex-free handles allow for comfortable grip, especially for users with arthritis.
- Adjustable handle height for customised comfort and posture.
- Easy folding mechanism to facilitate storage and transportation.
- Padded seat that flips up with a discreet pocket for smaller items and a storage bag supplied.
- Solid, puncture-proof 6" wheels for increased manoeuvrability, particularly indoors.

Code	Description	Seat Width	Total Height	SeatTo Floor Height	Unit Weight	SWL
WAF709600RE	Aspire Classic 6 Seat Walker - Red	355 mm	760 - 915 mm	528 mm	6 kg	130 kg
WAF709600BL	Aspire Classic 6 Seat Walker - Blue	355 mm	760 - 915 m m	528 mm	6 kg	130 kg

27 aspire

🗊 aspire



Aspire Deluxe Seat Walker / Rollator

- Deluxe, in-line rear wheels and brake mechanism increases the user's stepping zone.
- Robust aluminium frame and easy wipe-clean surfaces.
- Ultra-soft hand brakes with palm ball lock effortlessly and offer increased comfort for users with arthritis or wrist/hand joint pain
- Adjustable handle height for customised comfort and posture.





- Anatomically designed, latex-free handles allow for comfortable grip, especially for users with arthritis.
- Easy folding mechanism to facilitate storage and transportation.
- Padded seat that flips up with a discreet pocket for smaller items and a storage bag supplied.
- Solid, puncture-proof 8" wheels for increased manoeuvrability and peace of mind.

Code	Description	Seat Width	Total Height	Seat To Floor Height	Unit Weight	SWL
WAF705750RE	Aspire Deluxe Seat Walker - Red	355 mm	834 - 989 mm	587 mm	7 kg	130 kg



Aspire Mini Seat Walker / Rollator

- Robust aluminium frame and easy wipe-clean surfaces.
- Low 485mm seat height to suit shorter users.
- Ultra-soft hand brakes with palm ball lock effortlessly and offer increased comfort for users with arthritis or wrist/hand joint pain
- Anatomically designed, latex-free handles allow for comfortable grip, especially for users with arthritis.
- Adjustable handle height for customised comfort and posture.
- Easy folding mechanism to facilitate storage and transportation.
- Padded seat that flips up with a discreet pocket for smaller items and a storage bag supplied.
- Solid, puncture-proof 6" wheels for increased manoeuvrability and peace of mind.

Code	Description	Seat Width	Total Height	Seat To Floor Height	Unit Weight	SWL
WAF705800RE	Aspire Mini Seat Walker - Red	355 mm	743 - 905 mm	490 mm	6 kg	130 kg
WAF705800BL	Aspire Mini Seat Walker - Blue	355 mm	743 - 905 mm	490 mm	6 kg	130 kg



Aspire XL Seat Walker / Rollator

- Robust aluminium frame and easy wipe-clean surfaces.
- Wide padded seat and increased 180kg Safe Working Load to suit larger users.
- Ultra-soft hand brakes with palm ball lock effortlessly and offer increased comfort for users with arthritis or wrist/hand joint pain
- Anatomically designed, latex-free handles allow for comfortable grip, especially for users with arthritis.

🦪 aspire

🦪 aspire

🧊 aspire

- Adjustable handle height for customised comfort and posture.
- Easy folding mechanism to facilitate storage and transportation.
- Padded seat that flips up with a discreet pocket for smaller items and a storage bag supplied.
- Solid, puncture-proof 8" wheels for increased manoeuvrability and peace of mind.

Code	Description	Seat Width	Total Height	SeatTo Floor Height	Unit Weight	SWL
WAF750020SL	Aspire XL Seat Walker - Silver	495 mm	850 - 982 mm	560 mm	8 kg	180 kg



Forearm Walker - HD Gutter Frame

- Provides assistance for those who have difficulty walking
- Tutors the user into an upright gait
- Folding design for storage or transportation
- Height adjustable forearm supports to suit user requirements
- Adjustable hand grips to suit user forearm length
- Four swivel castors for manoeuvrability
- Powder coated mild steel construction
- Skis or glides available as an option as well as padded arm rest
- Large size designed for adults

Code	Description	Total Height	SWL
WAF699200	Forearm Walker - HD Gutter Frame	1000 - 1270 mm	180 kg



Aspire Manual Wheelchairs





Aspire Evoke 2 + Evoke 2 HD Wheelchair



- The Evoke 2 and Evoke 2 HD are designed for the user who wants a lightweight, highly adjustable manual wheelchair.
- Swing away and removable arm rests this facilitates user transfers and also makes for easier transport and storage.
- Swing-in / out, height adjustable and removable foot rests – this facilitates user transfers and allows optimal positioning for the changing needs.
- Strong, Nylon upholstery with a comfortable cushion and padded backrest
- For increased comfort, a padded overlay is available as an accessory.
- Depth adjustable seating with rear angle adjustability – for more accurate positioning, customisability and increased user comfort.
- Quick release 24" solid rear tyres the wheels are easily removed to reduce the weight of the chair for transport and storage.
- Crash Tested to AS/NZ3616.196

Code	Description	Seat Width	Overall Width	Seat Depth	Seat Height	Backrest Height	Armrest Height	S.W.L
MWS449500	Evoke 2 - 350 mm (14")	350 mm (14")	560 mm	400 - 475 mm	440 - 480 mm	430 - 490 mm	205 - 225 mm	140 kg
MWS449510	Evoke 2 - 400 mm (16")	400 mm (16")	610 mm	400 - 475 mm	440 - 480 mm	430 - 490 mm	205 - 225 mm	140 kg
MWS449520	Evoke 2 - 450 mm (18")	450 mm (18")	660 mm	400 - 475 mm	440 - 480 mm	430 - 490 mm	205 - 225 mm	140 kg
MWS449530	Evoke 2 - 500 mm (20")	500 mm (20")	710 mm	400 - 475 mm	440 - 480 mm	430 - 490 mm	205 - 225 mm	140 kg
MWS449540	Evoke 2 (HD) - 570 mm (22")	570 mm (22")	780 mm	400 - 475 mm	440 - 480 mm	430 - 490 mm	205 - 225 mm	180 kg
MWS449550	Evoke 2 (HD) - 610 mm (24")	610 mm (24")	820 mm	400 - 475 mm	440 - 480 mm	430 - 490 mm	205 - 225 mm	180 kg



Aspire Evoke 2 Junior Wheelchair

- The Evoke 2 Junior is designed to promote independent manoeuvring and handling by juniors, while making lifting and handling a breeze for parents and caregivers.
- Swing away and removable arm rests this facilitates user transfers and also makes for easier transport and storage.
- Swing-in / out, height adjustable and removable foot rests – this facilitates user transfers and allows optimal positioning for the changing needs.

🦪 aspire

- Strong, Nylon upholstery with a comfortable cushion and padded backrest.
- Depth adjustable seating with rear angle adjustability – for more accurate positioning, customisability and increased user comfort.
- Quick release 16" solid rear tyres the wheels are easily removed to reduce the weight of the chair for transport and storage.
- Crash Tested to AS/NZ3616.196

Code	Description	Seat Depth	Seat Width	Backrest Height	Seat To Floor Height	Unit Weight	SWL
MWS449560	Aspire Evoke 2 Junior - 300mm Wide	300 - 375 mm	300 mm	430 - 490 mm	360 - 410 mm	14 kg	140 kg

Aspire Tilt-In-Space Wheelchairs



Aspire Rehab RX Advanced Tilt-In-Space Wheelchair

- Designed for the more complex user, who requires assistance with positioning, repositioning and skin management.
- Tilt (26°) and backrest recline (30°) this facilitates weight re-distribution and the preservation of skin integrity. It also optimises user comfort, decreases user fatigue and assists with head control.
- Contoured, built-up backrest and seat cushion

 the backrest has a gentle lumbar build up
 and moderate lateral supports; together
 with the abductor and adductor build-up
 at the seat cushion, midline positioning is
 optimised.



- Advanced articulating legrests closely mimic the normal movement at the knee upon elevation, allowing the thigh to maintain contact with the seating surface, ensuring maximum support and skin contact with the cushion.
- Breathability support surfaces offer multi-layered foam with breathable multistretch covers - for pressure re-distribution, comfort and to enhance skin microclimate.
- Highly adjustable depth and width adjustable seat sizing. Height adjustable and removable armrests. Adjustable headrest.
- Crash Tested to AS/NZ3616.196

Code	Description	Seat Depth	Seat Width	Backrest Height	Seat To Floor Height	Unit Weight	SWL
MWS449710	Rehab RX - 410-460mm Wide	440 - 530 mm	410 - 460 mm	460 - 710 mm	400 - 475 mm	41.0 kg	180 kg
MWS449720	Rehab RX - 460-510mm Wide	440 - 530 mm	460 - 510 mm	460 - 710 mm	400 - 475 mm	43.0 kg	180 kg
MWS449740	Rehab RX - 510-560mm Wide	470 - 560 mm	510 - 560 mm	460 - 710 mm	400 - 475 mm	47.2 kg	180 kg



Aspire Rehab RX Junior Wheelchair

- Designed for the increased postural support and pressure care of children, via easy repositioning and advanced adjustability.
- Tilt (26°) and backrest recline (30°) this facilitates weight re-distribution and the preservation of skin integrity. It also optimises user comfort, decreases user fatigue and assists with head control.
- Highly adjustable depth and width adjustable seat sizing. Height adjustable and removable armrests. Adjustable headrest.
- Contoured, built-up backrest and seat cushion

 the backrest has a gentle lumbar build up
 and moderate lateral supports; together

with the abductor and adductor build-up at the seat cushion, midline positioning is optimised.

- Breathability support surfaces offer multilayered foam with breathable multi-stretch covers - for pressure re-distribution, comfort and to enhance skin microclimate.
- Advanced articulating legrests closely mimic the normal movement at the knee upon elevation, allowing the thigh to maintain contact with the seating surface, ensuring maximum support and skin contact with the cushion.
- Crash Tested to AS/NZ3616.196

Code	Description	Seat Depth	Seat Width	Backrest Height	Seat To Floor Height	Unit Weight	SWL
MWS449630	Rehab RX Junior - 320-370mm Wide	320 - 370 mm	320 - 370 mm	460 - 710 mm	400 - 475 mm	39.5 kg	180 kg



🥑 asoire



Aspire Rehab RS Classic Tilt-In-Space Wheelchair

- The Rehab RS is designed for a user who has limited ability to reposition themselves in a device and as a result, is at risk of pressure injuries.
- Tilt (26°) and backrest recline (30°) this facilitates weight re-distribution and the preservation of skin integrity. It also optimises user comfort, decreases user fatigue and assists with head control.
- Highly adjustable depth and width



adjustable seat sizing. Height adjustable and removable armrests. Adjustable headrest.

- Breathability support surfaces offer multilayered foam with breathable multi-stretch covers - for pressure re-distribution, comfort and to enhance skin microclimate.
- Heavy duty, padded, elevating leg rests supplied as standard.
- Crash Tested to AS/NZ3616.196

Code	Description	Seat Depth	Seat Width	Backrest Height	SeatTo Floor Height	Unit Weight	SWL
MWS449600	Rehab RS - 410-460mm Wide	440 - 530 mm	410 - 460 mm	460 - 710 mm	400 - 475 mm	40.7 kg	180 kg
MWS449610	Rehab RS - 460-510mm Wide	440 - 530 mm	460 - 510 mm	460 - 710 mm	400 - 475 mm	42.5 kg	180 kg
MWS449620	Rehab RS - 510-560mm Wide	470 - 560 mm	510 - 560 mm	460 - 710 mm	400 - 475 mm	46.7 kg	180 kg

Aspire Wheelchair Accessories



Code	Description
WCA848330	WChair Contoured Padded Overlay - 350mm - Aspire Evoke2
WCA848340	WChair Contoured Padded Overlay - 400mm - Aspire Evoke2
WCA848350	WChair Contoured Padded Overlay - 450mm - Aspire Evoke2
WCA848360	WChair Contoured Padded Overlay - 500mm - Aspire Evoke2
WCA848370	WChair Contoured Padded Overlay - 570mm - Aspire Evoke2
WCA848380	WChair Contoured Padded Overlay - 610mm - Aspire Evoke2
WCA848100	W'Chair Elevating Legrest - Standard - Left - 320/350/400/410mm - Aspire Assist2/Transit2/Evoke2/Rehab RS/Rehab RX
WCA848110	W'Chair Elevating Legrest - Standard - Right - 320/350/400/410mm - Aspire Assist2/Transit2/Evoke2/Rehab RS/ Rehab RX
WCA848120	W'Chair Elevating Legrest - Standard - Left - 450/460/500/510mm - Aspire Assist2/Transit2/Evoke2/Rehab RS/Rehab RX
WCA848130	W'Chair Elevating Legrest - Standard - Right - 450/460/500/510mm - Aspire Assist2/Transit2/Evoke2/Rehab RS/Rehab RX
WCA848140	W'Chair Elevating Legrest - Standard - LEFT- 570/610mm - Aspire Assist2/Evoke3
WCA848150	W'Chair Elevating Legrest - Standard - Right - 570/610mm - Aspire Assist2/Evoke2
WCA848160	W'Chair Elevating Legrest - Advanced - Left - 320/350/400/410mm - Aspire Assist2/Transit2/Evoke2/Rehab RS/Rehab RX
WCA848170	W'Chair Elevating Legrest - Advanced - Right - 320/350/400/410mm - Aspire Assist2/Transit2/Evoke2/Rehab RS/Rehab RX
WCA848180	W'Chair Elevating Legrest - Advanced - Left - 450/460/500/510mm - Aspire Assist2/Transit2/Evoke2/Rehab RS/Rehab RX
WCA848190	W'Chair Elevating Legrest - Advanced - Right - 450/460/500/510mm - Aspire Assist2/Transit2/Evoke2/Rehab RS/Rehab RX
WCA848220	W Chair Stump Support - Left - Silver - 320/350/400/450mm - Aspire Assist2/Transit2/Evoke2
WCA848230	W Chair Stump Support - Right - Silver - 320/350/400/450mm - Aspire Assist2/Transit2/Evoke2
WCA848240	W Chair Stump Support - Left - Silver - 500/550/570/610mm - Aspire Assist2/Transit2/Evoke2/Rehab RS/Rehab RX
WCA848250	W Chair Stump Support - Right - Silver - 500/550/570/610mm - Aspire Assist2/Transit2/Evoke2/Rehab RS/Rehab RX
WCA848270	W Chair AntiTip Bar - Aspire Evoke2/Assist2/Transit2/Rehab RS/Rehab RX
WCA848300	W'Chair Lap Belt - Metal Seatbelt Buckle - Aspire Assist2/Transit2/Evoke2/Rehab RS/Rehab RX
GSP492408	W'Chair Legrests - Swing Away - 350mm wide - Aspire Evoke2
GSP492414	W'Chair Legrests - Swing Away - 400mm wide - Aspire Evoke2
GSP492418	W'Chair Legrests - Swing Away - 450mm wide - Aspire Evoke2
GSP492422	W'Chair Legrests - Swing Away - 500mm wide - Aspire Evoke2
GSP492424	W'Chair Legrests - Swing Away - 570/610mm wide - Aspire Evoke2



Brochure Resources Online



Aspire Care Website



Product Application Guides Available on Request











