

# NSW HEALTH C963

**WALKERS & CRUTCHES CATALOGUE** 

**JANUARY 2020** 

VERSION: 2 | CREATED JANUARY 2020



Quality Equipment Enterprise Solutions Home & Community





### Aidacare & NSW Health

Aidacare are proud to be recognised as a NSW HealthShare / EnableNSW contracted service provider for a large range of products included in the recently announced Beds and Mattresses contract. Our commitment to all NSW HealthShare stakeholders including prescribing health care professionals and most importantly, the NSW HealthShare client, is to provide the very best service and product quality available.

For the healthcare prescriber, Aidacare offer industry leading support to assist in the provision of equipment for your clients.

- All contracted products are available for trial, anywhere in New South Wales, to enable assessment of suitability for your clients specific needs. Our equipment consultants can meet with you at any location, including client's private homes, to demonstrate and trial items.
- Dedicated NSW HealthShare display stock located in each of the Aidacare showrooms is available to enable prescribers to meet with clients on-site to trial and assess all required items in a professional environment.
- Many Aidacare showrooms even offer prescribers use of a private office and meeting room to utilise while meeting with clients or just to catch up on admin and reports.
- Specialist Complex Rehab & Seating staff to assist and support with prescription of customised and specialised equipment.
- Clinical education and support via the Aidacare Professional Training program and workshops.
- Service and repair of all contracted items with our workshop and mobile technicians.

General customer service and support including the booking of trials and any general advice and information required.

Our electronic catalogue is available for download at www.aidacare.com.au/globalassets/resources

For information on any of the extensive range of products available from Aidacare, please browse the complete online Aidacare catalogue at www.aidacare.com.au

For additional information, please contact your local Aidacare office directly on 1300 133 120 or email nswhealth@aidacare.com.au

EnableNSW provides equipment and services to people in NSW with chronic health conditions or disability to assist them with mobility, communication and self-care. We also provide assistance for people who need to travel significant distances for specialist medical care.

Aids and Equipment Program - The Aids and Equipment Program is a NSW Government initiative to assist eligible residents of NSW, who have a life-long or long-term disability, to live and participate within their community by providing appropriate equipment, aids and appliances.

Information on eligibility for the program, the range of equipment available, the lodgement of applications, mutual responsibilities and expectations, and steps for complaints and appeals is outlined in the Policy Directive - "Assistive Technology for Communication, Mobility, Respiratory Function & Self-Care". https://www1.health.nsw.gov.au/PDS/pages/doc.aspx?dn=PD2011\_027

People with a disability interested in seeking further information about the program should discuss their query with a healthcare provider, or Enable NSW on 1800 ENABLE (1800 362 253)

As part of HealthShare NSW, EnableNSW is responsible for the administration of NSW Health disability support and other assistance programs including:

- Aids and Equipment Program
- Continence Assistance
- Home Respiratory Program (HRP)
- Prosthetic Limb Service (formerly known as ALS)
- Specialised Equipment Essential for Discharge (SEED)

For more information, visit the EnableNSW website: http://www.enable.health.nsw.gov.au/

## State Wide Distribution Network



3A / 1 Moorebank Ave, Moorebank NSW 2170 Tel: (02) 8706 2300 Fax: (02) 9618 5111

SEVEN HILLS
2/4 Stanton Rd

Seven Hills NSW 2147 Tel: (02) 8706 2330 Fax: (02) 9624 4256

**BLACKTOWN** 

1/34 Campbell St Blacktown NSW 2178 Tel: (02) 9621 4077

WOLLONGONG
140 Church St
Wollongong NSW 2

Wollongong NSW 2500 Tel: (02) 8706 2359 **ALBURY** 

347 Urana Rd, Lavington NSW 2641 Tel: (02) 6040 3444 Fax: (02) 6040 3222

BALLINA

18 Endeavour Close, Ballina NSW 2478 Tel: (02) 6686 6644 Fax: (02) 6686 9383

COFFS HARBOUR

2/177 Orlando St, Coffs Harbour NSW 2450 Tel: (02) 6651 5333 Fax: (02) 6652 1272 T

PORT MACQUARIE

1/171 Lake Rd, Port Macquarie NSW 2444 Tel: (02) 6581 0018 Fax: (02) 6581 0038

OI OI

**ORANGE** 

Unit 1, 12-16 Ash St, Orange NSW 2800 Tel: 0437 134 295

İ

**TAMWORTH** 

3/248 Peel St, Tamworth NSW 2340 Tel: (02) 6761 2018 Fax: (02) 6761 2090



CALL 1300 133 120

Call this number to be automatically diverted to your local Aidacare office (for the price of a local call).



## Product Index



-					
-	200		$\sim$	h	~
v	ΙL	ıυ	u		63

Underarm Crutches - Red Dot.4Underarm Crutches - Bariatric.4Underarm Crutches - Wingnut.4Forarm Crutches.5Forearm Crutches Ergonomic.5	Kids Forearm Crutches.5Gutter Crutches.5Crutch - Cushion.5Crutch - Hand Grip.5Ferrules / Rubber Tips.5
Walking Sticks	
Folding Walking Stick	Walking Stick Palm Grip
Seat Walkers / Rollators & Accessories	
Aspire Mini Seat Walker	Aspire Vogue Lightweight Seat Walker
Forearm Walkers	
Aspire Vogue Forearm Seat Walker	
Walking Frames & Accessories	
Aspire Side Folding Walking Frame	Heavy Duty Folding Walking Frame







#### **Underarm Crutches - Red Dot**

- Patented clip lock adjustable handle and push button foot pieces for fast, precise measuring and fitting
- 22mm ferrules

- High density rubber cushions and handgrips providing long wear and greater comfort
- Lightweight aluminium frame, child weighing 1.1kg each, and youth adult tall weighing 1.25kg each

Code	Description	User Height	Crutch Height	S.W.L	Warranty
WAC692300	Underarm Crutches - Child	1220 - 1370 mm	787 - 940 mm	136 kg	1 Year
WAC692500	Underarm Crutches - Youth	1320 - 1570 mm	940 - 1140 mm	136 kg	1 Year
WAC692200	Underarm Crutches - Adult	1550 - 1750 mm	1120 - 1320 mm	136 kg	1 Year
WAC692400	Underarm Crutches - Adult Tall	1780 - 1980 mm	1320 - 1520 mm	136 kg	1 Year



#### **Underarm Crutches - Bariatric**

- Designed with strength to support the bariatric user
- Steel construction provides durability and support up to 249 kg
- Lightweight frame weighing 4.7kg
- 22mm ferrule

Code	Description	User Height	Crutch Height	S.W.L	Warranty
WAC688300	Underarm Crutches - Bariatric Adult	1575 - 1778 mm	1120 - 1320 mm	249 kg	1 Year
WAC688200	Underarm Crutches - Bariatric Adult Tall	1778 - 1981 mm	1320 - 1520 mm	249 kg	1 Year



## **Underarm Crutches - Wing Nut**

- Easy to use Wing-Nut height adjustment
- High density rubber cushions and handgrips providing long wear and greater comfort
- Lightweight aluminium frame
- 22mm ferrules

Code	Description	User Height	Crutch Height	S.W.L	Warranty
WAC692503	Underarm Crutches Youth - Small	1370 - 1575 mm	940 - 1145 mm	136 kg	1 Year
WAC692502	Underarm Crutches - Medium	1575 - 1778 mm	1120 - 1320 mm	136 kg	1 Year
WAC692501	Underarm Crutches Adult - Tall	1778 - 1980 mm	1350 - 1520 mm	136 kg	1 Year





#### Forearm Crutches

- Lightweight (600g each) and durable forearm crutches to suit a large range of users
- Cuff provides a secure fitting to the forearm

•	Adjustable height between cuff and handle & handle
	and ground

• 22mm ferrules

Code	Description	User Height	Crutch Height	S.W.L	Warranty
WAC692000	Forearm Crutches Youth - Small	1230 - 1640 mm	530 - 785 mm	180 kg	1 Year
WAC691400	Forearm Crutches Adult - Medium	1480 - 1860 mm	635 - 910 mm	180 kg	1 Year
WAC691800	Forearm Crutches Adult - Tall Large	1740 - 2040 mm	820 - 1075 mm	180 kg	1 Year



#### Forearm Crutches Ergonomic

- Ergonomic shaped handle distributes the pressure more evenly across the hand
- Lightweight (715g each)and durable forearm crutches to suit a large range of users
- Cuff provides a secure fitting to the forearm
- Adjustable height between cuff and handle & handle and ground
- 22mm ferrules

Code	Description	Total Height	Crutch Height	S.W.L	Warranty
WAC691500	Forearm Crutches Ergonomic - Medium	1585 - 1800 mm	580 - 835 mm	127 kg	1 Year
WAC691900	Forearm Crutches Ergonomic - Large	1706 - 2010 mm	680 - 935 mm	127 kg	1 Year



#### **Kids Forearm Crutches**

- Just for kids
- These light weight yellow Forearm Crutches are height adjustable

Code	Description	Handle Height	Ferrule	S.W.L	Warranty
WAC691600	Crutches Forearm Child 6 - 10 years	457 - 656 mm	18 mm	100 kg	1 Year



## **Gutter Crutches**

- Ideal for those with limited hand or wrist strength
- Lightweight tubular aluminium stem and padded vinyl covered gutter for maximum comfort
- Designed for users with arthritis or an impaired grip
- Forearm trough distributes a user's weight over a large support area
- Adjustable telescopic handle rotates through 360° and locks to suit individual forearm length

Code	Description	Gutter Height	Ferrule	Product Weight	S.W.L	Warranty
WAC694200	Gutter Crutches	980 - 1235 mm	25 mm	950g (each)	127 kg	1 Year



#### Crutch - Hand Grip

• Replacement Cushion handle to suit **Underarm Crutches** 

Code	Description
WAA686220	Crutch Grey Hand Grip Padding - Red Dot
WAA686225	Crutch Padding Grey Hand Grip - Wing Nut Crutches
WAA686230	Crutch Grey Hand Grip Bariatric



#### Ferrules/Rubber Tips

• High quality rubber ferrule tips suitable for use with walking frames, crutches, walking sticks and quad sticks

Code	Description
WAA686320	Rubber Ferrule 22 mm, Grey
WAA686330	Rubber Ferrule 25 mm, Grey



#### Crutch - Cushion

 Replacement Axilla Cushion for **Underarm Crutches** 

Code	Description
WAA686240	Crutch Grey Axilla Cushion - Red Dot
WAA686250	Crutch Grey Axilla Bariatric Cushion
WAA686270	Crutch Padding Hand Grip Yellow





#### T Handle Walking Stick

- Adjustable height to suit user
- Polished hardwood handle and matte black aluminium stick

Code	Description	Handle Height	Ferrule	Product Weight	S.W.L	Warranty
WAS743600	T Handle W/Stick	715 - 970 mm	18 mm	365 g	100kg	1Year



### Swan Neck Walking Stick with Strap

- Adjustable height to suit user
- Aluminium stick with unique foam swan handle for maximum comfort and strap for added security

Code	Description	Handle Height	Ferrule	Product Weight	S.W.L	Warranty
WAS743700	Swan Walking Stick with Strap	730 - 985 mm	18 mm	365 g	100 kg	1Year
WAS743800	T Handle Swan Walking Stick Adjustable	710 - 965 mm	22 mm	550 g	127 kg	1Year
WAS744750	Bariatric Walking Stick	737 - 966 mm	28 mm	400 g	225 kg	1Year



#### Pyramid Walking Stick

- Adjustable height to suit user
- Pyramid base for extra stability
- Soft foam hand grip

Code	Description	Handle Height	Base Measurement	Product Weight	S.W.L	Warranty
WAS743900	Walking Stick Pyramid	790 - 990 mm	240 x 300 mm	800 g	100 kg	1Year





### Walking Stick Palm Grip

- Unique, biometrically shaped handle designed to evenly spread pressure across the palm of the hand
- Very useful for those with arthritis or poor grip
- Right and left hand available

Code	Description	Handle Height	Ferrule	Product Weight	S.W.L	Warranty
WAS744300	Walking Stick Palm Right Grip	670 - 925 mm	18 mm	400 g	127 kg	1 Year
WAS744400	Walking Stick Palm Left Grip	670 - 925 mm	18 mm	400 g	127 kg	1 Year



### **Quad Base Walking Stick**

- Adjustable height to suit user
- Four point base for extra stability
- Soft foam hand grip

Code	Description	Handle Height	Base Measurement	Ferrule	Product Weight	S.W.L	Warranty
WAS744600	Walking Stick Quad Small Base	1170 - 1630 mm	200 x 150 mm	12.5 mm	1 kg	113 kg	1Year
WAS744500	Walking Stick Quad Large Base	1470 - 1930 mm	300 x 200 mm	15.6 mm	1.2 kg	113 kg	1 Year









Adjustable Handle Height



Robust Aluminium Frame



Padded Seat Flips up





Red Blue



Ultra-soft Hand Brakes with Palm Ball Lock for Increased Comfort, Especially for Users with Arthritis



Curved and Padded Backrest Offers Seated Support



Discreet Pocket Compartment

Allows the User to Store Smaller Objects Safely



Solid, Puncture-Proof 6" Wheels







Code	Description	Overall Height	Overall Length	Overall Width	Seat Depth	Seat Width	Seat Height	Backrest Height	Folded Width	Castor Size	Unit Weight
WAF705800RE	Mini Seat Walker - Red	743 - 905 mm	556 mm	615 mm	305 mm	355 mm	490 mm	690 mm	200 mm	6"	6 kg
WAF705800BL	Mini Seat Walker - Blue	743 - 905 mm	556 mm	615 mm	305 mm	355 mm	490 mm	690 mm	200 mm	6"	6 kg

## Seat Walkers / Rollators & Accessories



#### Aspire Classic 6" Seat Walker



Adjustable Handle Height



Robust Aluminium Frame



Padded Seat Flips up







Red

Blue

Anatomically Designed, Latex-Free Handles

Ultra-soft Hand Brakes with Palm Ball Lock for Increased Comfort, Especially for Users with Arthritis



Curved and Padded Backrest Offers Seated Support



Discreet Pocket Compartment

Allows the User to Store Smaller Objects Safely



Solid, Puncture-Proof 6" Wheels







Code	Description	Overall Height	Overall Length	Overall Width	Seat Depth	Seat Width	Seat Height	Backrest Height	Folded Width	Castor Size	Unit Weight
WAF709600RE	Classic 6" Seat Walker - Red	760 - 915 mm	624 mm	615 mm	305 mm	355 mm	528 mm	769 mm	200 mm	6"	6 kg
WAF709600BL	Classic 6" Seat Walker - Blue	760 - 915 mm	624 mm	615 mm	305 mm	355 mm	528 mm	769 mm	200 mm	6"	6 kg

#### Aspire Classic 8" Seat Walker



Adjustable Handle Height



Robust Aluminium Frame



Padded Seat Flips Up







Blue

Anatomically Designed, Latex-Free Handles

Ultra-soft Hand Brakes with Palm Ball Lock for Increased Comfort, Especially for Users with Arthritis



Curved and Padded Backrest Offers Seated Support



Discreet Pocket Compartment

Allows the User to Store Smaller Objects Safely



Solid, Puncture-Proof 8" Wheels







Code	Description	Overall Height	Overall Length	Overall Width	Seat Depth	Seat Width	Seat Height	Backrest Height	Folded Width	Castor Size	Unit Weight
WAF705700RE	Classic 8" Seat Walker - Red	817 - 972 mm	660 mm	628 mm	305 mm	355 mm	584 mm	815 mm	200 mm	8"	7 kg
WAF705700BL	Classic 8" Seat Walker - Blue	817 - 972 mm	660 mm	628 mm	305 mm	355 mm	584 mm	815 mm	200 mm	8"	7 kg



### Aspire Push Down Seat Walker





Anatomically Designed, Latex-Free Handles

Adjustable Handle Height



Padded Seat Flips up



Robust Aluminium Frame



Ideal for users with arthritis



Curved and Padded Backrest for Seated Support



Discreet Pocket Compartment for Storage

Allows the user to store smaller objects safely



Solid, Puncture-Proof 7" Wheels







Blue







Code	Description	Overall Height	Overall Length	Overall Width	Seat Depth	Seat Width	Seat Height	Backrest Height	Folded Width	Castor Size	Unit Weight
WAF705600RE	Push Down Seat Walker - Red	851 - 962 mm	656 mm	616 mm	305 mm	355 mm	580 mm	815 mm	255 mm	7″	7 kg
WAF705600BL	Push Down Seat Walker - Blue	851 - 962 mm	656 mm	616 mm	305 mm	355 mm	580 mm	815 mm	255 mm	7"	7 kg

### Aspire Superlite Adjustable Seat Walker



Compact Folding for Storage and Transport



Convenient Storage Bag Doubles as Seat

Height Adjustable Seat





Height Adjustable Handles



Tool Free Assembly

Swivelling 6" Front Wheels





1 Year Warranty



Code	Description	Overall Height	Overall Length		Seat Depth	Seat Width	Seat Height	Backrest Height		Castor Size	Unit Weight
WAF709650CH	Superlite Adjustable Seat Walker	740 - 910 mm	550 mm	600 mm	320 mm	350 mm	460 - 590 mm	700 mm	200 mm	6"	5 kg

## Seat Walkers / Rollators & Accessories



## Aspire XL Seat Walker





Adjustable Handle Height



Robust Aluminium Frame



Extra Large Padded Seat Flips Up



Wide Padded Seat Increased 180kg Safe Working Load



Ultra-soft Hand Brakes with Palm Ball Lock for Increased Comfort, Especially for Users with Arthritis



Curved and Padded Backrest for Seated Support



Discreet Pocket Compartment for Storage

Allows the User to Store Smaller Objects Safely



Solid, Puncture-Proof 8" Wheels



Silver



S.W.L 180kg



Code	Description		Overall Length	Overall Width	Seat Depth	Seat Width		Backrest Height	Folded Width	Castor Size	Unit Weight
WAF750020SL	XL Seat Walker	850 - 982 mm	697 mm	732 mm	340 mm	495 mm	560 mm	850 mm	230 mm	8"	8 kg

#### Aspire Deluxe Seat Walker



Anatomically Designed, Latex-Free Handles



Ultra-soft hand brakes with palm ball lock for increased comfort, especially for users with arthritis

Robust Aluminium Frame



Deluxe, in-line rear wheels and brake mechanism

Discreet Pocket Compartment



Allows the user to store smaller objects safely





Adjustable Handle Height

#### Curved and Padded Backrest



Padded Seat Flips up



Solid, Puncture-Proof 8" Wheels

Puncture-proof wheels for increased manoeuvrability over various surfaces





Red

Blue

S.W.L 130kg

Code	Description	Overall Height	Overall Length	Overall Width	Seat Depth	Seat Width	Seat Height	Backrest Height	Folded Width	Castor Size	Unit Weight
WAF705750RE	Deluxe Seat Walker - Red	834 - 989 mm	764 mm	625 mm	305 mm	355 mm	587 mm	836 mm	245 mm	8"	7 kg
WAF705750BL	Deluxe Seat Walker - Blue	834 - 989 mm	764 mm	625 mm	305 mm	355 mm	587 mm	836 mm	245 mm	8"	7 kg



### Aspire Vogue Lightweight Seat Walker



Advanced Locking Handbrakes





Cane Holder

Smooth Height Adjustment with AUTOLOCATE

Unbreakable TPE Super Soft Wheels



Luxury Look and Feel

New Wider Curved Backrest

For Comfort and Support



Removable Utility Bag



220mm Folding Width Minimises Storage Space. Locks when Folded







Red







Code	Description	Overall Height	Overall Length	Overall Width	Folded Width	Seat Depth	Seat Width	Seat Height	Wheel Size		Bag S.W.L
WAF705310RE	Vogue Walker - Red	785 - 920 mm	655 mm	620 mm	220 mm	240 mm	465 mm	540 mm	8"	6.1 kg	5 kg
WAF705310CH	Vogue Walker - Champagne	785 - 920 mm	655 mm	620 mm	220 mm	240 mm	465 mm	540 mm	8"	6.1 kg	5 kg

#### Drive Bariatric Rollator with Locking Brakes



Anatomically designed handles allow comfortable grip

Arthritic friendly looped cable brakes activate effortlessly

Padded seat flips up easily for access to storage basket or for folding Padded backrest

Fitted with solid 20cm puncture proof tyres for easy manoeuvrability and durability



Handle height can be adjusted between 88cm 100cm to suit the user

Foldable for easy storage and transportation

Includes a basket for storing essential items

S.W.L
227kg

Code	Description	Unit Dimensions	Seat Dimensions	Total Weight	S.W.L
WAF750100	Drive Bariatric Rollator w/locking Brakes	679 mm (W) x 774 mm (L)	457 mm (W) x 558 mm (H) x 330 mm (D)	12 kg	227 kg



#### Aspire Vogue Indoor Walker



Code	Description	Overall Height	Overall Length	Overall Width	Folded Width	Wheel Size	Total Weight	Bag S.W.L	Tray S.W.L	S.W.L
WAF705100RE	Vogue Walker - Indoor - Red/Black	860 - 970 mm	650 mm	540 mm	270 mm	6"	5.2kg	5kg	5kg	110kg
WAF705100CH	Vogue Walker - Indoor - Champagne/Black	860 - 970 mm	650 mm	540 mm	270 mm	6"	5.2kg	5kg	5kg	110kg

#### Seat Walker Accessories



Holder for Seat Walker

Attaches to A-Frame Seat Walkers. Top mount clamp to securely hold the walking stick in place WAA686925



Accessory to suit Aspire Deluxe Walkers



Accessory to suit Aspire Classic, Mini, Deluxe, Push-down & XL Seat Walker



WAA686920 - Aspire XL

IV Pole - Height Adjustable Accessory to suit Aspire VOGUE Walker & VOGUE Walker 2 WAA695650 - VOGUE Walker



Accessory to suit Aspire Vogue & VOGUE 2 Walkers & Rehab Wheelchairs



Accessory to suit Aspire VOGUE & VOGUE2 Walkers WAA695630





#### Aspire Meal Tray Walker



Code	Description	Overall Height	Overall Length	Overall Width	LowerTray Height	UpperTray Height	Tray Width	Castor Size	Total Weight
DLG272061	Meal Tray Walker	770 - 920 mm	560 mm	585 mm	430 mm	696 mm	450 mm	6"	8 kg



### Aspire Vogue Forearm Seat Walker





Advanced Locking Handbrakes



Smooth Height Adjustment with AUTOLOCATE

Unbreakable TPE Super Soft Wheels

User Comfort

Super-soft Wheels and Padded Forearm Surfaces Absorb Bumps and Jolts Creating a Smooth, Comfortable Ride



Smart Posture System

Height, Angle and Depth Adjustable Forearm Supports Achieve Optimal Centre of Gravity for Improved Balance and Weight Distribution



Locks when Folded

220mm Folding Width Minimises Storage Space



Rehab Focus

Easy to Push with a Small Turning Circle







Code	Description	Overall Height	Overall Length	Overall Width	Folded Width	Seat Depth	Seat Width	Seat Height	Wheel Size	Total Weight	Bag S.W.L
WAF705400	Vogue Forearm Walker - White	1020 - 1270 mm	660 mm	620 mm	225 mm	240 mm	460 mm	620 mm	8"	9 kg	5 kg



### Aspire Walking Tutor HD with Brakes



- Provides assistance for those who have difficulty walking
- · For users who have limited hand or wrist strength, can weight bear through forearms
- Unique lightweight (5kg) folding design for storage or transportation (300 x 640mm)
- · Adjustable hand grips to suit user forearm length
- · Height adjustable forearm supports to suit user requirements
- Skis available as an option instead of rear rubber tips

Code	Description	Depth	Height of Forearm Support	S.W.L	Wty
WAF699200	Walking Tutor Forearm Gutter Frame - HD / Brakes	840 mm	1030 - 1270 mm	180 kg	1Year



## Aspire Side Folding Walking Frame





## Walking Frames & Accessories Mobility Walking Frames & Aids







Wheels & Ski's model

#### Aluminium Non-Folding Walking Frame

- Provides assistance for those who have difficulty walking
- Lightweight aluminium frame ideal for indoor use
- Height adjustable to suit user requirements

<ul> <li>Supplied</li> </ul>	with	four	rubber	tips
------------------------------	------	------	--------	------

- Ferrule of 25mm
- · Lightweight at only 1.9kg

Code	Description	Width	Handle Ht.	S.W.L	Warranty
WAF704900	A/Walking Frame - Small	560 mm	665 - 740 mm	160 kg	1 Year
WAF704800	A/Walking Frame - Medium	560 mm	765 - 840 mm	160 kg	1 Year
WAF704700	A/Walking Frame - Large	560 mm	865 - 940 mm	160 kg	1 Year
WAF704900K	A/Walking Frame with Wheels & Ski's - Small	560 mm	665 - 740 mm	160 kg	1 Year
WAF704800K	A/Walking Frame with Wheels & Ski's - Medium	560 mm	765 - 840 mm	160 kg	1 Year
WAF704700K	A/Walking Frame with Wheels & Ski's - Large	560 mm	865 - 940 mm	160 kg	1 Year





- Provides assistance for those who have difficulty walking
- Folds flat for easy transportation and storage
- Lightweight aluminium frame ideal for indoor use
- Height adjustable to suit user requirements
- Supplied with four rubber tips
- Ferrule of 25mm
- · Lightweight at only 2kg



Code	Description	Width	Handle Ht.	S.W.L	Warranty
WAF706700	A/ Folding Walking Frame - Small	610 mm	710 - 787 mm	127 kg	1 Year
WAF706800	A/ Folding Walking Frame - Medium	610 mm	787 - 860 mm	127 kg	1 Year
WAF706900	A/ Folding Walking Frame - Large	610 mm	860 - 935 mm	127 kg	1 Year
WAF706700K	A/ Folding Walking Frame w/Wheels & Ski's - Small	610 mm	710 - 787 mm	127 kg	1 Year
WAF706800K	A/ Folding Walking Frame w/Wheels & Ski's - Medium	610 mm	787 - 860 mm	127 kg	1 Year
WAF706900K	A/ Folding Walking Frame w/Wheels & Ski's - Large	610 mm	860 - 935 mm	127 kg	1 Year



### Heavy Duty Folding Walking Frame

- Heavy Duty Frame for durability
- Folds flat for easy transportation and storage
- Extra wide frame width to suit larger users
- Height adjustable to suit user requirements
- Ferrule of 25mm
- Product weight at 8.3kg

Code	Description	Width	Handle Ht.	S.W.L	Warranty
WAF704745	Heavy Duty Folding Walking Frame	740 mm	810 - 990 mm	227 kg	1 Year
WAF704750	Heavy Duty Folding Walking Frame with 2 Front Wheels	740 mm	810 - 990 mm	227 kg	1 Year

#### Walking Frame Accessories





#### Glide White

Accessory to suit Aluminium Folding Walking Frames

WAA686690 - 19mm WAA686710 - 25mm WAA686730 - 28.6mm



Two Piece Ski - Grey (Pair)

Accessory to suitany rollator-style Walking Frames

WAA686540 - 25.4mm



Wheels for Walking Frame (Pair)

Accessory to suit Aluminium Walking Frames





1300 133 120 www.aidacare.com.au



:A40.3 090120