SWEP: SCRIPT FORM

K CARE SHOWER COMMODES



REFERRER DETAILS		CLIENT DETAILS	
Referrer Name:		First Name:	Last Name:
Position:		Client DOB:	
Phone:	Fax:	Aidacare Complex Rehab Representativ	e:
Email:			

PLEASE FAX / EMAIL THIS FORM TO 9384 2088 OR SWEP@AIDACARE.COM.AU



SHOWER COMMODE MAXI

- Removable Padded Arms and Swing-away Footrest
- 4 Stainless Steel Locking Castors
- Includes Bed Pan Carrier
- Seat Height: 590 mm
- Frame Width: Transit is 610 mm and Self Propelled is 660 mm
- Seat Width: Transit is 550 mm and Self Propelled is 600 mm
- Seat Depth: Transit is 440 mm and Self Propelled is 450 mm

Code	Description	Max User Weight.	Select
BTC062100SP1	Transit Bariatric (Open Front)	190 kg	
BTC062100SP2	Transit Bariatric (Closed Front)	190 kg	
BTC062200SP1	Bariatric (Open Front)	210 kg	
BTC062200SP2	Bariatric (Closed Front)	210 kg	



SHOWER COMMODE SELF PROPELLING STD & LOW

- Self Propelled Rear Wheels with Swing-up Arms and Swing - Away Foot Rests
- 2 Stainless Steel Locking Castors
- Stainless Steel Frame for Corrosion Resistances
- Includes Bed Pan Carrier
- Seat Height: Standard is 590 mm and Low is 500 mm
- Frame Width: 510 mmSeat Width: 440 mmSeat Depth: 430 mm

Code	Description	Max User Weight	Select
BTC077900SP1	Standard (Open Front)	120 kg	
BTC077900SP2	Standard (Closed Front)	120 kg	
BTC077900SP3	Low (Open Front)	120 kg	
BTC077900SP4	Low (Closed Front)	120 kg	



SHOWER COMMODE SELF PROPELLING - FOLDING

- Refer to features of Self Propelling Std & Low
- Able to Fold for Easy Storage and transfer
- Seat Height: Standard is 590 mm and Low is 500 mm
- Frame Width: 510 mm
- Seat Width: 440 mm
- Seat Depth: 430 mm

C	Code	Description	Max User Weight	Select
Е	3TC080500SP1	Folding (Open Front)	120 kg	
Е	3TC080500SP2	Folding (Closed Front)	120 kg	

Need more information on SWEP products?

Call us to request a complete SWEP Catalogue

1300 133 120



SHOWER COMMODE TRANSIT - STD & LOW

- Attendant Propelled with Swing-up Arms and Sliding Foot Plate
- 4 Stainless Steel Locking Castors
- Lower Height at seat front
- Seat Height: Standard is 590 mm and Low is 500 mm
- Frame Width: 510 mm
- Seat Width: 440 mm
- Seat Depth: 430 mm

Code	Description	Max User Weight	Select
BTC070600SP1	Low (Open Front)	120 kg	
BTC070600SP2	Low (Closed Front)	120 kg	
BTC070600SP3	Standard (Open Front)	120 kg	
BTC070600SP4	Standard (Closed Front)	120 kg	



SHOWER COMMODE TRANSIT FOLDING

- Refer to features of Self Propelling Std & Low
- Able to Fold for Easy Storage and transfer
- Seat Height: 590 mm
- Frame Width: 510 mm
- Seat Width: 440 mm
- Seat Depth: 430 mm

Code	Description	Max User Weight	Select
BTC071400SP1	Folding (Open Front)	120 kg	
BTC071400SP2	Folding (Closed Front)	120 kg	



SHOWER COMMODE SWING AWAY FOOT REST

- Refer to features of Self Propelling Std & Low
- Swing Away Foot Rests
- Seat Height: 590 mm
- Frame Width: 510 mm
- Seat Width: 440 mm
- Seat Depth: 430 mm

Code	Description	Max User Weight	Select
BTC070900SP1	Swing Away (Open Front)	120 kg	
BTC070900SP2	Swing Away (Closed Front)	120 kg	



SHOWER COMMODE TRANSIT TILT IN SPACE

- Gas Lift Tilt-in-Space with Swing-Up Arms
- Padded Calf Rest and Swing-away Foot Rests
- Removable Padded Arms and Swing-away Footrest
- 4 Stainless Steel Frame
- Includes Bed Pan Carrier
- Seat Height: 590 mm
- Frame Width: 510 mm
- Seat Width: 440 mm
- Seat Depth: 430 mm

Code	Description	Max User Weight	Select
BTC071100SP1	Tilt in Space (Open Front)	135 kg	
BTC071100SP2	Tilt in Space (Closed Front)	135 kg	

SWEP: SCRIPT FORM

K CARE SHOWER COMMODES - ACCESSORIES



Please select any accessories that are required. If none, please leave blank and sign at the bottom of this form.

Code	Description	Price	~
BTC085000	Hand Made Open Seat 320 mm		
BTC073900	Hand Made Open Seat 400 mm		
BTC074400	Hand Made Open Seat 440 mm		
BTC062800	Hand Made Open Seat 500 mm		
BTC062900	Hand Made Open Seat 550 mm		
BTC063000	Hand Made Open Seat 600 mm		
BTC063100	Hand Made Open Seat 650 mm		
BTC084900	Hand Made Closed Seat 320 mm		
BTC073800	Hand Made Closed Seat 400 mm		
BTC073700	Hand Made Closed Seat 440 mm		
BTC062400	Hand Made Closed Seat 500 mm		
BTC062500	Hand Made Closed Seat 550 mm		
BTC062600	Hand Made Closed Seat 600 mm		
BTC062700	Hand Made Closed Seat 650 mm		
BTC082100	Pressure Reducing Closed Seat 440 mm		
BTC082000	Pressure Reducing Open Seat 440 mm		
BTC074100	Hand Made Wedge Seat 440 mm		
BTC075900	Padded Tie-On Closed Seat 440 mm		
BTC076600	Bowl Carrier, Stainless Steel		
BTC076500	Bowl Carrier Junior, Stainless Steel		
BTC081600	Bowl With Lid And Handle		
BTC077000	Pan Carrier Junior		
BTC077100	Pan Carrier Stainless Steel		
BTT151800	Bed Pan		
BTT114800	Bed Pan (Slipper)		
BTC074600	IV Pole		
BTC075000	Oxygen Bottle Carrier		
BTC074700	Lock Down Arms		
BTC071600	Folding Safety Arms		
BTC072900	Drop Side Arms		
BTC079900	Stump Support Left		
BTC080200	Stump Support Right		
BTC075100	Pommel		
BTC078300	Anti Tip Bars		
BTC075300	Fixed Curved Padded Back		
BTC075200	Padded Tie On Back		
BTC078400	Brake Extension		
BTC079600	Contoured Headrest		
BTC078900	Footplate Extended		
BTC080300	Swing Away Legrests		
BTC080350	Swing Away Legrests For Ka114s55		
BTC078500	Elevating Legrest Left		
BTC078700	Elevating Legrest Right		
BTC071800	Pelvic Strap - Fitted		
BTC077200	Pelvic / Chest Strap		
BTC072300	Padded Arms - Fitted		
BTC072200	Padded Arms		
BTC072350	Padded Arms - Swing Up		
BTC072000	Swing Up Arms - Short		
BTC071900	Swing Up Arms S/Steel		
BTC073600	Gutter Supports - Fitted		

POPULAR ACCESSORIES



Hand Made Open Seat Available in various sizes



Hand Made Closed Seat Available in various sizes




Pan Carrier & Bed Pan



Anti Tip Bars



Elevating Legrest Right Available in Left



IV Pole

REFERRER SIGNATURE:	
DATE: / /	
Received by Aidacare:	
1	