

Medical Carts




Medical Carts

Our carts are designed to offer an attractive, practical and durable cart with multiple configurations to enable the efficient organisation and delivery of medication and equipment.

KEY FEATURES

- Three height options: Small, Medium, Large
- Modular design enables custom drawer combinations
- Integrated handle for ease of handling
- Large work area with slide out work surface are standard
- Choice of Breakaway Lock, Digital Lock, or Key Lock
- Smooth opening drawers with soft close
- Easy roll 125mm castors
- Robust panels to minimise dents and scratches
- Many options and accessories are available



Lock System Options



BREAK-AWAY LOCK with PLASTIC SEAL
The break-away lock is simple to operate and provides quick access to cart contents in any emergency situation. The seal can be broken easily and reapplied after restocking so the cart is ready for the next emergency.



KEY LOCK
The key lock system is a proven method of securing your medical cart and providing the necessary level of security.

Cart Colours



EMERGENCY
Crash Carts




PAEDIATRIC
Crash Carts




MEDICATION
Carts



MEDICATION
SACHET Carts




ANAESTHESIA
Carts




ISOLATION
Carts


Common Configurations*



LARGE 110cm / 43"



MEDIUM 100cm / 40"



SMALL 90cm / 37"

*Diagrams show standard configurations. Custom configurations are available upon request. Consult with your representative for standard cart handle configurations. Specifications, availability & components are subject to change without notice.



DIGITAL LOCK

A user friendly system where the cart is locked and unlocked using a preset code ensuring ease of access and the security required in clinical environments. Eliminate lost key situations. The digital lock is robust and waterproof. A master code and easy reset ensures restricted access.

5
YEAR WARRANTY
STRUCTURAL

2
YEAR WARRANTY
DIGITAL LOCK

1
YEAR WARRANTY
ACCESSORIES

Emergency Crash Carts

The Aidacare Emergency Crash Cart has sleek lines, and strong construction

Smooth manoeuvrability will make it the choice in emergency departments

The simple breakaway lock offers easy access, combined with essential security

Carts come standard with all essential accessories, to enhance the user experience

Customisable to suit a range of applications

Standard Inclusions: Draw Dividers, Defibrillator Shelf, Cardiac Chest Board, Oxygen Bottle Holder, Extended Sliding Surface, Waste Bin.



Code	Description	Drawers 3" (7cm)	Drawers 6" (16cm)	Drawers 9"(25cm)	Lock Type
WAT909410	Emergency Cart - Red - Medium - 100cm / 9H	2	2	1	Break-away
WAT909420	Emergency Cart - Red - Large - 110cm / 10H	3	2	1	Break-away

Optional accessories available, see page 7 for details.

Medication Sachet Carts

The Aidacare Medication Sachet Cart offers large capacity and flexible organisation of treatment supplies in 3" and 9" seamless drawers, featuring a unique integrated divider system

Includes expansive work area, with a slide out work surface as standard

The durable steel frame and unique spread castor position ensure stability, making the cart easy to move and control

Standard Inclusions: Drawer Dividers, Extended Sliding Surface & Waste Bin.



Code	Description	Drawers 3" (7cm)	Drawers 6" (16cm)	Drawers 9"(25cm)	Lock Type
WAT909520	Medication Sachet Cart - Small - 90cm / 8H	5	0	1	Digital
WAT909570	Medication Sachet Cart - Small - 90cm / 8H	5	0	1	Key
WAT909530	Medication Sachet Cart - Medium - 100cm / 9H	6	0	1	Digital
WAT909580	Medication Sachet Cart - Medium - 100cm / 9H	6	0	1	Key

Optional accessories available, see page 7 for details.

Medication Carts

The Aidacare Medication Cart is the powerful all-round solution for a multitude of applications including Treatment, Procedure, IV, Trauma, Dressing, Supply and General cart

Carts come standard with work surface and dividers for all drawers, this cart is ready for work

Fully customisable with a host of accessories and drawer configurations

Sleek lines, strong construction and smooth manoeuvrability

Standard Inclusions: Drawer Dividers, Extended Sliding Surface & Waste Bin.



Image shows Blue Medication Cart, also available in Cream

Cream Code	Blue Code	Description	Drawers 3" (7cm)	Drawers 6" (16cm)	Drawers 9"(25cm)	Lock Type
WAT909510	WAT909430	Medication Cart - Small - 90cm / 8H	1	2	1	Digital
WAT909560	WAT909470	Medication Cart - Small - 90cm / 8H	1	2	1	Key
WAT909540	WAT909450	Medication Cart - Medium - 100cm / 9H	2	2	1	Digital
WAT909590	WAT909490	Medication Cart - Medium - 100cm / 9H	2	2	1	Key
WAT909550	WAT909460	Medication Cart - Large - 110cm / 10H	1	3	1	Digital
WAT909600	WAT909500	Medication Cart - Large - 110cm / 10H	1	3	1	Key

Optional accessories available, see page 7 for details.

Paediatric Crash Carts

Designed to accommodate the Broselow System of colour coded drawers

Ensures doctors and technicians can quickly identify critical supplies in emergencies

Fitted standard with a keyless break-away lock and comprehensive range of standard accessories

Standard Inclusions: Defibrillator Shelf, Oxygen Bottle Holder, Extended Sliding Surface, Waste Bin.



Code	Description	Drawers 3" (7cm)	Drawers 6" (16cm)	Drawers 9"(25cm)	Lock Type
WAT909690	Paediatric Cart - Large - Multi - 110cm / 10H	8	1	0	Break-away

Optional accessories available, see page 7 for details.

Anaesthesia Carts

The Aidacare Anaesthesia Cart is the flexible solution to your workflow needs

Combined with its sleek lines, strong construction and easy manoeuvrability, it offers a fully customisable solution to meet individual needs

Packed with standard accessories, it allows for choice of cart height, draw depth, locking option and storage configuration to perfectly compliment your workflow needs

Standard Inclusions: Accessory Tower, Drawer Dividers, Extended Sliding Surface and Waste Bin.



Code	Description	Drawers 3" (7cm)	Drawers 6" (16cm)	Drawers 9"(25cm)	Lock Type
WAT909640	Anaesthesia Cart - Green - Medium - 100cm / 9H	2	2	1	Digital
WAT909670	Anaesthesia Cart - Green - Medium - 100cm / 9H	2	2	1	Break-away
WAT909650	Anaesthesia Cart - Green - Large - 110cm / 10H	1	3	1	Digital
WAT909680	Anaesthesia Cart - Green - Large - 110cm / 10H	1	3	1	Break-away

Optional accessories available, see page 7 for details.

Isolation Carts

The Aidacare Isolation Cart offers all the desirable features for high-risk infection control areas

From its smooth, brightly coloured exterior, to its configurable large-volume deep drawers

Capable of supporting heavier loads with its strong robust construction, this will be your cart of choice

Carts come standard with waste bin, extended surface and dividers for all drawers

Standard Inclusions: Drawer Dividers, Extended Sliding Surface & Waste Bin.



Code	Description	Drawers 3" (7cm)	Drawers 6" (16cm)	Drawers 9"(25cm)	Lock Type
WAT909620	Isolation Cart - Yellow - Medium - 100cm / 9H	0	0	3	No Lock

Optional accessories available, see page 7 for details.

Accessories and Common Parts




Powered Drive Wheel
WAT909700



Oxygen Tank Holder
WAT907100



Waste Bin with Lid
WAT907136



Plastic Seal Lock
WAT907690



CPR Board
WAT907120



Adjustable IV Pole
WAT907130



Sharps Container Holder - WAT907135
**Sharps container not included*



Small Tilt Bin Unit 6 Compartment
WAT907155



Medium Tilt Bin Unit 5 Compartment
WAT907150



Large Tilt Bin Unit 4 Compartment
WAT907140



Glove Box Container
WAT907134



Accessory Rail
WAT907133



Defibrillator Shelf
WAT907110



Tablet Holder
WAT907205




Laptop Holder
WAT907202



Grey Drawer Dividers
Short - WAT907330



Drawer Label Plastic Tag
WAT907105




Top Surface Clear Mat
WAT907131



Drawer Tab Fillers
Grey - WAT907325
Blue - WAT907326
Yellow - WAT907327
Green - WAT907328
Red - WAT907329



150mm Castor
Non-Braking - WAT907661
Braking - WAT907662
Tracking - WAT907663



Accessory Tower
WAT909858



Clear Sub-Dividers
All Drawers:
3" - WAT908840
6" & 9" - WAT909844



Cart labels
Emergency (Red) - WAT907400
Anaesthesia (Green) - WAT907380
Isolation (Yellow) - WAT907410
Medication (Blue) - WAT907420



Medical cart shown with Aidacare powered drive wheel.