

Caeli Combination

Mattress System

Combination mattress replacement system combining 125mm low air loss alternating cells with a 50mm protective foam underlay. Ideal for pressure injury prevention in the home and community and residential aged care environments.



Caeli

aidacare.com.au • sales@aidacare.com.au • 1300 133 120



Caeli Combination

Mattress System

Features & Benefits:

- **Combination mattress** replacement with 125mm alternating air cells and 50mm foam underlay.
- Low air loss cells assist with cooling, reducing moisture build up.
- **Dual Release CPR** valves for rapid deflation in an emergency.
- **Transport cap** keeps mattress inflated for up to 24 hours.
- Carry handles on both sides.
- Fluid resistant cover with two-way stretch material.
- Digital pump with easy-to-use interface.
- **Tie-Down Straps** securely anchor the mattress to the bed surface for improved patient transfers.





Code	Description	Length	Width	Height
PTM640100	Caeli Mattress Replacement - Single	2000mm	840mm	175mm
PTM640110	Caeli Mattress Replacement - King Single	2000mm	1080mm	175mm

